

# 5009

PRODUCT RANGE

wire pins... mini plates & screws... angled plates... cannulated screws... headless screws  
locking plates... intramedullary nails... spinal systems... dhs/dcs plates... elastic nails... wire  
locking plates... locking screws... external fixators... angled plates... trochanteric plates

quality & trust

*Implant Catalog*

[www.aysam.com.tr](http://www.aysam.com.tr)





**“We supply WORLD’s  
Orthopedic & Trauma  
Implants...”**





# Implant Products

We continue to add new products to our 5009 trauma products

01

## **SMALL/LARGE FRAGMENT LOCKING/NON LOCKING PLATES & SCREWS**

These products, components, treatment kits and tools are manufactured with medical interventions for the treatment of patients.

02

## **FEMUR/HUMERUS/TIBIA INTRAMEDULLARY NAILS**

All nail implants, according to the type, model, size (details are given in Aysam catalog) are used in the treatment of bone fractures, congenital or later acquired disorders of skeletal system. Surgeons select the right type, size of implants according to conditions of patients.

03

## **POSTERIOR SPINAL SYSTEMS, CERVICAL PLATES & SCREWS**

All Spinal System implants; type, model, size (details are given on Aysam catalog) and according to the surgeon's choice are used in the treatment of bone fractures, congenital or later acquired orthopedic, traumatological disorders of spinal column.

04

## **AYFIX DYNAMIC AXIAL, MINI RAIL, TUBULAR, HYBRID, CIRCULAR EXTERNAL FIXATORS**

External fixator implants and components are used for treatment of Orthopedic & Traumatologic disorders for fixation, distraction & compression purposes of different bones & joints in human body.

05

## **WIRES, PINS, STAPLES & DRILL BITS**

The wire-pin implants are designed and implemented to perform bone fixation in fractures. Fixation in corrective osteotomies, pseudoarthroses, arthrodesis, joint dislocations, replantation treatments, bone reconstructive & skeletal traction procedures. All type, model, size of implants are detailed in Aysam catalog.



## **OUR VISION IS TO BE THE LEADING COMPANY WORLD WIDE IN ORTHOPAEDIC TRAUMA WITH SUSTAINABLE GROWTH, BY FOLLOWING “CONTINUOUS IMPROVEMENT” AND “QUALITY IS IN DETAILS” PRINCIPLE.**

Since our establishment in 2001, we have continued to develop and produce with a focus on knowledge, technology, and human and animal health. While progressing with the awareness that our primary responsibility is to our patients and healthcare workers, we combine our knowledge and experience with modern engineering methods to help millions of patients heal and regain their mobility. We reach patients all over the world through our exports to over +100 countries across 6 continents, as well as through public hospitals, veterinary clinics, and distributors domestically. Emphasizing the importance of accessible healthcare for all, we ensure our products reach many aid organizations both domestically and internationally. Through our R&D and P&D activities, we enhance efficiency and contribute to the development and progress of our country, while continuing to support all projects endorsed by TÜBİTAK, universities, and other institutions with our products.

Aware that our responsibility extends beyond human and animal health, we adopt a sustainable production approach in all our production processes with the vision of “living well within planetary boundaries,” ensuring that we do not harm our country’s resources or the environment.

At Aysam, we prioritize Total Quality Management. Our quality management, based on principles such as customer satisfaction, leadership, collaboration, and participation, aims to meet our customers’ expectations, minimize errors, increase efficiency, and most importantly, continuously educate and develop our employees. We acknowledge that no process is flawless, from raw materials to all production and post-sales processes, and we embrace continuous improvement as our quality policy.

We host a total of 26,961 products, including 14,170 orthopedic trauma products, 285 total knee and hip arthroplasty, 741 sterilization containers and accessories, perforated sterilization baskets, 410 maxillofacial (implant) products, and 11,355 orthopedic surgical instruments and orthopedic surgical instrument sets. With our ongoing R&D and P&D efforts, we continue to expand our product portfolio each year.

We take pride in offering various solutions, from our doctors to our patients, with new developments in the field of orthopedic devices and medical technology, and we continue our work with the same dedication and enthusiasm as on our first day.

### **Quality Management**

In all our processes, being aware that we work in the field of human and animal health, we carry out production and post-production processes in compliance with CE, ISO 9001, ISO 13485 standards, and the 93/42/EEC Medical Devices Directive. We ensure the safe and effective production and distribution of our products, enabling the use of high-tech devices in the healthcare sector.

## QUALITY

We produce best quality products according to ISO 13485/MDD.

## TRUST

We always support our customers and distributors in the market.

## EXPERIENCE

Aysam has over 20 years experience on orthopaedic implant production.

## INNOVATION

Innovation is our indispensable condition to move forward.



f aysam.orthopedics

t AysamOrtho

in aysam

YouTube aysamorthopaedics

Instagram aysamorthopaedics

# INDEX

## SMALL FRAGMENT LOCKING & NON-LOCKING PLATES

- Small Fragment Locking Plates \_\_\_\_\_ 007
- Small Fragment Non-Locking Plates \_\_\_\_\_ 023

## SMALL FRAGMENT LOCKING & NON-LOCKING SCREWS

- Small Fragment Locking Screws \_\_\_\_\_ 031
- Small Fragment Non-Locking Screws \_\_\_\_\_ 035

## LARGE FRAGMENT LOCKING & NON-LOCKING PLATES

- Large Fragment Locking Plates \_\_\_\_\_ 039
- Large Fragment Non-Locking Plates \_\_\_\_\_ 054

## DHS - DCS PLATES & LAG-COMPRESSION SCREWS

- DHS-DCS Plates \_\_\_\_\_ 063
- DHS-DCS Lag & Compression Screws \_\_\_\_\_ 066

## LARGE FRAGMENT LOCKING & NON-LOCKING SCREWS

- Large Fragment Locking Screws \_\_\_\_\_ 067
- Large Fragment Non-Locking Screws \_\_\_\_\_ 073

## MINI FRAGMENT LOCKING PLATES & SCREWS

- 1.5 mm Mini Fragment Locking Plates and Screws \_\_\_\_\_ 077
- 2.0 mm Mini Fragment Locking Plates and Screws \_\_\_\_\_ 080
- 2.4 mm Mini Fragment Locking Plates and Screws \_\_\_\_\_ 083
- 2.7 mm Mini Fragment Locking Plates and Screws \_\_\_\_\_ 085

## CHEST WALL PLATES, SCREWS & SPLINTS

- Chest Wall Locking Plates \_\_\_\_\_ 89
- Chest Wall Locking Screws \_\_\_\_\_ 92
- Rib Intramedullary Splints \_\_\_\_\_ 92

## ANGLED BLADE PLATES & PEDIATRIC PLATES

- Angled Blade Plates \_\_\_\_\_ 93
- Pediatric Locking Plates \_\_\_\_\_ 96
- Pediatric Eight Plates \_\_\_\_\_ 96

## TROCHANTERIC PLATES & CABLE SYSTEMS

- Straight Cable Grip Plates \_\_\_\_\_ 97
- Trochanteric Cable Hook Plates \_\_\_\_\_ 98
- Trochanteric Cable Hook Grip Plates \_\_\_\_\_ 99
- Cable Locking Screws \_\_\_\_\_ 100
- Cable-Crimp Blocks \_\_\_\_\_ 100

## CANNULATED SCREWS

- Cannulated Screws \_\_\_\_\_ 101

## HEADLESS & HERBERT SCREWS

- Headless Screws \_\_\_\_\_ 109
- Herbert Screws \_\_\_\_\_ 111
- Interference Screws \_\_\_\_\_ 112

## INTRAMEDULLARY NAILS

### Humerus Intramedullary Nails

- Humerus Intramedullary Nails \_\_\_\_\_ 114

### Femur Intramedullary Nails

- Femur Intramedullary Nails ( ANTEGRADE ) \_\_\_\_\_ 115
- Femur Intramedullary Nails ( DFN/RETROGRADE ) \_\_\_\_\_ 116
- Femur Short Reconstruction Intramedullary Nails (PFN) \_\_\_\_\_ 116
- Femur Long Reconstruction Intramedullary Nails \_\_\_\_\_ 117
- APFNA Femur Short Intramedullary Nails \_\_\_\_\_ 118
- AROCOM Proximal Femur Antirotation  
Compression Intramedullary Nails \_\_\_\_\_ 119
- ANTIRO Proximal Femur Antirotation  
Compression Intramedullary Nails \_\_\_\_\_ 121

### Tibia Intramedullary Nails

- Tibia Intramedullary Nails \_\_\_\_\_ 123
- ANT Tibia Intramedullary Nails \_\_\_\_\_ 124

## ELASTIC NAILS \_\_\_\_\_ 125

## SPINAL SYSTEMS

- Posterior Spinal Systems \_\_\_\_\_ 127
- Cervical Plates & Screws \_\_\_\_\_ 132
- Anterior Thoracic-Lumbar Plates, Bolts & Screws \_\_\_\_\_ 134

## EXTERNAL FIXATORS

- Ayfix Dynamic Axial External Fixators \_\_\_\_\_ 135
- Mini Rail External Fixators \_\_\_\_\_ 148
- Tubular External Fixators \_\_\_\_\_ 150
- Hybrid External Fixators \_\_\_\_\_ 154
- Circular External Fixators \_\_\_\_\_ 159
- Pediatric Circular External Fixators \_\_\_\_\_ 166

## WIRES, PINS, STAPLES & DRILL BITS

- Kirschner Wires \_\_\_\_\_ 172
- Steinmann Pins \_\_\_\_\_ 173
- Schanz Screws \_\_\_\_\_ 174
- Rush Pins \_\_\_\_\_ 175
- Staples \_\_\_\_\_ 176
- Tibial Bolts \_\_\_\_\_ 177
- Knowless Hip Pins \_\_\_\_\_ 177
- Drill Bits \_\_\_\_\_ 178
- Gigli Saw \_\_\_\_\_ 178

"We supply WORLD's  
Orthopedic & Trauma  
Implants..."

# Small Fragment Plates

Locking/Non-Locking




## Small Fragment Locking Plates

### 3,5 mm Clavicle Locking Compression Plate



| SS       | TI       | HOLES | SIDE  | LENGTH (mm) |
|----------|----------|-------|-------|-------------|
| A1031906 | A1031806 | 6     | Left  | 70          |
| A1031907 | A1031807 | 7     | Left  | 84          |
| A1031908 | A1031808 | 8     | Left  | 98          |
| A1031909 | A1031809 | 9     | Left  | 112         |
| A1042906 | A1042806 | 6     | Right | 70          |
| A1042907 | A1042807 | 7     | Right | 84          |
| A1042908 | A1042808 | 8     | Right | 98          |
| A1042909 | A1042809 | 9     | Right | 112         |




**Intended use :** Temporary fixation of fractures, malunions, nonunions, osteotomies and pseudoarthroses of the clavicle, particularly for patients with osteopenic bone.

### 3,5 mm Clavicle Hook Locking Plate



| SS        | TI        | HOLES | SIDE  | LENGTH (mm) |
|-----------|-----------|-------|-------|-------------|
| A26009141 | A26009041 | 4     | Left  | 63          |
| A26009151 | A26009051 | 5     | Left  | 75          |
| A26009161 | A26009061 | 6     | Left  | 86          |
| A26009171 | A26009071 | 7     | Left  | 97          |
| A26009142 | A26009042 | 4     | Right | 63          |
| A26009152 | A26009052 | 5     | Right | 75          |
| A26009162 | A26009062 | 6     | Right | 86          |
| A26009172 | A26009072 | 7     | Right | 97          |




**Intended use :** Intended use: Temporary fixation of fractures, malunions, nonunions, osteotomies and pseudoarthroses of the clavicle, for the treatment of acromioclavicular dislocations, particularly for patients with osteopenic bone.

# Small Fragment Locking Plates

## 2,7/3,5 mm Distal Clavicle Locking Compression Plate



| SS       | TI       | HOLES | SIDE  | LENGTH (mm) |
|----------|----------|-------|-------|-------------|
| A5521903 | A1021803 | 3     | Left  | 66,5        |
| A5521904 | A1021804 | 4     | Left  | 78,7        |
| A5521905 | A1021805 | 5     | Left  | 91,1        |
| A5521906 | A1021806 | 6     | Left  | 103,2       |
| A5521907 | A1021807 | 7     | Left  | 113,8       |
| A5522903 | A1022803 | 3     | Right | 66,5        |
| A5522904 | A1022804 | 4     | Right | 78,7        |
| A5522905 | A1022805 | 5     | Right | 91,1        |
| A5522906 | A1022806 | 6     | Right | 103,2       |
| A5522907 | A1022807 | 7     | Right | 113,8       |




**Intended use :** Temporary fixation of fractures, malunions, nonunions, osteotomies and pseudoarthroses of the clavicle, for the treatment of acromioclavicular dislocations, particularly for patients with osteopenic bone.

## 3,5 mm Clavicle Reconstruction Locking Compression Plate



| SS       | TI       | HOLES | SIDE  | LENGTH (mm) |
|----------|----------|-------|-------|-------------|
| A1531806 | A1530806 | 6     | Left  | 94,3        |
| A1531807 | A1530807 | 7     | Left  | 107,8       |
| A1531813 | A1530813 | 8     | Left  | 121         |
| A1531809 | A1530809 | 9     | Left  | 133         |
| A1531906 | A1530906 | 6     | Right | 94,3        |
| A1531907 | A1530907 | 7     | Right | 107,8       |
| A1531913 | A1530913 | 8     | Right | 121         |
| A1531910 | A1530910 | 9     | Right | 133         |




**Intended use :** Temporary fixation of fractures, malunions, nonunions, osteotomies and pseudoarthroses of the clavicle, particularly for patients with osteopenic bone.

## Small Fragment Locking Plates

### 3,5 mm Proximal Humerus Locking Compression Plate 1



| SS       | TI       | HOLES | LENGTH (mm) |
|----------|----------|-------|-------------|
| A1620903 | A1620803 | 3     | 61,1        |
| A1620904 | A1620804 | 4     | 73,1        |
| A1620905 | A1620805 | 5     | 85,1        |
| A1620906 | A1620806 | 6     | 97,1        |
| A1620907 | A1620807 | 7     | 109,1       |
| A1620908 | A1620808 | 8     | 121,1       |




**Intended use :** Temporary fixation of fractures and fracture dislocations, malunions, nonunions, osteotomies and pseudoarthroses of the proximal humerus, particularly for adult patients with osteopenic bone.

### 3,5 mm Proximal Humerus Locking Compression Plate 2



| SS       | TI       | HOLES | LENGTH (mm) |
|----------|----------|-------|-------------|
| A1630903 | A1630803 | 3     | 87          |
| A1630904 | A1630804 | 4     | 99          |
| A1630905 | A1630805 | 5     | 111         |
| A1630906 | A1630806 | 6     | 123         |
| A1630907 | A1630807 | 7     | 135         |
| A1630908 | A1630808 | 8     | 147         |
| A1630909 | A1630809 | 9     | 159         |




**Intended use :** Temporary fixation of fractures and fracture dislocations, malunions, nonunions, osteotomies and pseudoarthroses of the proximal humerus, particularly for adult patients with osteopenic bone.

# Small Fragment Locking Plates

## 3,5 mm Proximal Humerus Locking Compression Plate 3



| SS       | TI       | HOLES | LENGTH (mm) |
|----------|----------|-------|-------------|
| A1640903 | A1640803 | 3     | 122         |
| A1640904 | A1640804 | 4     | 140         |
| A1640905 | A1640805 | 5     | 158         |
| A1640906 | A1640806 | 6     | 176         |
| A1640907 | A1640807 | 7     | 194         |




**Intended use :** Temporary fixation of fractures and fracture dislocations, malunions, nonunions, osteotomies and pseudoarthroses of the proximal humerus, particularly for adult patients with osteopenic bone.

## 3,5 mm Cloverleaf Locking Plate




| SS               | TI               | HOLES | LENGTH (mm) |
|------------------|------------------|-------|-------------|
| A101 02 270 0003 | A111 02 270 0003 | 3     | 88,1        |
| A101 02 270 0004 | A111 02 270 0004 | 4     | 104,1       |
| A101 02 270 0005 | A111 02 270 0005 | 5     | 120,1       |
| A101 02 270 0006 | A111 02 270 0006 | 6     | 136,1       |
| A101 02 270 0007 | A111 02 270 0007 | 7     | 152,1       |
| A101 02 270 0008 | A111 02 270 0008 | 8     | 168,1       |



**Intended use :** Temporary fixation of intra-articular fractures of the distal tibia and proximal humerus fractures in adult patients.

# Small Fragment Locking Plates

## 3,5 mm Distal Lateral Humerus Locking Plate




| SS       | TI       | HOLES | SIDE  | LENGTH (mm) |
|----------|----------|-------|-------|-------------|
| A5248105 | A5238105 | 5     | Left  | 56,5        |
| A5248107 | A5238107 | 7     | Left  | 83,5        |
| A5248109 | A5238109 | 9     | Left  | 108,5       |
| A5248111 | A5238111 | 11    | Left  | 133,5       |
| A5248113 | A5238113 | 13    | Left  | 158,5       |
| A5248115 | A5238115 | 15    | Left  | 183,5       |
| A5248205 | A5238205 | 5     | Right | 56,5        |
| A5248207 | A5238207 | 7     | Right | 83,5        |
| A5248209 | A5238209 | 9     | Right | 108,5       |
| A5248211 | A5238211 | 11    | Right | 133,5       |
| A5248213 | A5238213 | 13    | Right | 158,5       |
| A5248215 | A5238215 | 15    | Right | 183,5       |


**Intended use :** Temporary fixation of fractures, malunions, nonunions, osteotomies and pseudoarthroses of the distal humerus (lateral approach), particularly for adult patients with osteopenic bone.

## 2,7/3,5 mm Distal Humerus Dorsolateral Locking Plate

with lateral support



without lateral support



| SS       | TI       | HOLES | SIDE  | LENGTH (mm) |
|----------|----------|-------|-------|-------------|
| A9248105 | A9238105 | 5     | Left  | 103,9       |
| A9248107 | A9238107 | 7     | Left  | 126,9       |
| A9248109 | A9238109 | 9     | Left  | 152,9       |
| A9248111 | A9238111 | 11    | Left  | 178,9       |
| A9248113 | A9238113 | 13    | Left  | 204,9       |
| A9248205 | A9238205 | 5     | Right | 103,9       |
| A9248207 | A9238207 | 7     | Right | 126,9       |
| A9248209 | A9238209 | 9     | Right | 152,9       |
| A9248211 | A9238211 | 11    | Right | 178,9       |
| A9248213 | A9238213 | 13    | Right | 204,9       |


**Intended use :** Temporary fixation of fractures, malunions, nonunions, osteotomies and pseudoarthroses of the distal humerus (dorsolateral approach), particularly for adult patients with osteopenic bone.

## Small Fragment Locking Plates

### 3,5 mm Distal Medial Humerus Locking Plate 1



| SS       | TI       | HOLES | LENGTH (mm) |
|----------|----------|-------|-------------|
| A5249206 | A5249106 | 6     | 69,5        |
| A5249208 | A5249108 | 8     | 94          |
| A5249210 | A5249110 | 10    | 114         |
| A5249212 | A5249112 | 12    | 138,5       |
| A5249214 | A5249114 | 14    | 163         |
| A5249216 | A5249116 | 16    | 184         |




**Intended use :** Temporary fixation of fractures, malunions, nonunions, osteotomies and pseudoarthroses of the distal humerus (medial approach), particularly for adult patients with osteopenic bone.

### 2,7/3,5 mm Distal Medial Humerus Locking Plate 2



| SS          | TI          | HOLES | SIDE  | LENGTH (mm) |
|-------------|-------------|-------|-------|-------------|
| A1015279106 | A1115278106 | 6     | Left  | 84,8        |
| A1015279108 | A1115278108 | 8     | Left  | 109,3       |
| A1015279110 | A1115278110 | 10    | Left  | 133,8       |
| A1015279112 | A1115278112 | 12    | Left  | 158,3       |
| A1015279114 | A1115278114 | 14    | Left  | 182,8       |
| A1015279116 | A1115278116 | 16    | Left  | 207,3       |
| A1015279206 | A1115278206 | 6     | Right | 84,8        |
| A1015279208 | A1115278208 | 8     | Right | 109,3       |
| A1015279210 | A1115278210 | 10    | Right | 133,8       |
| A1015279212 | A1115278212 | 12    | Right | 158,3       |
| A1015279214 | A1115278214 | 14    | Right | 182,8       |
| A1015279216 | A1115278216 | 16    | Right | 207,3       |




**Intended use :** Temporary fixation of fractures, malunions, nonunions, osteotomies and pseudoarthroses of the distal humerus (medial approach), particularly for adult patients with osteopenic bone.

## Small Fragment Locking Plates

### 3,5 mm Y-Distal Humerus Locking Plate 1



| SS       | TI       | HOLES | SIDE  | LENGTH (mm) |
|----------|----------|-------|-------|-------------|
| A5319103 | A5318103 | 3     | Left  | 91,4        |
| A5319104 | A5318104 | 4     | Left  | 106,4       |
| A5319105 | A5318105 | 5     | Left  | 121,4       |
| A5319106 | A5318106 | 6     | Left  | 136,4       |
| A5319107 | A5318107 | 7     | Left  | 151,4       |
| A5319108 | A5318108 | 8     | Left  | 166,4       |
| A5319203 | A5318203 | 3     | Right | 91,4        |
| A5319204 | A5318204 | 4     | Right | 106,4       |
| A5319205 | A5318205 | 5     | Right | 121,4       |
| A5319206 | A5318206 | 6     | Right | 136,4       |
| A5319207 | A5318207 | 7     | Right | 151,4       |
| A5319208 | A5318208 | 8     | Right | 166,4       |




**Intended use :** Temporary fixation of fractures, malunions, nonunions, osteotomies and pseudoarthroses of the distal humerus, particularly for adult patients with osteopenic bone.

### 3,5 mm Olecranon Locking Compression Plate 1



| SS       | TI       | HOLES | SIDE  | LENGTH (mm) |
|----------|----------|-------|-------|-------------|
| A1651904 | A1651804 | 4     | Left  | 128,8       |
| A1651906 | A1651806 | 6     | Left  | 153,7       |
| A1651908 | A1651808 | 8     | Left  | 179,8       |
| A1662904 | A1662804 | 4     | Right | 128,8       |
| A1662906 | A1662806 | 6     | Right | 153,7       |
| A1662908 | A1662808 | 8     | Right | 179,8       |




**Intended use :** Temporary fixation of fractures, malunions, nonunions, osteotomies and pseudoarthroses of the proximal ulna, particularly for adult patients with osteopenic bone.

## Small Fragment Locking Plates

### 3,5 mm Olecranon Locking Compression Plate 2



| SS        | TI        | HOLES | SIDE  | LENGTH (mm) |
|-----------|-----------|-------|-------|-------------|
| A22024130 | A22014130 | 3     | Left  | 67,2        |
| A22024140 | A22014140 | 4     | Left  | 79,6        |
| A22024150 | A22014150 | 5     | Left  | 92,1        |
| A22024160 | A22014160 | 6     | Left  | 104,6       |
| A22024170 | A22014170 | 7     | Left  | 117         |
| A22024180 | A22014180 | 8     | Left  | 129,4       |
| A22024230 | A22014230 | 3     | Right | 67,2        |
| A22024240 | A22014240 | 4     | Right | 79,6        |
| A22024250 | A22014250 | 5     | Right | 92,1        |
| A22024260 | A22014260 | 6     | Right | 104,6       |
| A22024270 | A22014270 | 7     | Right | 117         |
| A22024280 | A22014280 | 8     | Right | 129,4       |




**Intended use :** Temporary fixation of fractures, malunions, nonunions, osteotomies and pseudoarthroses of the proximal ulna, particularly for adult patients with osteopenic bone.

### 3,5 mm Metaphysis Locking Compression Plate



| SS       | TI       | HOLES | LENGTH (mm) |
|----------|----------|-------|-------------|
| A1060906 | A1060806 | 6     | 86          |
| A1060907 | A1060807 | 7     | 99          |
| A1060908 | A1060808 | 8     | 112         |
| A1060909 | A1060809 | 9     | 125         |
| A1060911 | A1060811 | 11    | 151         |
| A1060913 | A1060813 | 13    | 177         |




**Intended use :** Temporary fixation of fractures, malunions, nonunions, osteotomies and pseudoarthroses of the various long bones, excluding the femur in adult patients, particularly for patients with osteopenic bone.

## Small Fragment Locking Plates

### 3,5 mm Locking Compression Hook Plate



| SS       | TI       | HOLES | LENGTH (mm) |
|----------|----------|-------|-------------|
| A1610905 | A1610805 | 5     | 64,4        |
| A1610906 | A1610806 | 6     | 72,9        |
| A1610907 | A1610807 | 7     | 81,4        |
| A1610908 | A1610808 | 8     | 89,9        |
| A1610909 | A1610809 | 9     | 98,4        |
| A1610910 | A1610810 | 10    | 106,9       |
| A1610911 | A1610811 | 11    | 115,4       |
| A1610912 | A1610812 | 12    | 123,9       |
| A1610913 | A1610813 | 13    | 132,4       |




**Intended use :** Temporary fixation of fractures, malunions, nonunions, osteotomies and pseudoarthroses of the ulna, radius and fibula, particularly for adult patients with osteopenic bone.

### 3,5 mm Locking Compression Plate




| SS       | TI       | HOLES | LENGTH (mm) |
|----------|----------|-------|-------------|
| A1840904 | A1820804 | 4     | 58          |
| A1840905 | A1820805 | 5     | 72          |
| A1840906 | A1820806 | 6     | 86          |
| A1840907 | A1820807 | 7     | 100         |
| A1840908 | A1820808 | 8     | 114         |
| A1840909 | A1820809 | 9     | 128         |
| A1840910 | A1820810 | 10    | 142         |
| A1840911 | A1820811 | 11    | 156         |
| A1840912 | A1820812 | 12    | 170         |



**Intended use :** Temporary fixation of fractures, malunions, nonunions, osteotomies and pseudoarthroses of the clavicle, scapula, olecranon, humerus, radius, ulna, pelvis, distal tibia, fibula, particularly in osteopenic bone for adult patients. These plates are also indicated for fracture fixation of diaphyseal and metaphyseal areas of long bones in pediatric patients.

# Small Fragment Locking Plates


## 1/3 Tubular Locking Plate



| SS       | TI       | HOLES | LENGTH (mm) |
|----------|----------|-------|-------------|
| A5051004 | A5050004 | 4     | 54          |
| A5051005 | A5050005 | 5     | 66          |
| A5051006 | A5050006 | 6     | 78          |
| A5051007 | A5050007 | 7     | 90          |
| A5051008 | A5050008 | 8     | 102         |
| A5051009 | A5050009 | 9     | 114         |
| A5051010 | A5050010 | 10    | 126         |
| A5051011 | A5050011 | 11    | 138         |
| A5051012 | A5050012 | 12    | 150         |
| A5051013 | A5050013 | 13    | 162         |

**Intended use :** Temporary fixation of fractures, malunions, nonunions, osteotomies and pseudoarthroses of small bones, including the ulna, radius and fibula, particularly for patients with osteopenic bone.

## 3,5 mm Reconstruction Locking Plate



| SS       | TI       | HOLES | LENGTH (mm) |
|----------|----------|-------|-------------|
| A1050904 | A1050804 | 4     | 48          |
| A1050905 | A1050805 | 5     | 60          |
| A1050906 | A1050806 | 6     | 72          |
| A1050907 | A1050807 | 7     | 84          |
| A1050908 | A1050808 | 8     | 96          |
| A1050909 | A1050809 | 9     | 108         |
| A1050910 | A1050810 | 10    | 120         |
| A1050911 | A1050811 | 11    | 132         |
| A1050912 | A1050812 | 12    | 144         |


**Intended use :** Temporary fixation of fractures, malunions, nonunions, osteotomies and pseudoarthroses of the clavicle, scapula, olecranon, humerus, radius, ulna, pelvis, distal tibia, fibula for adult patients. These plates are also indicated for fracture fixation of diaphyseal and metaphyseal areas of long bones in pediatric patients.

## Small Fragment Locking Plates

### 3,5 mm Ulna/Radius Anatomic Locking Plate



| SS       | TI       | HOLES | LENGTH (mm) |
|----------|----------|-------|-------------|
| A1822005 | A1821005 | 5     | 69          |
| A1822006 | A1821006 | 6     | 82          |
| A1822007 | A1821007 | 7     | 95,4        |
| A1822008 | A1821008 | 8     | 108,8       |
| A1822009 | A1821009 | 9     | 121,2       |




**Intended use :** Temporary fixation of fractures, malunions, nonunions, osteotomies and pseudoarthroses of the midshaft ulna and radius, particularly for patients with osteopenic bone.

### 3,5 mm Curved Reconstruction Locking Plate



| SS               | TI               | HOLES | LENGTH (mm) |
|------------------|------------------|-------|-------------|
| A101 02 221 0004 | A111 02 221 0004 | 4     | 61,4        |
| A101 02 221 0005 | A111 02 221 0005 | 5     | 77,8        |
| A101 02 221 0006 | A111 02 221 0006 | 6     | 93,7        |
| A101 02 221 0007 | A111 02 221 0007 | 7     | 109,4       |
| A101 02 221 0008 | A111 02 221 0008 | 8     | 125,0       |
| A101 02 221 0009 | A111 02 221 0009 | 9     | 140,2       |
| A101 02 221 0010 | A111 02 221 0010 | 10    | 155,0       |
| A101 02 221 0011 | A111 02 221 0011 | 11    | 169,9       |
| A101 02 221 0012 | A111 02 221 0012 | 12    | 184,7       |
| A101 02 221 0013 | A111 02 221 0013 | 13    | 198,0       |
| A101 02 221 0014 | A111 02 221 0014 | 14    | 212,0       |




**Intended use :** Temporary fixation of fractures, malunions, nonunions, osteotomies and pseudoarthroses of the pelvis, particularly for patients with osteopenic bone.

## Small Fragment Locking Plates

### 2,7/3,5 mm Distal Radius Locking Plate 1



| SS       | TI       | HOLES | SIDE  | LENGTH (mm) |
|----------|----------|-------|-------|-------------|
| A5159303 | A5159103 | 3     | Left  | 54          |
| A5159304 | A5159104 | 4     | Left  | 62          |
| A5159305 | A5159105 | 5     | Left  | 70          |
| A5159306 | A5159106 | 6     | Left  | 78          |
| A5159307 | A5159107 | 7     | Left  | 86          |
| A5159403 | A5159203 | 3     | Right | 54          |
| A5159404 | A5159204 | 4     | Right | 62          |
| A5159405 | A5159205 | 5     | Right | 70          |
| A5159406 | A5159206 | 6     | Right | 78          |
| A5159407 | A5159207 | 7     | Right | 86          |




**Intended use :** Temporary fixation of complex intra-articular and extra-articular fractures and osteotomies of the distal radius in adults, skeletally mature adolescents.

### 2,7/3,5 mm Distal Radius Locking Plate 2



| SS       | TI       | HOLES | SIDE  | LENGTH (mm) |
|----------|----------|-------|-------|-------------|
| A5169303 | A5169103 | 3     | Left  | 63,2        |
| A5169304 | A5169104 | 4     | Left  | 71,2        |
| A5169305 | A5169105 | 5     | Left  | 79,2        |
| A5169306 | A5169106 | 6     | Left  | 87,2        |
| A5169307 | A5169107 | 7     | Left  | 95,2        |
| A5169403 | A5169203 | 3     | Right | 63,2        |
| A5169404 | A5169204 | 4     | Right | 71,2        |
| A5169405 | A5169205 | 5     | Right | 79,2        |
| A5169406 | A5169206 | 6     | Right | 87,2        |
| A5169407 | A5169207 | 7     | Right | 95,2        |




**Intended use :** Temporary fixation of complex intra-articular and extra-articular fractures and osteotomies of the distal radius in adults, skeletally mature adolescents.

## Small Fragment Locking Plates

### 3,5 mm Locking Compression T-Plate, Right-Angled (3 Holes on Head)



| SS       | TI       | HOLES | LENGTH (mm) |
|----------|----------|-------|-------------|
| A1480903 | A1480803 | 3     | 49          |
| A1480904 | A1480804 | 4     | 61          |
| A1480905 | A1480805 | 5     | 73          |
| A1480906 | A1480806 | 6     | 85          |




**Intended use :** Temporary fixation of complex intra-articular and extra-articular fractures and osteotomies of the distal radius in adults, skeletally mature adolescents.

### 3,5 mm Locking Compression T-Plate Oblique-Angled



| SS       | TI       | HOLES | SIDE  | LENGTH (mm) |
|----------|----------|-------|-------|-------------|
| A1511903 | A1511803 | 3     | Left  | 53          |
| A1511904 | A1511804 | 4     | Left  | 66,5        |
| A1511905 | A1511805 | 5     | Left  | 80          |
| A1511906 | A1511806 | 6     | Left  | 93,5        |
| A1522903 | A1522803 | 3     | Right | 53          |
| A1522904 | A1522804 | 4     | Right | 66,5        |
| A1522905 | A1522805 | 5     | Right | 80          |
| A1522906 | A1522806 | 6     | Right | 93,5        |




**Intended use :** Temporary fixation of complex intra-articular and extra-articular fractures and osteotomies of the distal radius in adults, skeletally mature adolescents.

## Small Fragment Locking Plates

### 3,5 mm Locking Compression T-Plate, Right-Angled (4 Holes on Head)



| SS       | TI       | HOLES | LENGTH (mm) |
|----------|----------|-------|-------------|
| A1480930 | A1480830 | 3     | 46          |
| A1480940 | A1480840 | 4     | 57          |
| A1480950 | A1480850 | 5     | 68          |
| A1480960 | A1480860 | 6     | 80          |




**Intended use :** Temporary fixation of complex intra-articular and extra-articular fractures and osteotomies of the distal radius in adults, skeletally mature adolescents.

### 2,7 mm Radius Head Locking Plate



| SS       | TI       | HOLES | LENGTH (mm) |
|----------|----------|-------|-------------|
| A5159603 | A5159503 | 3     | 33          |
| A5159605 | A5159505 | 5     | 47          |




**Intended use :** Temporary fixation of extra-articular and intra-articular fractures of the proximal radius and multifragmented radial neck fractures in adults, skeletally mature adolescents.

## Small Fragment Locking Plates

### 2,7 mm Hallux-Valgus Locking Plate



| TI               | HOLES | SIDE       | LENGTH (mm) |
|------------------|-------|------------|-------------|
| A312 02 062 0402 | 2     | 90° Left   | 33          |
| A312 02 062 0403 | 3     | 90° Left   | 35          |
| A312 02 062 0404 | 4     | 90° Left   | 37          |
| A312 02 062 0202 | 2     | 110° Left  | 33          |
| A312 02 062 0203 | 3     | 110° Left  | 35          |
| A312 02 062 0204 | 4     | 110° Left  | 37          |
| A312 02 061 0302 | 2     | 90° Right  | 33          |
| A312 02 061 0303 | 3     | 90° Right  | 35          |
| A312 02 061 0304 | 4     | 90° Right  | 37          |
| A312 02 061 0102 | 2     | 110° Right | 33          |
| A312 02 061 0103 | 3     | 110° Right | 35          |
| A312 02 061 0104 | 4     | 110° Right | 37          |




**Intended use :** Temporary fixation of hallux valgus osteotomies of adult patients.

### Attachment Locking Plate Small



| TI               | TYPE  |
|------------------|-------|
| A110 02 190 0006 | Small |




**Intended use :** For augmented temporary fixation of fractures, used as an attachment of 3.5 mm locking compression plates.



## Small Fragment Non-Locking Plates


### 3,5 mm Clavicle Plate



**Intended use :** Temporary fixation of fractures, malunions, nonunions, osteotomies and pseudoarthroses of the clavicle.

| SS       | TI       | HOLES | SIDE  | LENGTH (mm) |
|----------|----------|-------|-------|-------------|
| A5089106 | A5088106 | 6     | Left  | 70          |
| A5089107 | A5088107 | 7     | Left  | 84          |
| A5089108 | A5088108 | 8     | Left  | 98          |
| A5089109 | A5088109 | 9     | Left  | 112         |
| A5089206 | A5088206 | 6     | Right | 70          |
| A5089207 | A5088207 | 7     | Right | 84          |
| A5089208 | A5088208 | 8     | Right | 98          |
| A5089209 | A5088209 | 9     | Right | 112         |

### 3,5 mm Clavicle Hook Plate



**Intended use :** Temporary fixation of fractures, malunions, nonunions, osteotomies and pseudoarthroses of the clavicle, for the treatment of acromioclavicular dislocations


| SS       | TI       | HOLES | SIDE  | LENGTH (mm) |
|----------|----------|-------|-------|-------------|
| A5269104 | A5268104 | 4     | Left  | 63          |
| A5269105 | A5268105 | 5     | Left  | 75          |
| A5269106 | A5268106 | 6     | Left  | 86          |
| A5269107 | A5268107 | 7     | Left  | 97          |
| A5269204 | A5268204 | 4     | Right | 63          |
| A5269205 | A5268205 | 5     | Right | 75          |
| A5269206 | A5268206 | 6     | Right | 86          |
| A5269207 | A5268207 | 7     | Right | 97          |

## Small Fragment Non-Locking Plates

### 3,5 mm Cloverleaf Plate



| SS               | TI               | HOLES | LENGTH (mm) |
|------------------|------------------|-------|-------------|
| A100 02 270 0003 | A110 02 270 0003 | 3     | 88,1        |
| A100 02 270 0004 | A110 02 270 0004 | 4     | 104,1       |
| A100 02 270 0005 | A110 02 270 0005 | 5     | 120,1       |
| A100 02 270 0006 | A110 02 270 0006 | 6     | 136,1       |
| A100 02 270 0007 | A110 02 270 0007 | 7     | 152,1       |
| A100 02 270 0008 | A110 02 270 0008 | 8     | 168,1       |




**Intended use :** Temporary fixation of intra-articular fractures of the distal tibia and proximal humerus fractures in adult patients.

### 3,5 mm Y-Distal Humerus Plate 1




| SS       | TI       | HOLES | SIDE  | LENGTH (mm) |
|----------|----------|-------|-------|-------------|
| A5219103 | A5218103 | 3     | Left  | 91,4        |
| A5219104 | A5218104 | 4     | Left  | 106,4       |
| A5219105 | A5218105 | 5     | Left  | 121,4       |
| A5219106 | A5218106 | 6     | Left  | 136,4       |
| A5219107 | A5218107 | 7     | Left  | 151,4       |
| A5219108 | A5218108 | 8     | Left  | 166,4       |
| A5219203 | A5218203 | 3     | Right | 91,4        |
| A5219204 | A5218204 | 4     | Right | 106,4       |
| A5219205 | A5218205 | 5     | Right | 121,4       |
| A5219206 | A5218206 | 6     | Right | 136,4       |
| A5219207 | A5218207 | 7     | Right | 151,4       |
| A5219208 | A5218208 | 8     | Right | 166,4       |



**Intended use :** Temporary fixation of fractures, malunions, nonunions, osteotomies and pseudoarthroses of the distal humerus in adult patients.

## Small Fragment Non-Locking Plates


### 3,5 mm Y-Distal Humerus Plate 2



| SS               | TI               | HOLES | LENGTH (mm) |
|------------------|------------------|-------|-------------|
| A100 02 224 1205 | A110 02 224 1205 | 5     | 137         |
| A100 02 224 1206 | A110 02 224 1206 | 6     | 153         |
| A100 02 224 1207 | A110 02 224 1207 | 7     | 169         |
| A100 02 224 1209 | A110 02 224 1209 | 9     | 201         |

**Intended use :** Temporary fixation of fractures, malunions, nonunions, osteotomies and pseudoarthroses of the distal humerus in adult patients.

### 3,5 mm Ulnar Olecranon Plate



| SS               | TI        | HOLES | SIDE  | LENGTH (mm) |
|------------------|-----------|-------|-------|-------------|
| A100 02 224 1104 | A22014040 | 4     | Left  | 79,6        |
| A100 02 224 1105 | A22014050 | 5     | Left  | 92,1        |
| A100 02 224 1103 | A22014060 | 6     | Left  | 104,6       |
| A100 02 224 1004 | A22014030 | 4     | Right | 79,6        |
| A100 02 224 1005 | A22014070 | 5     | Right | 92,1        |
| A100 02 224 1003 | A22014080 | 6     | Right | 104,6       |


**Intended use :** Temporary fixation of fractures, malunions, nonunions, osteotomies and pseudoarthroses of the proximal ulna in adult patients.

# Small Fragment Non-Locking Plates

## 3,5 mm Distal Medial Humerus Reconstruction Plate



| TI       | HOLES | LENGTH (mm) |
|----------|-------|-------------|
| A5249010 | 10    | 69,6        |
| A5249012 | 12    | 85          |
| A5249014 | 14    | 100,4       |




**Intended use :** Temporary fixation of fractures, malunions, nonunions, osteotomies and pseudoarthroses of the distal humerus (medial approach) in adults.

## 3,5 mm Y-Distal Humerus Reconstruction Plate




| SS                | TI                | LENGTH (mm) |
|-------------------|-------------------|-------------|
| A 100 02 224 0000 | A 110 02 224 0000 | 88,6        |



**Intended use :** Temporary fixation of fractures, malunions, nonunions, osteotomies and pseudoarthroses of the distal humerus in adults.

## Small Fragment Non-Locking Plates


### 1/3 Tubular Plate



**Intended use :** Fractures, osteotomies, nonunions and malunions of small bones, including the ulna, radius and fibula.

| SS               | TI               | HOLES | LENGTH (mm) |
|------------------|------------------|-------|-------------|
| A100 02 181 0004 | A110 02 181 0004 | 4     | 54          |
| A100 02 181 0005 | A110 02 181 0005 | 5     | 66          |
| A100 02 181 0006 | A110 02 181 0006 | 6     | 78          |
| A100 02 181 0007 | A110 02 181 0007 | 7     | 90          |
| A100 02 181 0008 | A110 02 181 0008 | 8     | 102         |
| A100 02 181 0009 | A110 02 181 0009 | 9     | 114         |
| A100 02 181 0010 | A110 02 181 0010 | 10    | 126         |
| A100 02 181 0011 | A110 02 181 0011 | 11    | 138         |
| A100 02 181 0012 | A110 02 181 0012 | 12    | 150         |

### Quarter 1/4 Tubular Plate



**Intended use :** Temporary fixation of fractures, malunions, nonunions, osteotomies and pseudoarthroses of small bones, including the ulna, radius and fibula, metacarpal and metatarsal bones.

| SS               | TI               | HOLES | LENGTH (mm) |
|------------------|------------------|-------|-------------|
| A100 02 182 0003 | A110 02 170 0003 | 3     | 23          |
| A100 02 182 0004 | A110 02 170 0004 | 4     | 31          |
| A100 02 182 0005 | A110 02 170 0005 | 5     | 39          |
| A100 02 182 0006 | A110 02 170 0006 | 6     | 47          |
| A100 02 182 0007 | A110 02 170 0007 | 7     | 55          |
| A100 02 182 0008 | A110 02 170 0008 | 8     | 63          |
| A100 02 182 0009 | A110 02 170 0009 | 9     | 71          |
| A100 02 182 0010 | A110 02 170 0010 | 10    | 79          |

## Small Fragment Non-Locking Plates

### 3,5 mm Thin Narrow Plate DCP



| SS               | TI               | HOLES | LENGTH (mm) |
|------------------|------------------|-------|-------------|
| A100 02 171 0004 | A110 02 171 0004 | 4     | 54          |
| A100 02 171 0005 | A110 02 171 0005 | 5     | 66          |
| A100 02 171 0006 | A110 02 171 0006 | 6     | 78          |
| A100 02 171 0007 | A110 02 171 0007 | 7     | 90          |
| A100 02 171 0008 | A110 02 171 0008 | 8     | 102         |
| A100 02 171 0009 | A110 02 171 0009 | 9     | 114         |
| A100 02 171 0010 | A110 02 171 0010 | 10    | 126         |
| A100 02 171 0012 | A110 02 171 0012 | 12    | 150         |




**Intended use :** Temporary fixation of fractures, malunions, nonunions, osteotomies and pseudoarthroses of the clavicle, scapula, olecranon, humerus, radius, ulna, pelvis, distal tibia, fibula. These plates are also indicated for fracture fixation of diaphyseal and metaphyseal areas of long bones in pediatric patients.

### 3,5 mm Straight Reconstruction Plate




| SS               | TI               | HOLES | LENGTH (mm) |
|------------------|------------------|-------|-------------|
| A100 02 220 0004 | A110 02 220 0004 | 4     | 48          |
| A100 02 220 0005 | A110 02 220 0005 | 5     | 60          |
| A100 02 220 0006 | A110 02 220 0006 | 6     | 72          |
| A100 02 220 0007 | A110 02 220 0007 | 7     | 84          |
| A100 02 220 0008 | A110 02 220 0008 | 8     | 96          |
| A100 02 220 0009 | A110 02 220 0009 | 9     | 108         |
| A100 02 220 0010 | A110 02 220 0010 | 10    | 120         |
| A100 02 220 0011 | A110 02 220 0011 | 11    | 132         |
| A100 02 220 0012 | A110 02 220 0012 | 12    | 144         |



**Intended use :** Temporary fixation of fractures, malunions, nonunions, osteotomies and pseudoarthroses of the clavicle, scapula, olecranon, humerus, radius, ulna, pelvis, distal tibia, fibula for adult patients. These plates are also indicated for fracture fixation of diaphyseal and metaphyseal areas of long bones in pediatric patients.

## Small Fragment Non-Locking Plates


### 3,5 mm Curved Reconstruction Plate



**Intended use :** Temporary fixation of fractures, malunions, nonunions, osteotomies and pseudoarthroses of the pelvis

| SS               | Ti               | HOLES | LENGTH (mm) |
|------------------|------------------|-------|-------------|
| A100 02 221 0002 | A110 02 221 0002 | 2     | 29,1        |
| A100 02 221 0003 | A110 02 221 0003 | 3     | 45,1        |
| A100 02 221 0004 | A110 02 221 0004 | 4     | 60,8        |
| A100 02 221 0005 | A110 02 221 0005 | 5     | 76,9        |
| A100 02 221 0006 | A110 02 221 0006 | 6     | 92,7        |
| A100 02 221 0007 | A110 02 221 0007 | 7     | 108,2       |
| A100 02 221 0008 | A110 02 221 0008 | 8     | 123,6       |
| A100 02 221 0009 | A110 02 221 0009 | 9     | 138,7       |
| A100 02 221 0010 | A110 02 221 0010 | 10    | 153,9       |
| A100 02 221 0011 | A110 02 221 0011 | 11    | 168,4       |
| A100 02 221 0012 | A110 02 221 0012 | 12    | 183,2       |

### 3,5 mm T-Plate Right-Angled (3 Holes On Head)



**Intended use :** Temporary fixation of complex intra-articular and extra-articular fractures and osteotomies of the distal radius in adults, skeletally mature adolescents.


| SS              | Ti              | HOLES | LENGTH (mm) |
|-----------------|-----------------|-------|-------------|
| A300 02 079 003 | A310 02 079 003 | 3     | 49          |
| A300 02 079 004 | A310 02 079 004 | 4     | 61          |
| A300 02 079 005 | A310 02 079 005 | 5     | 73          |
| A300 02 079 006 | A310 02 079 006 | 6     | 85          |

## Small Fragment Non-Locking Plates

### 3,5 mm T-Plate Right-Angled (4 Holes On Head)



| SS               | TI               | HOLES | LENGTH (mm) |
|------------------|------------------|-------|-------------|
| A300 02 080 0003 | A310 02 080 0003 | 3     | 46          |
| A300 02 080 0004 | A310 02 080 0004 | 4     | 57          |
| A300 02 080 0005 | A310 02 080 0005 | 5     | 68          |
| A300 02 080 0006 | A310 02 080 0006 | 6     | 79          |




**Intended use :** Temporary fixation of complex intra-articular and extra-articular fractures and osteotomies of the distal radius in adults, skeletally mature adolescents.

### 3,5 mm Compression T-Plate Oblique-Angled



| SS               | TI               | HOLES | SIDE  | LENGTH (mm) |
|------------------|------------------|-------|-------|-------------|
| A300 02 090 0003 | A310 02 090 0003 | 3     | Left  | 52,5        |
| A300 02 090 0004 | A310 02 090 0004 | 4     | Left  | 66          |
| A300 02 090 0005 | A310 02 090 0005 | 5     | Left  | 79,5        |
| A300 02 090 0006 | A310 02 090 0006 | 6     | Left  | 93          |
| A300 02 091 0003 | A310 02 091 0003 | 3     | Right | 52,5        |
| A300 02 091 0004 | A310 02 091 0004 | 4     | Right | 66          |
| A300 02 091 0005 | A310 02 091 0005 | 5     | Right | 79,5        |
| A300 02 091 0006 | A310 02 091 0006 | 6     | Right | 93          |



**Intended use :** Temporary fixation of complex intra-articular and extra-articular fractures and osteotomies of the distal radius in adults, skeletally mature adolescents.

"We supply WORLD's  
Orthopedic & Trauma  
Implants..."



# Small Fragment Screws

Locking/Non-Locking



## Small Fragment Locking Screws

### 3,5 mm Self-Tapping Locking Screw



| SS       | TI       | DIA. x LENGTH (mm) | SS       | TI       | DIA. x LENGTH (mm) |
|----------|----------|--------------------|----------|----------|--------------------|
| A1993510 | A1983510 | 3,5X10             | A1993530 | A1983530 | 3,5X30             |
| A1993512 | A1983512 | 3,5X12             | A1993532 | A1983532 | 3,5X32             |
| A1993514 | A1983514 | 3,5X14             | A1993534 | A1983534 | 3,5X34             |
| A1993516 | A1983516 | 3,5X16             | A1993536 | A1983536 | 3,5X36             |
| A1993518 | A1983518 | 3,5X18             | A1993538 | A1983538 | 3,5X38             |
| A1993520 | A1983520 | 3,5X20             | A1993540 | A1983540 | 3,5X40             |
| A1993522 | A1983522 | 3,5X22             | A1993542 | A1983542 | 3,5X42             |
| A1993524 | A1983524 | 3,5X24             | A1993544 | A1983544 | 3,5X44             |
| A1993526 | A1983526 | 3,5X26             | A1993546 | A1983546 | 3,5X46             |
| A1993528 | A1983528 | 3,5X28             | A1993548 | A1983548 | 3,5X48             |
|          |          |                    | A1993550 | A1983550 | 3,5X50             |
|          |          |                    | A1993555 | A1983555 | 3,5X55             |
|          |          |                    | A1993560 | A1983560 | 3,5X60             |



### 3,5 mm Low-Profile Locking Screw



| TI       | DIA. x LENGTH (mm) | TI       | DIA. x LENGTH (mm) |
|----------|--------------------|----------|--------------------|
| A1973510 | 3,5X10             | A1973520 | 3,5X20             |
| A1973512 | 3,5X12             | A1973522 | 3,5X22             |
| A1973514 | 3,5X14             | A1973524 | 3,5X24             |
| A1973516 | 3,5X16             | A1973526 | 3,5X26             |
| A1973518 | 3,5X18             | A1973528 | 3,5X28             |



### 3,5 mm Locking Cancellous Screw Full Threaded



| SS        | TI        | DIA. x LENGTH (mm) | SS        | TI        | DIA. x LENGTH (mm) |
|-----------|-----------|--------------------|-----------|-----------|--------------------|
| A33803510 | A33703510 | 3,5X10             | A33803530 | A33703530 | 3,5X30             |
| A33803512 | A33703512 | 3,5X12             | A33803532 | A33703532 | 3,5X32             |
| A33803514 | A33703514 | 3,5X14             | A33803534 | A33703534 | 3,5X34             |
| A33803516 | A33703516 | 3,5X16             | A33803536 | A33703536 | 3,5X36             |
| A33803518 | A33703518 | 3,5X18             | A33803538 | A33703538 | 3,5X38             |
| A33803520 | A33703520 | 3,5X20             | A33803540 | A33703540 | 3,5X40             |
| A33803522 | A33703522 | 3,5X22             | A33803542 | A33703542 | 3,5X42             |
| A33803524 | A33703524 | 3,5X24             | A33803544 | A33703544 | 3,5X44             |
| A33803526 | A33703526 | 3,5X26             | A33803546 | A33703546 | 3,5X46             |
| A33803528 | A33703528 | 3,5X28             | A33803548 | A33703548 | 3,5X48             |
|           |           |                    | A33803550 | A33703550 | 3,5X50             |
|           |           |                    | A33803555 | A33703555 | 3,5X55             |
|           |           |                    | A33803560 | A33703560 | 3,5X60             |



## Small Fragment Locking Screws

### 3,5 mm Locking Cancellous Screw Partially Threaded



| SS        | TI        | DIA. x LENGTH (mm) | SS        | TI        | DIA. x LENGTH (mm) |
|-----------|-----------|--------------------|-----------|-----------|--------------------|
| A33813520 | A33713520 | 3,5X20             | A33813538 | A33713538 | 3,5X38             |
| A33813522 | A33713522 | 3,5X22             | A33813540 | A33713540 | 3,5X40             |
| A33813524 | A33713524 | 3,5X24             | A33813542 | A33713542 | 3,5X42             |
| A33813526 | A33713526 | 3,5X26             | A33813544 | A33713544 | 3,5X44             |
| A33813528 | A33713528 | 3,5X28             | A33813546 | A33713546 | 3,5X46             |
| A33813530 | A33713530 | 3,5X30             | A33813548 | A33713548 | 3,5X48             |
| A33813532 | A33713532 | 3,5X32             | A33813550 | A33713550 | 3,5X50             |
| A33813534 | A33713534 | 3,5X34             | A33813555 | A33713555 | 3,5X55             |
| A33813536 | A33713536 | 3,5X36             | A33813560 | A33713560 | 3,5X60             |




### 3,5 mm Cable Locking Screw



TI  
A900 03 020 0005

### 3,5 mm Locking Screw Plug



TI  
A900 03 020 0009

### 2,4 mm Locking Screw (5,5 mm Head Diameter)



| SS       | TI       | DIA. x LENGTH (mm) | SS       | TI       | DIA. x LENGTH (mm) |
|----------|----------|--------------------|----------|----------|--------------------|
| A1992306 | A2065014 | 2,4X6              | A1992324 | A2065032 | 2,4X24             |
| A1992308 | A2065016 | 2,4X8              | A1992326 | A2065034 | 2,4X26             |
| A1992310 | A2065018 | 2,4X10             | A1992328 | A2065036 | 2,4X28             |
| A1992312 | A2065020 | 2,4X12             | A1992330 | A2065038 | 2,4X30             |
| A1992314 | A2065022 | 2,4X14             | A1992332 | A2065040 | 2,4X32             |
| A1992316 | A2065024 | 2,4X16             | A1992334 | A2065042 | 2,4X34             |
| A1992318 | A2065026 | 2,4X18             | A1992336 | A2065044 | 2,4X36             |
| A1992320 | A2065028 | 2,4X20             | A1992338 | A2065046 | 2,4X38             |
| A1992322 | A2065030 | 2,4X22             | A1992340 | A2065048 | 2,4X40             |



# Small Fragment Locking Screws

## 2,7 mm Mini Locking Screw



| TI       | DIA. x LENGTH (mm) | TI       | DIA. x LENGTH (mm) |
|----------|--------------------|----------|--------------------|
| A1982710 | 2,7X10             | A1982726 | 2,7X26             |
| A1982712 | 2,7X12             | A1982728 | 2,7X28             |
| A1982714 | 2,7X14             | A1982730 | 2,7X30             |
| A1982716 | 2,7X16             | A1982732 | 2,7X32             |
| A1982718 | 2,7X18             | A1982734 | 2,7X34             |
| A1982720 | 2,7X20             | A1982736 | 2,7X36             |
| A1982722 | 2,7X22             | A1982738 | 2,7X38             |
| A1982724 | 2,7X24             | A1982740 | 2,7X40             |

## 2,7 mm Mini Locking Screw (5,5 mm Head Diameter)



| SS       | TI       | DIA. x LENGTH (mm) | SS       | TI       | DIA. x LENGTH (mm) |
|----------|----------|--------------------|----------|----------|--------------------|
| A1992706 | A2065050 | 2,7X6              | A1992720 | A2065085 | 2,7X20             |
| A1992708 | A2065055 | 2,7X8              | A1992722 | A2065090 | 2,7X22             |
| A1992710 | A2065060 | 2,7X10             | A1992724 | A1992005 | 2,7X24             |
| A1992712 | A2065065 | 2,7X12             | A1992726 | A1992007 | 2,7X26             |
| A1992714 | A2065070 | 2,7X14             | A1992728 | A1992009 | 2,7X28             |
| A1992716 | A2065075 | 2,7X16             | A1992730 | A1992011 | 2,7X30             |
| A1992718 | A2065080 | 2,7X18             | A1992732 | A1992013 | 2,7X32             |
|          |          |                    | A1992734 | A1992015 | 2,7X34             |
|          |          |                    | A1992736 | A1992017 | 2,7X36             |
|          |          |                    | A1992738 | A1992019 | 2,7X38             |
|          |          |                    | A1992740 | A1992021 | 2,7X40             |

## 3,5 mm Self-Drilling Locking Screw



| SS       | TI       | DIA. x LENGTH (mm) | SS       | TI       | DIA. x LENGTH (mm) |
|----------|----------|--------------------|----------|----------|--------------------|
| A2393510 | A2383510 | 3,5X10             | A2393530 | A2383530 | 3,5X30             |
| A2393512 | A2383512 | 3,5X12             | A2393532 | A2383532 | 3,5X32             |
| A2393514 | A2383514 | 3,5X14             | A2393534 | A2383534 | 3,5X34             |
| A2393516 | A2383516 | 3,5X16             | A2393536 | A2383536 | 3,5X36             |
| A2393518 | A2383518 | 3,5X18             | A2393538 | A2383538 | 3,5X38             |
| A2393520 | A2383520 | 3,5X20             | A2393540 | A2383540 | 3,5X40             |
| A2393522 | A2383522 | 3,5X22             | A2393542 | A2383542 | 3,5X42             |
| A2393524 | A2383524 | 3,5X24             | A2393544 | A2383544 | 3,5X44             |
| A2393526 | A2383526 | 3,5X26             | A2393546 | A2383546 | 3,5X46             |
| A2393528 | A2383528 | 3,5X28             | A2393548 | A2383548 | 3,5X48             |
|          |          |                    | A2393550 | A2383550 | 3,5X50             |
|          |          |                    | A2393555 | A2383555 | 3,5X55             |
|          |          |                    | A2393560 | A2383560 | 3,5X60             |

## Small Fragment Locking/Non-Locking Screws

### 3,5 mm Self-Tapping Locking Cannulated Screw



| SS       | TI       | DIA. x LENGTH (mm) | SS       | TI       | DIA. x LENGTH (mm) |
|----------|----------|--------------------|----------|----------|--------------------|
| A1953510 | A1963510 | 3,5X10             | A1863532 | A1963532 | 3,5X32             |
| A1953512 | A1963512 | 3,5X12             | A1863534 | A1963534 | 3,5X34             |
| A1953514 | A1963514 | 3,5X14             | A1863536 | A1963536 | 3,5X36             |
| A1953516 | A1963516 | 3,5X16             | A1863538 | A1963538 | 3,5X38             |
| A1953518 | A1963518 | 3,5X18             | A1953540 | A1963540 | 3,5X40             |
| A1953520 | A1963520 | 3,5X20             | A1863518 | A1963565 | 3,5X42             |
| A1953522 | A1963522 | 3,5X22             | A1863520 | A1963570 | 3,5X44             |
| A1953524 | A1963524 | 3,5X24             | A1863522 | A1963546 | 3,5X46             |
| A1953526 | A1963526 | 3,5X26             | A1863524 | A1963548 | 3,5X48             |
| A1953528 | A1963528 | 3,5X28             | A1953550 | A1963550 | 3,5X50             |
| A1953530 | A1963530 | 3,5X30             | A1953555 | A1963555 | 3,5X55             |
|          |          |                    | A1953560 | A1963560 | 3,5X60             |

### 2,4 mm Mini Locking Screw



| TI               | DIA. x LENGTH (mm) | TI               | DIA. x LENGTH (mm) |
|------------------|--------------------|------------------|--------------------|
| A110 03 213 0060 | 2,4X6              | A100 03 216 0015 | 2,4X15             |
| A100 03 216 0070 | 2,4X7              | A110 03 213 0160 | 2,4X16             |
| A110 03 213 0080 | 2,4X8              | A110 03 213 0180 | 2,4X18             |
| A100 03 216 0090 | 2,4X9              | A110 03 213 0200 | 2,4X20             |
| A110 03 213 0100 | 2,4X10             | A110 03 213 0220 | 2,4X22             |
| A100 03216 0011  | 2,4X11             | A110 03 213 0240 | 2,4X24             |
| A110 03 213 0120 | 2,4X12             | A110 03 213 0260 | 2,4X26             |
| A100 03 216 0013 | 2,4X13             | A110 03 213 0280 | 2,4X28             |
| A110 03 213 0140 | 2,4X14             | A110 03 213 0300 | 2,4X30             |

### 2,4 mm Mini Screw



| TI              | DIA. x LENGTH (mm) | TI              | DIA. x LENGTH (mm) |
|-----------------|--------------------|-----------------|--------------------|
| A1982306        | 2,4X6              | A100 03216 0025 | 2,4X15             |
| A100 03216 0017 | 2,4X7              | A1982316        | 2,4X16             |
| A1982308        | 2,4X8              | A1982318        | 2,4X18             |
| A100 03216 0019 | 2,4X9              | A1982320        | 2,4X20             |
| A1982310        | 2,4X10             | A1982322        | 2,4X22             |
| A100 03216 0021 | 2,4X11             | A1982324        | 2,4X24             |
| A1982312        | 2,4X12             | A1982326        | 2,4X26             |
| A100 03216 0023 | 2,4X13             | A1982328        | 2,4X28             |
| A1982314        | 2,4X14             | A1982330        | 2,4X30             |

## Small Fragment Non-Locking Screws

### 3,5 mm Cortical Screw (Pitch 1,25 mm)



| SS               | TI               | DIA. x LENGTH (mm) |
|------------------|------------------|--------------------|
| A100 03 218 0100 | A110 03 218 0100 | 3,5X10             |
| A100 03 218 0120 | A110 03 218 0120 | 3,5X12             |
| A100 03 218 0140 | A110 03 218 0140 | 3,5X14             |
| A100 03 218 0160 | A110 03 218 0160 | 3,5X16             |
| A100 03 218 0180 | A110 03 218 0180 | 3,5X18             |
| A100 03 218 0200 | A110 03 218 0200 | 3,5X20             |
| A100 03 218 0220 | A110 03 218 0220 | 3,5X22             |
| A100 03 218 0240 | A110 03 218 0240 | 3,5X24             |
| A100 03 218 0260 | A110 03 218 0260 | 3,5X26             |
| A100 03 218 0280 | A110 03 218 0280 | 3,5X28             |

| SS               | TI               | DIA. x LENGTH (mm) |
|------------------|------------------|--------------------|
| A100 03 218 0300 | A110 03 218 0300 | 3,5X30             |
| A100 03 218 0320 | A110 03 218 0320 | 3,5X32             |
| A100 03 218 0340 | A110 03 218 0340 | 3,5X34             |
| A100 03 218 0360 | A110 03 218 0360 | 3,5X36             |
| A100 03 218 0380 | A110 03 218 0380 | 3,5X38             |
| A100 03 218 0400 | A110 03 218 0400 | 3,5X40             |
| A100 03 218 0420 | A110 03 218 0420 | 3,5X42             |
| A100 03 218 0440 | A110 03 218 0440 | 3,5X44             |
| A100 03 218 0460 | A110 03 218 0460 | 3,5X46             |
| A100 03 218 0480 | A110 03 218 0480 | 3,5X48             |
| A100 03 218 0500 | A110 03 218 0500 | 3,5X50             |
| A100 03 218 0520 | A110 03 218 0520 | 3,5X52             |
| A100 03 218 0540 | A110 03 218 0540 | 3,5X54             |
| A100 03 218 0560 | A110 03 218 0560 | 3,5X56             |
| A100 03 218 0580 | A110 03 218 0580 | 3,5X58             |
| A100 03 218 0600 | A110 03 218 0600 | 3,5X60             |



### 3,5 mm Cancellous Bone Screw Full Threaded



| SS               | TI               | DIA. x LENGTH (mm) |
|------------------|------------------|--------------------|
| A100 03 219 0100 | A110 03 219 0100 | 3,5X10             |
| A100 03 219 0120 | A110 03 219 0120 | 3,5X12             |
| A100 03 219 0140 | A110 03 219 0140 | 3,5X14             |
| A100 03 219 0160 | A110 03 219 0160 | 3,5X16             |
| A100 03 219 0180 | A110 03 219 0180 | 3,5X18             |
| A100 03 219 0200 | A110 03 219 0200 | 3,5X20             |
| A100 03 219 0220 | A110 03 219 0220 | 3,5X22             |
| A100 03 219 0240 | A110 03 219 0240 | 3,5X24             |
| A100 03 219 0260 | A110 03 219 0260 | 3,5X26             |
| A100 03 219 0280 | A110 03 219 0280 | 3,5X28             |
| A100 03 219 0300 | A110 03 219 0300 | 3,5X30             |

| SS               | TI               | DIA. x LENGTH (mm) |
|------------------|------------------|--------------------|
| A100 03 219 0320 | A110 03 219 0320 | 3,5X32             |
| A100 03 219 0340 | A110 03 219 0340 | 3,5X34             |
| A100 03 219 0360 | A110 03 219 0360 | 3,5X36             |
| A100 03 219 0380 | A110 03 219 0380 | 3,5X38             |
| A100 03 219 0400 | A110 03 219 0400 | 3,5X40             |
| A100 03 219 0420 | A6183542         | 3,5X42             |
| A100 03 219 0440 | A110 03 219 0450 | 3,5X44             |
| A100 03 219 0460 | A6183546         | 3,5X46             |
| A100 03 219 0480 | A6183548         | 3,5X48             |
| A100 03 219 0500 | A110 03 219 0500 | 3,5X50             |
| A100 03 219 0520 | A6183552         | 3,5X52             |
| A100 03 219 0540 | A110 03 219 0540 | 3,5X54             |
| A100 03 219 0560 | A6183556         | 3,5X56             |
| A100 03 219 0580 | A110 03 219 0580 | 3,5X58             |
| A100 03 219 0600 | A6183560         | 3,5X60             |



### 3,5 mm Cancellous Bone Screw Partially Threaded



| SS       | TI       | DIA. x LENGTH (mm) |
|----------|----------|--------------------|
| A6353520 | A6153520 | 3,5X20             |
| A6353522 | A6153522 | 3,5X22             |
| A6353524 | A6153524 | 3,5X24             |
| A6353526 | A6153526 | 3,5X26             |
| A6353528 | A6153528 | 3,5X28             |
| A6353530 | A6153530 | 3,5X30             |
| A6353532 | A6153532 | 3,5X32             |

| SS       | TI       | DIA. x LENGTH (mm) |
|----------|----------|--------------------|
| A6353534 | A6153534 | 3,5X34             |
| A6353536 | A6153536 | 3,5X36             |
| A6353538 | A6153538 | 3,5X38             |
| A6353540 | A6153540 | 3,5X40             |
| A6353542 | A6153542 | 3,5X42             |
| A6353544 | A6153544 | 3,5X44             |
| A6353546 | A6153546 | 3,5X46             |
| A6353548 | A6153548 | 3,5X48             |
| A6353550 | A6153550 | 3,5X50             |



## Small Fragment Non-Locking Screws

### 4,0 mm Cancellous Bone Screw Full Threaded



| SS               | TI               | DIA. x LENGTH (mm) |
|------------------|------------------|--------------------|
| A100 03 220 0100 | A110 03 220 0100 | 4,0X10             |
| A100 03 220 0120 | A110 03 220 0120 | 4,0X12             |
| A100 03 220 0140 | A110 03 220 0140 | 4,0X14             |
| A100 03 220 0160 | A110 03 220 0160 | 4,0X16             |
| A100 03 220 0180 | A110 03 220 0180 | 4,0X18             |



| SS               | TI               | DIA. x LENGTH (mm) |
|------------------|------------------|--------------------|
| A100 03 220 0200 | A110 03 220 0200 | 4,0X20             |
| A100 03 220 0220 | A110 03 220 0220 | 4,0X22             |
| A100 03 220 0240 | A110 03 220 0240 | 4,0X24             |
| A100 03 220 0260 | A110 03 220 0260 | 4,0X26             |
| A100 03 220 0280 | A110 03 220 0280 | 4,0X28             |
| A100 03 220 0300 | A110 03 220 0300 | 4,0X30             |
| A100 03 220 0350 | A110 03 220 0350 | 4,0X35             |
| A100 03 220 0400 | A110 03 220 0400 | 4,0X40             |
| A100 03 220 0450 | A110 03 220 0450 | 4,0X45             |
| A100 03 220 0500 | A110 03 220 0500 | 4,0X50             |
| A100 03 220 0550 | A110 03 220 0550 | 4,0X55             |
| A100 03 220 0600 | A110 03 220 0600 | 4,0X60             |

### 4,0 mm Cancellous Bone Screw Partially Threaded



| SS               | TI               | DIA. x LENGTH (mm) |
|------------------|------------------|--------------------|
| A100 03 021 0200 | A110 03 021 0200 | 4,0X20             |
| A100 03 021 0220 | A110 03 021 0220 | 4,0X22             |
| A100 03 021 0240 | A110 03 021 0240 | 4,0X24             |
| A100 03 021 0260 | A110 03 021 0260 | 4,0X26             |
| A100 03 021 0280 | A110 03 021 0280 | 4,0X28             |
| A100 03 021 0300 | A110 03 021 0300 | 4,0X30             |
| A100 03 021 0350 | A110 03 021 0350 | 4,0X35             |
| A100 03 021 0400 | A110 03 021 0400 | 4,0X40             |
| A100 03 021 0450 | A110 03 021 0450 | 4,0X45             |
| A100 03 021 0500 | A110 03 021 0500 | 4,0X50             |
| A100 03 021 0550 | A110 03 021 0550 | 4,0X55             |
| A100 03 021 0600 | A110 03 021 0600 | 4,0X60             |



### 2,7 mm Cortical Screw



| SS               | TI               | DIA. x LENGTH (mm) |
|------------------|------------------|--------------------|
| A100 03 217 0060 | A110 03 222 0110 | 2,7X6              |
| A100 03 217 0080 | A110 03 222 0112 | 2,7X8              |
| A100 03 217 0100 | A110 03 222 0140 | 2,7X10             |
| A100 03 217 0120 | A110 03 222 0160 | 2,7X12             |
| A100 03 217 0140 | A110 03 222 0180 | 2,7X14             |
| A100 03 217 0160 | A110 03 222 0200 | 2,7X16             |

| SS               | TI               | DIA. x LENGTH (mm) |
|------------------|------------------|--------------------|
| A100 03 217 0180 | A110 03 222 0220 | 2,7X18             |
| A100 03 217 0200 | A110 03 222 0240 | 2,7X20             |
| A100 03 217 0220 | A110 03 222 0260 | 2,7X22             |
| A100 03 217 0240 | A110 03 222 0280 | 2,7X24             |
| A100 03 217 0260 | A110 03 222 0300 | 2,7X26             |
| A100 03 217 0280 | A110 03 222 0320 | 2,7X28             |
| A100 03 217 0300 | A110 03 222 0340 | 2,7X30             |
| A100 03 217 0320 | A110 03 222 0360 | 2,7X32             |
| A100 03 217 0340 | A110 03 222 0380 | 2,7X34             |
| A100 03 217 0360 | A110 03 222 0400 | 2,7X36             |
| A100 03 217 0380 | A110 03 222 0420 | 2,7X38             |
| A100 03 217 0400 | A110 03 222 0440 | 2,7X40             |



Small Fragment Non-Locking Screws

2,4 mm Cortical Screw



| SS               | TI               | DIA. x LENGTH<br>(mm) | SS               | TI               | DIA. x LENGTH<br>(mm) |
|------------------|------------------|-----------------------|------------------|------------------|-----------------------|
| A100 03 213 0060 | A110 03 222 0460 | 2,4X6                 | A100 03 213 0140 | A110 03 222 0540 | 2,4X14                |
| A100 03 213 0080 | A110 03 222 0480 | 2,4X8                 | A100 03 213 0160 | A110 03 222 0560 | 2,4X16                |
| A100 03 213 0100 | A110 03 222 0500 | 2,4X10                | A100 03 213 0180 | A110 03 222 0580 | 2,4X18                |
| A100 03 213 0120 | A110 03 222 0520 | 2,4X12                | A100 03 213 0200 | A110 03 222 0600 | 2,4X20                |
|                  |                  |                       | A100 03 213 0220 | A110 03 222 0620 | 2,4X22                |
|                  |                  |                       | A100 03 213 0240 | A110 03 222 0640 | 2,4X24                |
|                  |                  |                       | A100 03 213 0260 | A110 03 222 0650 | 2,4X26                |
|                  |                  |                       | A100 03 213 0280 | A110 03 222 0700 | 2,4X28                |
|                  |                  |                       | A100 03 213 0300 | A110 03 223 0100 | 2,4X30                |



2,0 mm Cortical Screw



| SS               | DIA. x LENGTH<br>(mm) | SS               | DIA. x LENGTH<br>(mm) |
|------------------|-----------------------|------------------|-----------------------|
| A100 03 216 0060 | 2,0X6                 | A100 03 216 0220 | 2,0X22                |
| A100 03 216 0080 | 2,0X8                 | A100 03 216 0240 | 2,0X24                |
| A100 03 216 0100 | 2,0X10                | A100 03 216 0260 | 2,0X26                |
| A100 03 216 0120 | 2,0X12                | A100 03 216 0280 | 2,0X28                |
| A100 03 216 0140 | 2,0X14                | A100 03 216 0300 | 2,0X30                |
| A100 03 216 0160 | 2,0X16                | A100 03 216 0340 | 2,0X34                |
| A100 03 216 0180 | 2,0X18                | A100 03 216 0360 | 2,0X36                |
| A100 03 216 0200 | 2,0X20                | A100 03 216 0380 | 2,0X38                |



1,5 mm Cortical Screw



| SS               | DIA. x LENGTH<br>(mm) | SS               | DIA. x LENGTH<br>(mm) |
|------------------|-----------------------|------------------|-----------------------|
| A100 03 215 0060 | 1,5X6                 | A100 03 215 0120 | 1,5X12                |
| A100 03 215 0070 | 1,5X7                 | A100 03 215 0140 | 1,5X14                |
| A100 03 215 0080 | 1,5X8                 | A100 03 215 0160 | 1,5X16                |
| A100 03 215 0090 | 1,5X9                 | A100 03 215 0180 | 1,5X18                |
| A100 03 215 0100 | 1,5X10                | A100 03 215 0200 | 1,5X20                |



"We supply WORLD's  
Orthopedic & Trauma  
Implants..."


# Large Fragment Plates

Locking/Non-Locking



# Large Fragment Locking Plates


## 4,5/5,0 mm Proximal Femur Locking Compression Plate 1



**Intended use :** Temporary fixation of fractures of the trochanteric region, fractures of the proximal femur, metastatic fractures of the proximal femur, osteotomies of the proximal femur and fixation of nonunions or malunions of proximal femur in adult patients.

| SS       | TI       | HOLES | SIDE  | LENGTH (mm) |
|----------|----------|-------|-------|-------------|
| A1120903 | A1120803 | 3     | Left  | 91,0        |
| A1120905 | A1120805 | 5     | Left  | 126,8       |
| A1120907 | A1120807 | 7     | Left  | 162,8       |
| A1120909 | A1120809 | 9     | Left  | 198,8       |
| A1120911 | A1120811 | 11    | Left  | 234,8       |
| A1120913 | A1120813 | 13    | Left  | 270,8       |
| A1130903 | A1130803 | 3     | Right | 91,0        |
| A1130905 | A1130805 | 5     | Right | 126,8       |
| A1130907 | A1130807 | 7     | Right | 162,8       |
| A1130909 | A1130809 | 9     | Right | 198,8       |
| A1130911 | A1130811 | 11    | Right | 234,8       |
| A1130913 | A1130813 | 13    | Right | 270,8       |

## 4,5/5,0/7,3 mm Proximal Femur Locking Compression Plate 2



**Intended use :** Temporary fixation of fractures of the trochanteric region, fractures of the proximal femur, metastatic fractures of the proximal femur, osteotomies of the proximal femur and fixation of nonunions or malunions of proximal femur in adult patients.

| SS       | TI       | HOLES | SIDE  | LENGTH (mm) |
|----------|----------|-------|-------|-------------|
| A1120902 | A1120802 | 2     | Left  | 119         |
| A1120904 | A1120804 | 4     | Left  | 159         |
| A1120906 | A1120806 | 6     | Left  | 195         |
| A1120908 | A1120808 | 8     | Left  | 231         |
| A1120910 | A1120810 | 10    | Left  | 267         |
| A1130902 | A1130802 | 2     | Right | 119         |
| A1130904 | A1130804 | 4     | Right | 159         |
| A1130906 | A1130806 | 6     | Right | 195         |
| A1130908 | A1130808 | 8     | Right | 231         |
| A1130910 | A1130914 | 10    | Right | 267         |

# Large Fragment Locking Plates


## Proximal Femur Compression Plate 7,3 mm Cannulated Locking Screw



| SS       | TI       | DIA. x LENGTH (mm) |
|----------|----------|--------------------|
| A6696050 | A6666050 | 7,3X50             |
| A6696055 | A6666055 | 7,3X55             |
| A6696060 | A6666060 | 7,3X60             |
| A6696065 | A6666065 | 7,3X65             |
| A6696070 | A6666070 | 7,3X70             |
| A6696075 | A6666075 | 7,3X75             |
| A6696080 | A6666080 | 7,3X80             |
| A6696085 | A6666085 | 7,3X85             |
| A6696090 | A6666090 | 7,3X90             |
| A6696095 | A6666095 | 7,3X95             |
| A6696100 | A6666100 | 7,3X100            |
| A6696105 | A6666105 | 7,3X105            |
| A6696110 | A6666110 | 7,3X110            |
| A6696115 | A6666115 | 7,3X115            |
| A6696120 | A6666120 | 7,3X120            |

**Intended use :** Temporary fixation of proximal femur fractures. Used only together with Proximal Femur Locking Compression Plate 2

## 4,5/5,0 mm Locking Plate Broad Curved



| SS       | TI       | HOLES | LENGTH (mm) |
|----------|----------|-------|-------------|
| A1871909 | A1871809 | 9     | 195,4       |
| A1871910 | A1871810 | 10    | 212,7       |
| A1871911 | A1872811 | 11    | 230,5       |
| A1871912 | A1872812 | 12    | 248,4       |
| A1871913 | A1872813 | 13    | 265,7       |
| A1871914 | A1872814 | 14    | 282,9       |
| A1871915 | A1872815 | 15    | 300,7       |
| A1871916 | A1872816 | 16    | 318,5       |


**Intended use :** Temporary fixation of fractures of femur shaft, fixation of periprosthetic fractures, fixation of nonunions or malunions of femur shaft in adult patients.

# Large Fragment Locking Plates

## 4,5/5,0 mm Distal Lateral Femur Locking Plate 1



| SS       | TI       | HOLES | SIDE  | LENGTH (mm) |
|----------|----------|-------|-------|-------------|
| A1441905 | A1441805 | 5     | Left  | 156,4       |
| A1441907 | A1441807 | 7     | Left  | 196,4       |
| A1441909 | A1441809 | 9     | Left  | 236,4       |
| A1441911 | A1441811 | 11    | Left  | 276,4       |
| A1441913 | A1441813 | 13    | Left  | 316,4       |
| A1452905 | A1452805 | 5     | Right | 156,4       |
| A1452907 | A1452807 | 7     | Right | 196,4       |
| A1452909 | A1452809 | 9     | Right | 236,4       |
| A1452911 | A1452811 | 11    | Right | 276,4       |
| A1452913 | A1452813 | 13    | Right | 316,4       |




**Intended use :** Temporary fixation of fractures, osteotomies, nonunions or malunions of the distal femur (intercondylar, supracondylar, intra-extra articular) in adult patients. Those plates are used by lateral approach.

## 4,5/5,0 mm Distal Lateral Femur Locking Plate 2



| SS       | TI       | HOLES | SIDE  | LENGTH (mm) |
|----------|----------|-------|-------|-------------|
| A6729106 | A6728106 | 6     | Left  | 158         |
| A6729107 | A6728107 | 7     | Left  | 174         |
| A6729108 | A6728108 | 8     | Left  | 191         |
| A6729109 | A6728109 | 9     | Left  | 206         |
| A6729110 | A6728110 | 10    | Left  | 223         |
| A6729112 | A6728111 | 11    | Left  | 250         |
| A6729206 | A6728206 | 6     | Right | 158         |
| A6729207 | A6728207 | 7     | Right | 174         |
| A6729208 | A6728208 | 8     | Right | 191         |
| A6729209 | A6728209 | 9     | Right | 206         |
| A6729210 | A6728210 | 10    | Right | 223         |
| A6729211 | A6728211 | 11    | Right | 250         |




**Intended use :** Temporary fixation of fractures, osteotomies, nonunions or malunions of the distal femur (intercondylar, supracondylar, intra-extra articular) in adult patients. Those plates are used by lateral approach.

# Large Fragment Locking Plates

## 4,5/5,0 mm Distal Medial Femur Locking Plate



| SS       | TI       | HOLES | SIDE  | LENGTH (mm) |
|----------|----------|-------|-------|-------------|
| A6739104 | A6738104 | 4     | Left  | 115,9       |
| A6739106 | A6738106 | 6     | Left  | 147,9       |
| A6739108 | A6738108 | 8     | Left  | 179,9       |
| A6739110 | A6738110 | 10    | Left  | 211,9       |
| A6739204 | A6738204 | 4     | Right | 115,9       |
| A6739206 | A6738206 | 6     | Right | 147,9       |
| A6739208 | A6738208 | 8     | Right | 179,9       |
| A6739210 | A6738210 | 10    | Right | 211,9       |




**Intended use :** Temporary fixation of fractures, osteotomies, nonunions or malunions of the distal femur (intercondylar, supracondylar, intra-extra articular) in adult patients. Those plates are used by medial approach.

## 4,5/5,0 mm Condylar Femur Locking Plate




| SS       | TI       | HOLES | SIDE  | LENGTH (mm) |
|----------|----------|-------|-------|-------------|
| A6749106 | A6748106 | 6     | Left  | 149,1       |
| A6749108 | A6748108 | 8     | Left  | 181,1       |
| A6749110 | A6748110 | 10    | Left  | 213,1       |
| A6749112 | A6748112 | 12    | Left  | 245,1       |
| A6749206 | A6748206 | 6     | Right | 149,1       |
| A6749208 | A6748208 | 8     | Right | 181,1       |
| A6749210 | A6748210 | 10    | Right | 213,1       |
| A6749212 | A6748212 | 12    | Right | 245,1       |



**Intended use :** Temporary fixation of fractures, osteotomies, nonunions or malunions of the distal femur (intercondylar, supracondylar, intra-extra articular) in adult patients. Those plates are used by lateral approach.

# Large Fragment Locking Plates


## 4,5/5,0 mm Proximal Lateral Tibia Locking Plate 1



| SS       | TI       | HOLES | SIDE  | LENGTH (mm) |
|----------|----------|-------|-------|-------------|
| A1461905 | A1461805 | 5     | Left  | 143,3       |
| A1461907 | A1461807 | 7     | Left  | 183,3       |
| A1461909 | A1461809 | 9     | Left  | 223,3       |
| A1461911 | A1461811 | 11    | Left  | 263,3       |
| A1461913 | A1461813 | 13    | Left  | 303,3       |
| A1472905 | A1472805 | 5     | Right | 143,3       |
| A1472907 | A1472807 | 7     | Right | 183,3       |
| A1472909 | A1472809 | 9     | Right | 223,3       |
| A1472911 | A1472811 | 11    | Right | 263,3       |
| A1472913 | A1472813 | 13    | Right | 303,3       |

**Intended use :** Temporary fixation of fractures, osteotomies, nonunions or malunions of the proximal tibia in adult patients. Those plates are used by medial approach.

## 4,5/5,0 mm Proximal Lateral Tibia Locking Plate 2




| SS       | TI       | HOLES | SIDE  | LENGTH (mm) |
|----------|----------|-------|-------|-------------|
| A9461905 | A9461805 | 5     | Left  | 102         |
| A9461907 | A9461807 | 7     | Left  | 138         |
| A9461909 | A9461809 | 9     | Left  | 174         |
| A9461911 | A9461811 | 11    | Left  | 210         |
| A9461913 | A9461813 | 13    | Left  | 246         |
| A9472905 | A9472805 | 5     | Right | 102         |
| A9472907 | A9472807 | 7     | Right | 138         |
| A9472909 | A9472809 | 9     | Right | 174         |
| A9472911 | A9472811 | 11    | Right | 210         |
| A9472913 | A9472813 | 13    | Right | 246         |

**Intended use :** Temporary fixation of fractures, osteotomies, nonunions or malunions of the proximal tibia in adult patients. Those plates are used by lateral approach.

# Large Fragment Locking Plates


## 4,5/5,0 mm Proximal Medial Tibia Locking Plate



**Intended use :** Temporary fixation of fractures, osteotomies, nonunions or malunions of the proximal tibia in adult patients. Those plates are used by medial approach.

| SS       | TI       | HOLES | SIDE  | LENGTH (mm) |
|----------|----------|-------|-------|-------------|
| A5599104 | A5598104 | 4     | Left  | 90,6        |
| A5599106 | A5598106 | 6     | Left  | 122,6       |
| A5599108 | A5598108 | 8     | Left  | 154,6       |
| A5599110 | A5598110 | 10    | Left  | 186,6       |
| A5599112 | A5598112 | 12    | Left  | 218,6       |
| A5599204 | A5598204 | 4     | Right | 90,6        |
| A5599206 | A5598206 | 6     | Right | 122,6       |
| A5599208 | A5598208 | 8     | Right | 154,6       |
| A5599210 | A5598210 | 10    | Right | 186,6       |
| A5599212 | A5598212 | 12    | Right | 218,6       |

## 4,5/5,0 mm Locking Compression L-Plate



**Intended use :** Temporary fixation of fractures, osteotomies, nonunions or malunions of the proximal tibia in adult patients.


| SS       | TI       | HOLES | SIDE  | LENGTH (mm) |
|----------|----------|-------|-------|-------------|
| A1561903 | A1561803 | 3     | Left  | 65,1        |
| A1561904 | A1561804 | 4     | Left  | 81,1        |
| A1561905 | A1561805 | 5     | Left  | 97,1        |
| A1561906 | A1561806 | 6     | Left  | 113,1       |
| A1561907 | A1561807 | 7     | Left  | 129,1       |
| A1561908 | A1561808 | 8     | Left  | 145,1       |
| A1572903 | A1572803 | 3     | Right | 65,1        |
| A1572904 | A1572804 | 4     | Right | 81,1        |
| A1572905 | A1572805 | 5     | Right | 97,1        |
| A1572906 | A1572806 | 6     | Right | 113,1       |
| A1572907 | A1572807 | 7     | Right | 129,1       |
| A1572908 | A1572808 | 8     | Right | 145,1       |

# Large Fragment Locking Plates

## 4,5/5,0 mm Locking Compression T-Plate



| SS       | TI       | HOLES | LENGTH (mm) |
|----------|----------|-------|-------------|
| A1490904 | A1490804 | 4     | 81,7        |
| A1490905 | A1490805 | 5     | 97,7        |
| A1490906 | A1490806 | 6     | 113,7       |
| A1490907 | A1490807 | 7     | 129,7       |
| A1490908 | A1490808 | 8     | 145,7       |




**Intended use :** Temporary fixation of fractures, osteotomies, nonunions or malunions of the proximal tibia in adult patients.

## 3,5/4,5/5,0 mm Metaphysis Locking Compression Plate




| SS       | TI       | HOLES | LENGTH (mm) |
|----------|----------|-------|-------------|
| A1850903 | A1850803 | 3     | 118         |
| A1850904 | A1850804 | 4     | 136         |
| A1850905 | A1850805 | 5     | 154         |
| A1850906 | A1850806 | 6     | 172         |
| A1850907 | A1850807 | 7     | 190         |
| A1850908 | A1850808 | 8     | 208         |
| A1850909 | A1850809 | 9     | 226         |
| A1850911 | A1850811 | 11    | 262         |



**Intended use :** Temporary fixation of extra-articular fractures, osteotomies, nonunions or malunions of the metaphyseal area of distal tibia, distal & proximal humerus that extend into the shaft.

# Large Fragment Locking Plates

## 4,5/5,0 mm Narrow Locking Compression Plate




| SS       | TI       | HOLES | LENGTH (mm) |
|----------|----------|-------|-------------|
| A1821105 | A1820905 | 5     | 110,6       |
| A1821106 | A1820906 | 6     | 128,6       |
| A1821107 | A1820907 | 7     | 146,6       |
| A1821108 | A1820908 | 8     | 164,6       |
| A1821109 | A1820909 | 9     | 182,6       |
| A1821110 | A1820910 | 10    | 200,6       |
| A1821111 | A1820911 | 11    | 218,6       |
| A1821112 | A1820912 | 12    | 236,6       |
| A1830816 | A1820913 | 13    | 254,6       |
| A1821114 | A1820914 | 14    | 272,6       |

**Intended use :** Temporary fixation of fractures, osteotomies, nonunions and malunions, pseudoarthroses of the humerus, tibia, femur in adult patients and adolescent patients.

## Attachment Locking Plate Medium




| TI               | TYPE   | LENGTH (mm) |
|------------------|--------|-------------|
| A110 02 190 0007 | Medium | 34,0        |



**Intended use :** For augmented Fixation of fractures, used as an attachment of 4,5/5,0 mm Narrow Locking Compression Plates.

## 4,5/5,0 mm Broad Locking Compression Plate



| SS       | TI       | HOLES | LENGTH (mm) |
|----------|----------|-------|-------------|
| A1830905 | A1830805 | 5     | 109         |
| A1830906 | A1830806 | 6     | 127         |
| A1830907 | A1830807 | 7     | 145         |
| A1830908 | A1830808 | 8     | 163         |
| A1830909 | A1830809 | 9     | 181         |
| A1830910 | A1830810 | 10    | 199         |
| A1830911 | A1830811 | 11    | 217         |
| A1830912 | A1830812 | 12    | 235         |
| A1830913 | A1830813 | 13    | 253         |
| A1830914 | A1830814 | 14    | 271         |

**Intended use :** Temporary fixation of fractures, osteotomies, nonunions and malunions of the femur. Those plates are used by lateral approach.

# Large Fragment Locking Plates


## Attachment Locking Plate Large



| TI               | TYPE  | LENGTH (mm) |
|------------------|-------|-------------|
| A110 02 190 0008 | Large | 41,8        |

**Intended use :** For augmented Fixation of fractures, used as an attachment of 4,5/5,0 mm Broad Locking Compression Plates.


## 4,5/5,0 mm Distal Tibia Antrolateral Locking Plate 1



| SS       | TI       | HOLES | SIDE  | LENGTH (mm) |
|----------|----------|-------|-------|-------------|
| A1261905 | A1261805 | 5     | Left  | 108,3       |
| A1261907 | A1261807 | 7     | Left  | 140,3       |
| A1261909 | A1261809 | 9     | Left  | 172,3       |
| A1261911 | A1261811 | 11    | Left  | 204,3       |
| A1261913 | A1261813 | 13    | Left  | 236,3       |
| A1261915 | A1261815 | 15    | Left  | 268,3       |
| A1272905 | A1272805 | 5     | Right | 108,3       |
| A1272907 | A1272807 | 7     | Right | 140,3       |
| A1272909 | A1272809 | 9     | Right | 172,3       |
| A1272911 | A1272811 | 11    | Right | 204,3       |
| A1272913 | A1272813 | 13    | Right | 236,3       |
| A1272915 | A1272815 | 15    | Right | 268,3       |

**Intended use :** Temporary fixation of fractures, osteotomies, nonunions and malunions of the distal tibia. Those plates are used by anterolateral approach.

## 3,5/4,0 mm Distal Tibia Antrolateral Locking Plate 2




| SS       | TI       | HOLES | SIDE  | LENGTH (mm) |
|----------|----------|-------|-------|-------------|
| A1431905 | A1431805 | 5     | Left  | 89,8        |
| A1431907 | A1431807 | 7     | Left  | 119,8       |
| A1431909 | A1431809 | 9     | Left  | 149,8       |
| A1431911 | A1431811 | 11    | Left  | 179,8       |
| A1431913 | A1431813 | 13    | Left  | 209,8       |
| A1431915 | A1431815 | 15    | Left  | 239,8       |
| A1431917 | A1431817 | 17    | Left  | 269,8       |
| A1432905 | A1432805 | 5     | Right | 89,8        |
| A1432907 | A1432807 | 7     | Right | 119,8       |
| A1432909 | A1432809 | 9     | Right | 149,8       |
| A1432911 | A1432811 | 11    | Right | 179,8       |
| A1432913 | A1432813 | 13    | Right | 209,8       |
| A1432915 | A1432815 | 15    | Right | 239,8       |
| A1432917 | A1432817 | 17    | Right | 269,8       |

**Intended use :** Temporary fixation of fractures, osteotomies, nonunions and malunions of the distal tibia. Those plates are used by anterolateral approach.

# Large Fragment Locking Plates


## 3,5/4,0 mm Distal Medial Tibia Locking Compression Plate 1



| SS       | TI       | HOLES | SIDE  | LENGTH (mm) |
|----------|----------|-------|-------|-------------|
| A1381904 | A1381804 | 4     | Left  | 116         |
| A1381906 | A1381806 | 6     | Left  | 142         |
| A1381908 | A1381808 | 8     | Left  | 168         |
| A1381910 | A1381810 | 10    | Left  | 194         |
| A1381912 | A1381812 | 12    | Left  | 220         |
| A1381914 | A1381814 | 14    | Left  | 246         |
| A1392904 | A1392804 | 4     | Right | 116         |
| A1392906 | A1392806 | 6     | Right | 142         |
| A1392908 | A1392808 | 8     | Right | 168         |
| A1392910 | A1392810 | 10    | Right | 194         |
| A1392912 | A1392812 | 12    | Right | 220         |
| A1392914 | A1392814 | 14    | Right | 246         |

**Intended use :** Temporary fixation of fractures, osteotomies, nonunions and malunions of the distal tibia. Those plates are used by medial approach.

## 3,5/4,0 mm Distal Medial Tibia Locking Compression Plate 2




| SS       | TI       | HOLES | SIDE  | LENGTH (mm) |
|----------|----------|-------|-------|-------------|
| A1401906 | A1401806 | 6     | Left  | 137,9       |
| A1401908 | A1401808 | 8     | Left  | 161,9       |
| A1401910 | A1401810 | 10    | Left  | 185,9       |
| A1401912 | A1401812 | 12    | Left  | 209,9       |
| A1401914 | A1401814 | 14    | Left  | 233,9       |
| A1412906 | A1412806 | 6     | Right | 137,9       |
| A1412908 | A1412808 | 8     | Right | 161,9       |
| A1412910 | A1412810 | 10    | Right | 185,9       |
| A1412912 | A1412812 | 12    | Right | 209,9       |
| A1412914 | A1412814 | 14    | Right | 233,9       |

**Intended use :** Temporary fixation of fractures, osteotomies, nonunions and malunions of the distal tibia. Those plates are used by medial approach.

# Large Fragment Locking Plates

## 3,5/4,0 mm Distal Tibia Locking Pilon Plate




**Intended use :** Temporary fixation of fractures, osteotomies, nonunions and malunions of the distal articular fractures of the tibia (pilon fractures) with emphasis on complex articular reconstruction.

| SS       | TI       | HOLES | LENGTH (mm) |
|----------|----------|-------|-------------|
| A1411816 | A1401816 | 16    | 111,5       |
| A1411818 | A1401818 | 18    | 142,5       |
| A1411820 | A1401820 | 20    | 173,5       |
| A1411822 | A1401822 | 22    | 204,5       |
| A1411824 | A1401824 | 24    | 235,5       |
| A1411826 | A1401826 | 26    | 266,5       |

## 3,5/4,0 mm Distal Fibula Locking Plate 1




**Intended use :** Temporary fixation of fractures, osteotomies and nonunions of the metaphyseal and diaphyseal region of the distal fibula.



| SS       | TI       | HOLES | SIDE  | LENGTH (mm) |
|----------|----------|-------|-------|-------------|
| A5097103 | A5096103 | 3     | Left  | 71          |
| A5097104 | A5096104 | 4     | Left  | 81          |
| A5097105 | A5096105 | 5     | Left  | 91          |
| A5097106 | A5096106 | 6     | Left  | 101         |
| A5097107 | A5096107 | 7     | Left  | 111         |
| A5097108 | A5096108 | 8     | Left  | 121         |
| A5097109 | A5096109 | 9     | Left  | 131         |
| A5097110 | A5096110 | 10    | Left  | 141         |
| A5097111 | A5096111 | 11    | Left  | 151         |
| A5097112 | A5096112 | 12    | Left  | 161         |
| A5097203 | A5096203 | 3     | Right | 71          |
| A5097204 | A5096204 | 4     | Right | 81          |
| A5097205 | A5096205 | 5     | Right | 91          |
| A5097206 | A5096206 | 6     | Right | 101         |
| A5097207 | A5096207 | 7     | Right | 111         |
| A5097208 | A5096208 | 8     | Right | 121         |
| A5097209 | A5096209 | 9     | Right | 131         |
| A5097210 | A5096210 | 10    | Right | 141         |
| A5097211 | A5096211 | 11    | Right | 151         |
| A5097212 | A5096212 | 12    | Right | 161         |

# Large Fragment Locking Plates


## 2,7/3,5/4,0 mm Distal Fibula Locking Plate 2



**Intended use :** Temporary fixation of fractures, osteotomies and nonunions of the metaphyseal and diaphyseal region of the distal fibula.

| SS        | TI        | HOLES | SIDE  | LENGTH (mm) |
|-----------|-----------|-------|-------|-------------|
| A14319913 | A14318813 | 3     | Left  | 75,5        |
| A14319914 | A14318814 | 4     | Left  | 87,5        |
| A14319905 | A14318805 | 5     | Left  | 99,5        |
| A14319906 | A14318806 | 6     | Left  | 111,5       |
| A14319907 | A14318807 | 7     | Left  | 123,5       |
| A14319908 | A14318808 | 8     | Left  | 135,5       |
| A14319909 | A14318809 | 9     | Left  | 147,5       |
| A14319910 | A14318810 | 10    | Left  | 159,5       |
| A14319911 | A14318811 | 11    | Left  | 171,5       |
| A14319912 | A14318812 | 12    | Left  | 183,5       |
| A14329913 | A14328813 | 3     | Right | 75,5        |
| A14329914 | A14328814 | 4     | Right | 87,5        |
| A14329905 | A14328805 | 5     | Right | 99,5        |
| A14329906 | A14328806 | 6     | Right | 111,5       |
| A14329907 | A14328807 | 7     | Right | 123,5       |
| A14329908 | A14328808 | 8     | Right | 135,5       |
| A14329909 | A14328809 | 9     | Right | 147,5       |
| A14329910 | A14328810 | 10    | Right | 159,5       |
| A14329911 | A14328811 | 11    | Right | 171,5       |
| A14329912 | A14328812 | 12    | Right | 183,5       |

## 4,0/4,5/5,0 mm Reconstruction Locking Compression Plate




**Intended use :** Temporary fixation of fractures, osteotomies, nonunions and malunions of large bones, including pelvis, femur, tibia & humerus in Adult patients.

| SS       | TI       | HOLES | LENGTH (mm) |
|----------|----------|-------|-------------|
| A1810905 | A1810805 | 5     | 86          |
| A1810906 | A1810806 | 6     | 102         |
| A1810907 | A1810807 | 7     | 118         |
| A1810908 | A1810808 | 8     | 134         |
| A1810909 | A1810809 | 9     | 150         |
| A1810910 | A1810810 | 10    | 166         |
| A1810912 | A1810812 | 12    | 198         |

# Large Fragment Locking Plates


## 4,0/4,5/5,0 mm Curved Reconstruction Locking Plate



**Intended use :** Temporary fixation of fractures, osteotomies, malunions and non-unions of pelvis.

| SS               | TI               | HOLES | LENGTH (mm) |
|------------------|------------------|-------|-------------|
| A101 02 228 0002 | A111 02 228 0002 | 2     | 30          |
| A101 02 228 0003 | A111 02 228 0003 | 3     | 45          |
| A101 02 228 0004 | A111 02 228 0004 | 4     | 61          |
| A101 02 228 0005 | A111 02 228 0005 | 5     | 78          |
| A101 02 228 0006 | A111 02 228 0006 | 6     | 94          |
| A101 02 228 0007 | A111 02 228 0007 | 7     | 109         |
| A101 02 228 0008 | A111 02 228 0008 | 8     | 125         |
| A101 02 228 0009 | A111 02 228 0009 | 9     | 140         |
| A101 02 228 0010 | A111 02 228 0010 | 10    | 155         |
| A101 02 228 0011 | A111 02 228 0011 | 11    | 169         |
| A101 02 228 0012 | A111 02 228 0012 | 12    | 184         |


## 3,5 mm Calcaneal Locking Plate 1



**Intended use :** Temporary fixation of fractures, osteotomies, malunions and nonunions of calcaneus in adult & skeletally mature adolescents.

| SS       | TI       | SIDE  | LENGTH (mm) |
|----------|----------|-------|-------------|
| A5989100 | A5988100 | Left  | 62,6        |
| A5989101 | A5988101 | Right | 62,6        |

## 3,5 mm Calcaneal Locking Plate 2




**Intended use :** Temporary fixation of fractures, osteotomies, malunions and nonunions of calcaneus in adult & skeletally mature adolescents.

| SS       | TI       | SIDE  | LENGTH (mm) |
|----------|----------|-------|-------------|
| A5999100 | A5998100 | Left  | 60,5        |
| A5999101 | A5998101 | Right | 60,5        |

# Large Fragment Locking Plates


## 3,5 mm Calcaneal Locking Plate 3



| SS       | TI       | SIDE  | LENGTH (mm) |
|----------|----------|-------|-------------|
| A5979100 | A5978100 | Left  | 45,7        |
| A5979101 | A5978101 | Right | 45,7        |

**Intended use :** Temporary fixation of fractures, osteotomies, malunions and nonunions of calcaneus in adult & skeletally mature adolescents.


## 3,5 mm Calcaneal Locking Plate 4



| SS       | TI       | SIDE  | LENGTH (mm) |
|----------|----------|-------|-------------|
| A5969100 | A5968100 | Left  | 65,6        |
| A5969101 | A5968101 | Right | 65,6        |

**Intended use :** Temporary fixation of fractures, osteotomies, malunions and nonunions of calcaneus in adult & skeletally mature adolescents.

## 3,5 mm Calcaneal Locking Plate 5




| SS       | TI       | SIDE  | LENGTH (mm) |
|----------|----------|-------|-------------|
| A5959100 | A5958100 | Left  | 59,1        |
| A5959101 | A5958101 | Right | 59,1        |

**Intended use :** Temporary fixation of fractures, osteotomies, malunions and nonunions of calcaneus in adult & skeletally mature adolescents.



# Large Fragment Non-Locking Plates


## 4,5 mm Proximal Femur Compression Plate 1



**Intended use :** Temporary fixation of fractures of the trochanteric region, fractures of the proximal femur, metastatic fractures of the proximal femur, osteotomies of the proximal femur and fixation of nonunions or malunions of proximal femur in adult patients.

| SS               | TI               | HOLES | SIDE  | LENGTH (mm) |
|------------------|------------------|-------|-------|-------------|
| A100 02 252 0005 | A110 02 252 0005 | 5     | Left  | 126,8       |
| A100 02 252 0007 | A110 02 252 0007 | 7     | Left  | 162,8       |
| A100 02 252 0009 | A110 02 252 0009 | 9     | Left  | 198,8       |
| A100 02 252 0011 | A110 02 252 0011 | 11    | Left  | 234,8       |
| A100 02 252 0013 | A110 02 252 0013 | 13    | Left  | 270,8       |
| A100 02 253 0005 | A110 02 253 0005 | 5     | Right | 126,8       |
| A100 02 253 0007 | A110 02 253 0007 | 7     | Right | 162,8       |
| A100 02 253 0009 | A110 02 253 0009 | 9     | Right | 198,8       |
| A100 02 253 0011 | A110 02 253 0011 | 11    | Right | 234,8       |
| A100 02 253 0013 | A110 02 253 0013 | 13    | Right | 270,8       |

## 4,5 mm Distal Lateral Femur Plate 1




**Intended use :** Temporary fixation of fractures, osteotomies, nonunions or malunions of the distal femur (intercondylar, supracondylar, intra-extra articular) in adult patients. Those plates are used by lateral approach.

| SS               | TI               | HOLES | SIDE  | LENGTH (mm) |
|------------------|------------------|-------|-------|-------------|
| A100 02 245 0005 | A110 02 245 0005 | 5     | Left  | 156,4       |
| A100 02 245 0007 | A110 02 245 0007 | 7     | Left  | 196,4       |
| A100 02 245 0009 | A110 02 245 0009 | 9     | Left  | 236,4       |
| A100 02 245 0011 | A110 02 245 0011 | 11    | Left  | 276,4       |
| A100 02 245 0013 | A110 02 245 0013 | 13    | Left  | 306,4       |
| A100 02 246 0005 | A110 02 246 0005 | 5     | Right | 156,4       |
| A100 02 246 0007 | A110 02 246 0007 | 7     | Right | 196,4       |
| A100 02 246 0009 | A110 02 246 0009 | 9     | Right | 236,4       |
| A100 02 246 0011 | A110 02 246 0011 | 11    | Right | 276,4       |
| A100 02 246 0013 | A110 02 246 0013 | 13    | Right | 306,4       |

## Large Fragment Non-Locking Plates


### 4,5 mm Distal Lateral Femur Plate 2



**Intended use :** Temporary fixation of fractures, osteotomies, nonunions or malunions of the distal femur (intercondylar, supracondylar, intra-extra articular) in adult patients. Those plates are used by lateral approach.

| SS       | TI       | HOLES | SIDE  | LENGTH (mm) |
|----------|----------|-------|-------|-------------|
| A5729106 | A5728106 | 6     | Left  | 158         |
| A5729107 | A5728107 | 7     | Left  | 174         |
| A5729108 | A5728108 | 8     | Left  | 191         |
| A5729109 | A5728109 | 9     | Left  | 206         |
| A5729110 | A5728110 | 10    | Left  | 223         |
| A5729112 | A5728112 | 12    | Left  | 250         |
| A5729206 | A5728206 | 6     | Right | 158         |
| A5729207 | A5728207 | 7     | Right | 174         |
| A5729208 | A5728208 | 8     | Right | 191         |
| A5729209 | A5728209 | 9     | Right | 206         |
| A5729210 | A5728210 | 10    | Right | 223         |
| A5729212 | A5728212 | 12    | Right | 250         |

### 4,5/6,5 mm Femur Condylar Buttress Plate



**Intended use :** Temporary fixation of fractures, osteotomies, nonunions or malunions of the distal femur (intercondylar, supracondylar, intra-extra articular) in adult patients. Those plates are used by lateral approach.


| SS               | TI               | HOLES | SIDE  | LENGTH (mm) |
|------------------|------------------|-------|-------|-------------|
| A100 02 185 0005 | A110 02 185 0005 | 5     | Left  | 133         |
| A100 02 185 0006 | A110 02 185 0006 | 6     | Left  | 149         |
| A100 02 185 0007 | A110 02 185 0007 | 7     | Left  | 165         |
| A100 02 185 0008 | A110 02 185 0008 | 8     | Left  | 181         |
| A100 02 185 0009 | A110 02 185 0009 | 9     | Left  | 197         |
| A100 02 185 0010 | A110 02 185 0010 | 10    | Left  | 213         |
| A100 02 185 0011 | A110 02 185 0011 | 11    | Left  | 229         |
| A100 02 185 0012 | A110 02 185 0012 | 12    | Left  | 245         |
| A100 02 185 0013 | A110 02 185 0013 | 13    | Left  | 264         |
| A100 02 185 0014 | A110 02 185 0014 | 14    | Left  | 277         |
| A100 02 185 0015 | A110 02 185 0015 | 15    | Left  | 293         |
| A100 02 186 0005 | A110 02 186 0005 | 5     | Right | 133         |
| A100 02 186 0006 | A110 02 186 0006 | 6     | Right | 149         |
| A100 02 186 0007 | A110 02 186 0007 | 7     | Right | 165         |
| A100 02 186 0008 | A110 02 186 0008 | 8     | Right | 181         |
| A100 02 186 0009 | A110 02 186 0009 | 9     | Right | 197         |
| A100 02 186 0010 | A110 02 186 0010 | 10    | Right | 213         |
| A100 02 186 0011 | A110 02 186 0011 | 11    | Right | 229         |
| A100 02 186 0012 | A110 02 186 0012 | 12    | Right | 245         |
| A100 02 186 0013 | A110 02 186 0013 | 13    | Right | 264         |
| A100 02 186 0014 | A110 02 186 0014 | 14    | Right | 277         |
| A100 02 186 0015 | A110 02 186 0015 | 15    | Right | 293         |

## Large Fragment Non-Locking Plates

### 4,5 mm Distal Medial Femur Buttress Plate



| SS       | TI       | HOLES | SIDE  | LENGTH (mm) |
|----------|----------|-------|-------|-------------|
| A5739104 | A5738104 | 4     | Left  | 115,9       |
| A5739106 | A5738106 | 6     | Left  | 147,9       |
| A5739108 | A5738108 | 8     | Left  | 179,9       |
| A5739204 | A5738204 | 4     | Right | 115,9       |
| A5739206 | A5738206 | 6     | Right | 147,9       |
| A5739208 | A5738208 | 8     | Right | 179,9       |




**Intended use :** Temporary fixation of fractures, osteotomies, nonunions or malunions of the distal femur (intercondylar, supracondylar, intra-extra articular) in adult patients. Those plates are used by medial approach.

### 4,5/6,5 mm Broad Compression Plate



| SS               | TI               | HOLES | LENGTH (mm) |
|------------------|------------------|-------|-------------|
| A100 02 161 0004 | A110 02 161 0004 | 4     | 71          |
| A100 02 161 0005 | A110 02 161 0005 | 5     | 87          |
| A100 02 161 0006 | A110 02 161 0006 | 6     | 103         |
| A100 02 161 0007 | A110 02 161 0007 | 7     | 119         |
| A100 02 161 0008 | A110 02 161 0008 | 8     | 135         |
| A100 02 161 0009 | A110 02 161 0009 | 9     | 151         |
| A100 02 161 0010 | A110 02 161 0010 | 10    | 167         |
| A100 02 161 0011 | A110 02 161 0011 | 11    | 183         |
| A100 02 161 0012 | A110 02 161 0012 | 12    | 199         |
| A100 02 161 0013 | A110 02 161 0013 | 13    | 215         |
| A100 02 161 0014 | A110 02 161 0014 | 14    | 231         |
| A100 02 161 0015 | A110 02 161 0015 | 15    | 247         |
| A100 02 161 0016 | A110 02 161 0016 | 16    | 263         |
| A100 02 161 0017 | A110 02 161 0017 | 17    | 279         |
| A100 02 161 0018 | A110 02 161 0018 | 18    | 295         |




**Intended use :** Temporary fixation of fractures, malunions, nonunions, osteotomies and pseudoarthroses of the femur.

## Large Fragment Non-Locking Plates

### 4,5/6,5 mm Narrow Compression Plate



| SS               | TI               | HOLES | LENGTH (mm) |
|------------------|------------------|-------|-------------|
| A100 02 151 0002 | A110 02 151 0002 | 2     | 39          |
| A100 02 151 0003 | A110 02 151 0003 | 3     | 55          |
| A100 02 151 0004 | A110 02 151 0004 | 4     | 71          |
| A100 02 151 0005 | A110 02 151 0005 | 5     | 87          |
| A100 02 151 0006 | A110 02 151 0006 | 6     | 103         |
| A100 02 151 0007 | A110 02 151 0007 | 7     | 119         |
| A100 02 151 0008 | A110 02 151 0008 | 8     | 135         |
| A100 02 151 0009 | A110 02 151 0009 | 9     | 151         |
| A100 02 151 0010 | A110 02 151 0010 | 10    | 167         |
| A100 02 151 0011 | A110 02 151 0011 | 11    | 183         |
| A100 02 151 0012 | A110 02 151 0012 | 12    | 199         |
| A100 02 151 0013 | A110 02 151 0013 | 13    | 215         |
| A100 02 151 0014 | A110 02 151 0014 | 14    | 231         |
| A100 02 151 0015 | A110 02 151 0015 | 15    | 247         |
| A100 02 151 0016 | A110 02 151 0016 | 16    | 263         |




**Intended use :** Temporary fixation of fractures, malunions, nonunions, osteotomies and pseudoarthroses of the humerus and tibia.

### 4,5 mm Straight Reconstruction Plate



| SS               | TI               | HOLES | LENGTH (mm) |
|------------------|------------------|-------|-------------|
| A100 02 277 0004 | A110 02 277 0004 | 4     | 71          |
| A100 02 277 0005 | A110 02 277 0005 | 5     | 87          |
| A100 02 277 0006 | A110 02 277 0006 | 6     | 103         |
| A100 02 277 0007 | A110 02 277 0007 | 7     | 119         |
| A100 02 277 0008 | A110 02 277 0008 | 8     | 135         |
| A100 02 277 0009 | A110 02 277 0009 | 9     | 151         |
| A100 02 277 0010 | A110 02 277 0010 | 10    | 167         |
| A100 02 227 0011 | A110 02 227 0011 | 11    | 183         |
| A100 02 227 0012 | A110 02 227 0012 | 12    | 199         |
| A100 02 227 0013 | A110 02 227 0013 | 13    | 215         |
| A100 02 227 0014 | A110 02 227 0014 | 14    | 231         |
| A100 02 227 0015 | A110 02 227 0015 | 15    | 247         |
| A100 02 227 0016 | A110 02 227 0016 | 16    | 263         |



**Intended use :** Temporary fixation of fractures, osteotomies, nonunions and malunions of large bones, including pelvis, femur, tibia & humerus in Adult patients.

# Large Fragment Non-Locking Plates


## 4,5 mm Curved Reconstruction Plate



**Intended use :** Temporary fixation of fractures, malunions, nonunions, osteotomies of the pelvis.

| SS               | TI               | HOLES | LENGTH (mm) |
|------------------|------------------|-------|-------------|
| A100 02 228 0004 | A110 02 228 0004 | 4     | 61,4        |
| A100 02 228 0005 | A110 02 228 0005 | 5     | 77,3        |
| A100 02 228 0006 | A110 02 228 0006 | 6     | 93,2        |
| A100 02 228 0007 | A110 02 228 0007 | 7     | 108,9       |
| A100 02 228 0008 | A110 02 228 0008 | 8     | 124,7       |
| A100 02 228 0009 | A110 02 228 0009 | 9     | 139,8       |
| A100 02 228 0010 | A110 02 228 0010 | 10    | 154,9       |
| A100 02 228 0011 | A110 02 228 0011 | 11    | 169,6       |
| A100 02 228 0012 | A110 02 228 0012 | 12    | 184         |
| A100 02 228 0013 | A110 02 228 0013 | 13    | 198         |
| A100 02 228 0014 | A110 02 228 0014 | 14    | 211,6       |
| A100 02 228 0015 | A110 02 228 0015 | 15    | 224,7       |
| A100 02 228 0016 | A110 02 228 0016 | 16    | 237,5       |

## 4,5 mm Semi Tubular Plate



**Intended use :** Temporary fixation of fractures, malunions, nonunions, osteotomies of the fibula, ulna, radius.


| SS               | TI               | HOLES | LENGTH (mm) |
|------------------|------------------|-------|-------------|
| A100 02 180 0004 | A110 02 180 0004 | 4     | 71          |
| A100 02 180 0005 | A110 02 180 0005 | 5     | 87          |
| A100 02 180 0006 | A110 02 180 0006 | 6     | 103         |
| A100 02 180 0007 | A110 02 180 0007 | 7     | 119         |
| A100 02 180 0008 | A110 02 180 0008 | 8     | 135         |
| A100 02 180 0009 | A110 02 180 0009 | 9     | 151         |
| A100 02 180 0010 | A110 02 180 0010 | 10    | 167         |
| A100 02 180 0011 | A110 02 180 0011 | 11    | 183         |
| A100 02 180 0012 | A110 02 180 0012 | 12    | 231         |

## Large Fragment Non-Locking Plates

### 4,5 mm Spoon Plate



| SS               | HOLES | LENGTH (mm) |
|------------------|-------|-------------|
| A100 02 260 0005 | 5     | 101,4       |
| A100 02 260 0006 | 6     | 120,6       |




**Intended use :** Temporary fixation of fractures, malunions, nonunions, osteotomies of proximal and distal tibia.

### 4,5/6,5 mm Cobra Head Plate



| SS               | HOLES | LENGTH (mm) |
|------------------|-------|-------------|
| A200 02 180 0006 | 6     | 143         |
| A200 02 180 0007 | 7     | 159         |
| A200 02 180 0008 | 8     | 175         |
| A200 02 180 0009 | 9     | 191         |
| A200 02 180 0010 | 10    | 207         |
| A200 02 180 0011 | 11    | 223         |
| A200 02 180 0012 | 12    | 239         |
| A200 02 180 0013 | 13    | 255         |




**Intended use :** Temporary fixation for the purpose of hip arthrodesis.

# Large Fragment Non-Locking Plates

## 4,5/6,5 mm T-Buttress Plate



| SS               | TI               | HOLES | LENGTH (mm) |
|------------------|------------------|-------|-------------|
| A101 02 190 0003 | A111 02 190 0003 | 3     | 63,5        |
| A101 02 190 0004 | A111 02 190 0004 | 4     | 79,5        |
| A101 02 190 0005 | A111 02 190 0005 | 5     | 95,5        |
| A101 02 190 0006 | A111 02 190 0006 | 6     | 111,5       |
| A101 02 190 0007 | A111 02 190 0007 | 7     | 127,5       |
| A101 02 190 0008 | A111 02 190 0008 | 8     | 143,5       |




**Intended use :** Temporary fixation of fractures, malunions, nonunions, osteotomies of the proximal tibia in adult patients. Those plates are used both by medial and lateral approach.

## 4,5/6,5 mm L-Buttress Plate




| SS               | TI               | HOLES | SIDE  | LENGTH (mm) |
|------------------|------------------|-------|-------|-------------|
| A100 02 203 0003 | A110 02 203 0003 | 3     | Left  | 65,1        |
| A100 02 203 0004 | A110 02 203 0004 | 4     | Left  | 81,1        |
| A100 02 203 0005 | A110 02 203 0005 | 5     | Left  | 97,1        |
| A100 02 203 0006 | A110 02 203 0006 | 6     | Left  | 113,1       |
| A100 02 203 0007 | A110 02 203 0007 | 7     | Left  | 129,1       |
| A100 02 203 0008 | A110 02 203 0008 | 8     | Left  | 145,1       |
| A100 02 204 0003 | A110 02 204 0003 | 3     | Right | 65,1        |
| A100 02 204 0004 | A110 02 204 0004 | 4     | Right | 81,1        |
| A100 02 204 0005 | A110 02 204 0005 | 5     | Right | 97,1        |
| A100 02 204 0006 | A110 02 204 0006 | 6     | Right | 113,1       |
| A100 02 204 0007 | A110 02 204 0007 | 7     | Right | 129,1       |
| A100 02 204 0008 | A110 02 204 0008 | 8     | Right | 145,1       |



**Intended use :** Temporary fixation of fractures, malunions, nonunions, osteotomies of the proximal tibia in adult patients. Those plates are used both by medial and lateral approach.

# Large Fragment Non-Locking Plates


## 4,5 mm Proximal Lateral Tibia Plate



| SS               | TI               | HOLES | SIDE  | LENGTH (mm) |
|------------------|------------------|-------|-------|-------------|
| A100 02 256 0005 | A110 02 256 0005 | 5     | Left  | 141,8       |
| A100 02 256 0007 | A110 02 256 0007 | 7     | Left  | 181,8       |
| A100 02 256 0009 | A110 02 256 0009 | 9     | Left  | 221,8       |
| A100 02 256 0011 | A110 02 256 0011 | 11    | Left  | 261,8       |
| A100 02 256 0013 | A110 02 256 0013 | 13    | Left  | 301,8       |
| A100 02 257 0005 | A110 02 257 0005 | 5     | Right | 141,8       |
| A100 02 257 0007 | A110 02 257 0007 | 7     | Right | 181,8       |
| A100 02 257 0009 | A110 02 257 0009 | 9     | Right | 221,8       |
| A100 02 257 0011 | A110 02 257 0011 | 11    | Right | 261,8       |
| A100 02 257 0013 | A110 02 257 0013 | 13    | Right | 301,8       |

**Intended use :** Temporary fixation of fractures, malunions, nonunions, osteotomies and pseudoarthroses of the proximal tibia. Those plates are used by lateral approach.

## 4,5 mm Distal Tibia Lateral Plate



| SS               | TI               | HOLES | SIDE  | LENGTH (mm) |
|------------------|------------------|-------|-------|-------------|
| A100 02 240 0005 | A110 02 240 0005 | 5     | Left  | 108,3       |
| A100 02 240 0007 | A110 02 240 0007 | 7     | Left  | 140,3       |
| A100 02 240 0009 | A110 02 240 0009 | 9     | Left  | 172,3       |
| A100 02 240 0011 | A110 02 240 0011 | 11    | Left  | 204,3       |
| A100 02 240 0013 | A110 02 240 0013 | 13    | Left  | 236,3       |
| A100 02 241 0005 | A110 02 241 0005 | 5     | Right | 108,3       |
| A100 02 241 0007 | A110 02 241 0007 | 7     | Right | 140,3       |
| A100 02 241 0009 | A110 02 241 0009 | 9     | Right | 172,3       |
| A100 02 241 0011 | A110 02 241 0011 | 11    | Right | 204,3       |
| A100 02 241 0013 | A110 02 241 0013 | 13    | Right | 236,3       |


**Intended use :** Temporary fixation of fractures, malunions, nonunions, osteotomies and pseudoarthroses of the distal tibia.

# Large Fragment Non-Locking Plates

## 3,5 mm Distal Tibia Medial Plate 2



| SS               | TI               | HOLES | SIDE  | LENGTH (mm) |
|------------------|------------------|-------|-------|-------------|
| A100 02 242 0005 | A110 02 242 0005 | 5     | Left  | 125,9       |
| A100 02 242 0007 | A110 02 242 0007 | 7     | Left  | 149,9       |
| A100 02 242 0009 | A110 02 242 0009 | 9     | Left  | 173,9       |
| A100 02 242 0011 | A110 02 242 0011 | 11    | Left  | 197,9       |
| A100 02 242 0013 | A110 02 242 0013 | 13    | Left  | 221,9       |
| A100 02 243 0005 | A110 02 243 0005 | 5     | Right | 125,9       |
| A100 02 243 0007 | A110 02 243 0007 | 7     | Right | 149,9       |
| A100 02 243 0009 | A110 02 243 0009 | 9     | Right | 173,9       |
| A100 02 243 0011 | A110 02 243 0011 | 11    | Right | 197,9       |
| A100 02 243 0013 | A110 02 243 0013 | 13    | Right | 221,9       |




**Intended use :** Temporary fixation of fractures, malunions, nonunions, osteotomies and pseudoarthroses of the distal tibia. Those plates are used by medial approach.

## 3,5 mm Distal Fibula Plate 1



| SS       | TI       | HOLES | SIDE  | LENGTH (mm) |
|----------|----------|-------|-------|-------------|
| A5099102 | A5098102 | 2     | Left  | 61          |
| A5099104 | A5098104 | 4     | Left  | 81          |
| A5099105 | A5098105 | 5     | Left  | 91          |
| A5099106 | A5098106 | 6     | Left  | 101         |
| A5099107 | A5098107 | 7     | Left  | 111         |
| A5099108 | A5098108 | 8     | Left  | 121         |
| A5099109 | A5098109 | 9     | Left  | 131         |
| A5099110 | A5098110 | 10    | Left  | 141         |
| A5099202 | A5098202 | 2     | Right | 61          |
| A5099204 | A5098204 | 4     | Right | 81          |
| A5099205 | A5098205 | 5     | Right | 91          |
| A5099206 | A5098206 | 6     | Right | 101         |
| A5099207 | A5098207 | 7     | Right | 111         |
| A5099208 | A5098208 | 8     | Right | 121         |
| A5099209 | A5098209 | 9     | Right | 131         |
| A5099210 | A5098210 | 10    | Right | 141         |



**Intended use :** Temporary fixation of fractures, malunions, nonunions, osteotomies and pseudoarthroses of the distal fibula.

"We supply WORLD's  
Trauma Implants..."

# Dhs/Dcs Locking/Non-Locking Plates




**aysam**

# DHS-DCS Plates & Lag-Compression Screws

## DHS Locking Plate 135°




| SS       | TI       | HOLES | LENGTH (mm) |
|----------|----------|-------|-------------|
| A5839102 | A5838102 | 2     | 91,3        |
| A5839103 | A5838103 | 3     | 107,3       |
| A5839104 | A5838104 | 4     | 123,3       |
| A5839105 | A5838105 | 5     | 139,3       |
| A5839106 | A5838106 | 6     | 153,3       |
| A5839107 | A5838107 | 7     | 171,3       |
| A5839108 | A5838108 | 8     | 187,3       |
| A5839109 | A5838109 | 9     | 203,3       |
| A5839110 | A5838110 | 10    | 219,3       |
| A5839111 | A5838111 | 11    | 235,3       |
| A5839112 | A5838112 | 12    | 251,3       |




**Intended use :** Temporary fixation of fractures, malunions, nonunions, osteotomies and pseudoarthroses of the proximal femur (intertrochanteric fractures, subtrochanteric fractures, basilar neck fractures).

## DCS Locking Plate 95°



| SS       | TI       | HOLES | LENGTH (mm) |
|----------|----------|-------|-------------|
| A5829102 | A5828102 | 2     | 86,8        |
| A5829103 | A5828103 | 3     | 102,8       |
| A5829104 | A5828104 | 4     | 118,8       |
| A5829105 | A5828105 | 5     | 134,8       |
| A5829106 | A5828106 | 6     | 150,8       |
| A5829107 | A5828107 | 7     | 166,8       |
| A5829108 | A5828108 | 8     | 182,8       |
| A5829109 | A5828109 | 9     | 198,8       |
| A5829110 | A5828110 | 10    | 214,8       |
| A5829111 | A5828111 | 11    | 230,8       |
| A5829112 | A5828112 | 12    | 246,8       |




**Intended use :** Temporary fixation of fractures, malunions, nonunions, osteotomies and pseudoarthroses of the distal femur (Supracondylar, intercondylar fractures).

# DHS-DCS Plates & Lag-Compression Screws

## DHS Plate 135°



| SS               | TI               | HOLES | LENGTH (mm) |
|------------------|------------------|-------|-------------|
| A100 02 310 0002 | A110 02 310 0002 | 2     | 91,4        |
| A100 02 310 0003 | A110 02 310 0003 | 3     | 107,4       |
| A100 02 310 0004 | A110 02 310 0004 | 4     | 123,4       |
| A100 02 310 0005 | A110 02 310 0005 | 5     | 139,4       |
| A100 02 310 0006 | A110 02 310 0006 | 6     | 155,4       |
| A100 02 310 0007 | A110 02 310 0007 | 7     | 171,4       |
| A100 02 310 0008 | A110 02 310 0008 | 8     | 187,4       |
| A100 02 310 0009 | A110 02 310 0009 | 9     | 203,4       |
| A100 02 310 0010 | A110 02 310 0010 | 10    | 219,4       |
| A100 02 310 0011 | A110 02 310 0011 | 11    | 235,4       |
| A100 02 310 0012 | A110 02 310 0012 | 12    | 251,4       |




**Intended use :** Temporary fixation of fractures, malunions, nonunions, osteotomies and pseudoarthroses of the proximal femur (intertrochanteric fractures, subtrochanteric fractures, basilar neck fractures).

## DCS Plate 95°




| SS               | TI               | HOLES | LENGTH (mm) |
|------------------|------------------|-------|-------------|
| A100 02 314 0002 | A110 02 314 0002 | 2     | 86,7        |
| A100 02 314 0003 | A110 02 314 0003 | 3     | 102,7       |
| A100 02 314 0004 | A110 02 314 0004 | 4     | 118,7       |
| A100 02 314 0005 | A110 02 314 0005 | 5     | 134,7       |
| A100 02 314 0006 | A110 02 314 0006 | 6     | 150,7       |
| A100 02 314 0007 | A110 02 314 0007 | 7     | 166,7       |
| A100 02 314 0008 | A110 02 314 0008 | 8     | 182,7       |
| A100 02 314 0009 | A110 02 314 0009 | 9     | 198,7       |
| A100 02 314 0010 | A110 02 314 0010 | 10    | 214,7       |
| A100 02 314 0011 | A110 02 314 0011 | 11    | 230,7       |
| A100 02 314 0012 | A110 02 314 0012 | 12    | 246,7       |



**Intended use :** Temporary fixation of fractures, malunions, nonunions, osteotomies and pseudoarthroses of the distal femur (Supracondylar, intercondylar fractures).

# DHS-DCS Plates & Lag-Compression Screws

## DHS/DCS Lag Screw



**Intended use :** Temporary fixation and compression at fracture site.  
Used together with compression screw, DHS/DCS plates.

| SS               | TI               | LENGTH (mm) |
|------------------|------------------|-------------|
| A100 03 011 0500 | A110 03 011 0500 | 50          |
| A100 03 011 0550 | A110 03 011 0550 | 55          |
| A100 03 011 0600 | A110 03 011 0600 | 60          |
| A100 03 011 0650 | A110 03 011 0650 | 65          |
| A100 03 011 0700 | A110 03 011 0700 | 70          |
| A100 03 011 0750 | A110 03 011 0750 | 75          |
| A100 03 011 0800 | A110 03 011 0800 | 80          |
| A100 03 011 0850 | A110 03 011 0850 | 85          |
| A100 03 011 0900 | A110 03 011 0900 | 90          |
| A100 03 011 0950 | A110 03 011 0950 | 95 mm       |
| A100 03 011 1000 | A110 03 011 1000 | 100 mm      |
| A100 03 011 1050 | A110 03 011 1050 | 105 mm      |
| A100 03 011 1100 | A110 03 011 1100 | 110 mm      |
| A100 03 011 1150 | A110 03 011 1150 | 115 mm      |
| A100 03 011 1200 | A110 03 011 1200 | 120 mm      |
| A100 03 011 1250 | A110 03 011 1250 | 125 mm      |

## DHS/DCS Compression Screw



**Intended use :** Compression at fracture site.  
Used together with lag screw, DHS/DCS plates.

| SS               | TI       | LENGTH (mm) |
|------------------|----------|-------------|
| A100 03 011 1051 | A5878050 | 28          |

"We supply WORLD's  
Orthopedic & Trauma  
Implants..."


# Large Fragment Screws

Locking/Non-Locking



# Large Fragment Locking Screws

## 5,0 mm Self-Tapping Locking Screw



| SS       | TI       | DIA. x LENGTH (mm) |
|----------|----------|--------------------|
| A2095014 | A2085014 | 5,0X14             |
| A2095016 | A2085016 | 5,0X16             |
| A2095018 | A2085018 | 5,0X18             |
| A2095020 | A2085020 | 5,0X20             |
| A2095022 | A2085022 | 5,0X22             |
| A2095024 | A2085024 | 5,0X24             |
| A2095026 | A2085026 | 5,0X26             |
| A2095028 | A2085028 | 5,0X28             |



| SS       | TI       | DIA. x LENGTH (mm) |
|----------|----------|--------------------|
| A2095030 | A2085030 | 5,0X30             |
| A2095032 | A2085032 | 5,0X32             |
| A2095034 | A2085034 | 5,0X34             |
| A2095036 | A2085036 | 5,0X36             |
| A2095038 | A2085038 | 5,0X38             |
| A2095040 | A2085040 | 5,0X40             |
| A2095042 | A2085042 | 5,0X42             |
| A2095044 | A2085044 | 5,0X44             |
| A2095046 | A2085046 | 5,0X46             |
| A2095048 | A2085048 | 5,0X48             |
| A2095050 | A2085050 | 5,0X50             |
| A2095055 | A2085055 | 5,0X55             |
| A2095060 | A2085060 | 5,0X60             |
| A2095065 | A2085065 | 5,0X65             |
| A2095070 | A2085070 | 5,0X70             |
| A2095075 | A2085075 | 5,0X75             |
| A2095080 | A2085080 | 5,0X80             |
| A2095085 | A2085085 | 5,0X85             |
| A2095090 | A2085090 | 5,0X90             |

## 5,0 mm Locking Cancellous Screw Full Threaded



| SS        | TI        | DIA. x LENGTH (mm) |
|-----------|-----------|--------------------|
| A33802035 | A33702035 | 5,0X35             |
| A33802040 | A33702040 | 5,0X40             |
| A33802045 | A33702045 | 5,0X45             |
| A33802050 | A33702050 | 5,0X50             |
| A33802055 | A33702055 | 5,0X55             |
| A33802060 | A33702060 | 5,0X60             |
| A33802065 | A33702065 | 5,0X65             |
| A33802070 | A33702070 | 5,0X70             |
| A33802075 | A33702075 | 5,0X75             |
| A33802080 | A33702080 | 5,0X80             |
| A33802085 | A33702085 | 5,0X85             |
| A33802090 | A33702090 | 5,0X90             |
| A33802095 | A33702095 | 5,0X95             |
| A33802100 | A33702100 | 5,0X100            |
| A33802105 | A33702105 | 5,0X105            |
| A33802110 | A33702110 | 5,0X110            |


## Large Fragment Locking Screws

### 5,0 mm Locking Cancellous Screw Partially Threaded



| SS        | TI        | DIA. x LENGTH (mm) |
|-----------|-----------|--------------------|
| A33801035 | A33701035 | 5,0X35             |
| A33801040 | A33701040 | 5,0X40             |
| A33801045 | A33701045 | 5,0X45             |
| A33801050 | A33701050 | 5,0X50             |
| A33801055 | A33701055 | 5,0X55             |
| A33801060 | A33701060 | 5,0X60             |
| A33801065 | A33701065 | 5,0X65             |
| A33801070 | A33701070 | 5,0X70             |
| A33801075 | A33701075 | 5,0X75             |
| A33801080 | A33701080 | 5,0X80             |
| A33801085 | A33701085 | 5,0X85             |
| A33801090 | A33701090 | 5,0X90             |
| A33801095 | A33701095 | 5,0X95             |
| A33801100 | A33701100 | 5,0X100            |
| A33801105 | A33701105 | 5,0X105            |
| A33801110 | A33701110 | 5,0X110            |

### 5,0 mm Self-Drilling Locking Screw



| SS       | TI       | DIA. x LENGTH (mm) |
|----------|----------|--------------------|
| A2495018 | A2485018 | 5,0X18             |
| A2495020 | A2485020 | 5,0X20             |
| A2495022 | A2485022 | 5,0X22             |
| A2495024 | A2485024 | 5,0X24             |
| A2495026 | A2485026 | 5,0X26             |
| A2495028 | A2485028 | 5,0X28             |
| A2495030 | A2485030 | 5,0X30             |


| SS       | TI       | DIA. x LENGTH (mm) |
|----------|----------|--------------------|
| A2495032 | A2485032 | 5,0X32             |
| A2495034 | A2485034 | 5,0X34             |
| A2495036 | A2485036 | 5,0X36             |
| A2495038 | A2485038 | 5,0X38             |
| A2495040 | A2485040 | 5,0X40             |
| A2495042 | A2485042 | 5,0X42             |
| A2495044 | A2485044 | 5,0X44             |
| A2495046 | A2485046 | 5,0X46             |
| A2495048 | A2485048 | 5,0X48             |
| A2495050 | A2485050 | 5,0X50             |
| A2495055 | A2485055 | 5,0X55             |
| A2495060 | A2485060 | 5,0X60             |
| A2495065 | A2485065 | 5,0X65             |
| A2495070 | A2485070 | 5,0X70             |
| A2495075 | A2485075 | 5,0X75             |
| A2495080 | A2485080 | 5,0X80             |
| A2495085 | A2485085 | 5,0X85             |
| A2495090 | A2485090 | 5,0X90             |

# Large Fragment Locking Screws

## 5,0 mm Cable Locking Screw




**TI**  
A900 03 020 0008



## Locking Screw Plug




**TI**  
A900 03 020 0009



## 6,5 mm Locking Cancellous Screw Full Threaded



| SS        | TI        | DIA. x LENGTH<br>(mm) |
|-----------|-----------|-----------------------|
| A34865035 | A34765035 | 6,5X35                |
| A34865040 | A34765040 | 6,5X40                |
| A34865045 | A34765045 | 6,5X45                |
| A34865050 | A34765050 | 6,5X50                |
| A34865055 | A34765055 | 6,5X55                |
| A34865060 | A34765060 | 6,5X60                |
| A34865065 | A34765065 | 6,5X65                |
| A34865070 | A34765070 | 6,5X70                |
| A34865075 | A34765075 | 6,5X75                |
| A34865080 | A34765080 | 6,5X80                |
| A34865085 | A34765085 | 6,5X85                |
| A34865090 | A34765090 | 6,5X90                |
| A34865095 | A34765095 | 6,5X95                |
| A34865100 | A34765100 | 6,5X100               |
| A34865105 | A34765105 | 6,5X105               |
| A34865110 | A34765110 | 6,5X110               |



\* DIA. = DIAMETER (mm) \*\* LENGTH = LENGTH (mm)

## Large Fragment Locking Screws

### 6,5 mm Locking Cancellous Screw 16 mm Threaded



| SS        | Ti        | DIA. x LENGTH (mm) |
|-----------|-----------|--------------------|
| A34867035 | A34767035 | 6,5X35             |
| A34867040 | A34767040 | 6,5X40             |
| A34867045 | A34767045 | 6,5X45             |
| A34867050 | A34767050 | 6,5X50             |
| A34867055 | A34767055 | 6,5X55             |
| A34867060 | A34767060 | 6,5X60             |
| A34867065 | A34767065 | 6,5X65             |
| A34867070 | A34767070 | 6,5X70             |
| A34867075 | A34767075 | 6,5X75             |
| A34867080 | A34767080 | 6,5X80             |
| A34867085 | A34767085 | 6,5X85             |
| A34867090 | A34767090 | 6,5X90             |
| A34867095 | A34767095 | 6,5X95             |
| A34867100 | A34767100 | 6,5X100            |
| A34867105 | A34767105 | 6,5X105            |
| A34867110 | A34767110 | 6,5X110            |

### 6,5 mm Locking Cancellous Screw 32 mm Threaded



| SS        | Ti        | DIA. x LENGTH (mm) |
|-----------|-----------|--------------------|
| A34866045 | A34766045 | 6,5X45             |
| A34866050 | A34766050 | 6,5X50             |
| A34866055 | A34766055 | 6,5X55             |
| A34866060 | A34766060 | 6,5X60             |
| A34866065 | A34766065 | 6,5X65             |
| A34866070 | A34766070 | 6,5X70             |
| A34866075 | A34766075 | 6,5X75             |
| A34866080 | A34766080 | 6,5X80             |
| A34866085 | A34766085 | 6,5X85             |
| A34866090 | A34766090 | 6,5X90             |
| A34866095 | A34766095 | 6,5X95             |
| A34866100 | A34766100 | 6,5X100            |
| A34866105 | A34766105 | 6,5X105            |
| A34866110 | A34766110 | 6,5X110            |

# Large Fragment Locking Screws

## 5,0 mm Self-Tapping Locking Cannulated Screw



| SS       | TI       | DIA. x LENGTH (mm) | SS       | TI       | DIA. x LENGTH (mm) |
|----------|----------|--------------------|----------|----------|--------------------|
| A2045014 | A2055014 | 5,0X14             | A2045032 | A2055032 | 5,0X32             |
| A2045016 | A2055016 | 5,0X16             | A2045034 | A2055034 | 5,0X34             |
| A2045018 | A2055018 | 5,0X18             | A2045036 | A2055036 | 5,0X36             |
| A2045020 | A2055020 | 5,0X20             | A2045038 | A2055038 | 5,0X38             |
| A2045022 | A2055022 | 5,0X22             | A2045040 | A2055040 | 5,0X40             |
| A2045024 | A2055024 | 5,0X24             | A2045042 | A2055042 | 5,0X42             |
| A2045026 | A2055026 | 5,0X26             | A2045044 | A2055044 | 5,0X44             |
| A2045028 | A2055028 | 5,0X28             | A2045046 | A2055046 | 5,0X46             |
| A2045030 | A2055030 | 5,0X30             | A2045048 | A2055048 | 5,0X48             |
|          |          |                    | A2045050 | A2055050 | 5,0X50             |
|          |          |                    | A2045055 | A2055055 | 5,0X55             |
|          |          |                    | A2045060 | A2055060 | 5,0X60             |
|          |          |                    | A2045065 | A2055065 | 5,0X65             |
|          |          |                    | A2045070 | A2055070 | 5,0X70             |
|          |          |                    | A2045075 | A2055075 | 5,0X75             |
|          |          |                    | A2045080 | A2055080 | 5,0X80             |
|          |          |                    | A2045085 | A2055085 | 5,0X85             |
|          |          |                    | A2045090 | A2055090 | 5,0X90             |



## 4,5 mm Self-Tapping Cortical Screw



| SS               | TI               | DIA. x LENGTH (mm) | SS               | TI               | DIA. x LENGTH (mm) |
|------------------|------------------|--------------------|------------------|------------------|--------------------|
| A100 03 223 0140 | A110 03 223 0140 | 4,5X14             | A100 03 223 0400 | A110 03 223 0400 | 4,5X40             |
| A100 03 223 0160 | A110 03 223 0160 | 4,5X16             | A100 03 223 0420 | A110 03 223 0420 | 4,5X42             |
| A100 03 223 0180 | A110 03 223 0180 | 4,5X18             | A100 03 223 0440 | A110 03 223 0440 | 4,5X44             |
| A100 03 223 0200 | A110 03 223 0200 | 4,5X20             | A100 03 223 0460 | A110 03 223 0460 | 4,5X46             |
| A100 03 223 0220 | A110 03 223 0220 | 4,5X22             | A100 03 223 0480 | A110 03 223 0480 | 4,5X48             |
| A100 03 223 0240 | A110 03 223 0240 | 4,5X24             | A100 03 223 0500 | A110 03 223 0500 | 4,5X50             |
| A100 03 223 0260 | A110 03 223 0260 | 4,5X26             | A100 03 223 0520 | A110 03 223 0520 | 4,5X52             |
| A100 03 223 0280 | A110 03 223 0280 | 4,5X28             | A100 03 223 0540 | A110 03 223 0540 | 4,5X54             |
| A100 03 223 0300 | A110 03 223 0300 | 4,5X30             | A100 03 223 0560 | A110 03 223 0560 | 4,5X56             |
| A100 03 223 0320 | A110 03 223 0320 | 4,5X32             | A100 03 223 0580 | A110 03 223 0580 | 4,5X58             |
| A100 03 223 0340 | A110 03 223 0340 | 4,5X34             | A100 03 223 0600 | A110 03 223 0600 | 4,5X60             |
| A100 03 223 0360 | A110 03 223 0360 | 4,5X36             | A100 03 223 0650 | A110 03 223 0650 | 4,5X65             |
| A100 03 223 0380 | A110 03 223 0380 | 4,5X38             | A100 03 223 0700 | A110 03 223 0700 | 4,5X70             |



# Large Fragment Non-Locking Screws

## 6,5 mm Cancellous Screw 16 mm Threaded



| SS               | Ti               | DIA. x LENGTH (mm) |
|------------------|------------------|--------------------|
| A100 03 225 0500 | A110 03 225 0500 | 6,5X50             |
| A100 03 225 0550 | A110 03 225 0550 | 6,5X55             |
| A100 03 225 0600 | A110 03 225 0600 | 6,5X60             |
| A100 03 225 0650 | A110 03 225 0650 | 6,5X65             |
| A100 03 225 0700 | A110 03 225 0700 | 6,5X70             |
| A100 03 225 0750 | A110 03 225 0750 | 6,5X75             |
| A100 03 225 0800 | A110 03 225 0800 | 6,5X80             |
| A100 03 225 0850 | A110 03 225 0850 | 6,5X85             |
| A100 03 225 0900 | A110 03 225 0900 | 6,5X90             |
| A100 03 225 0950 | A110 03 225 0950 | 6,5X95             |
| A100 03 225 1000 | A110 03 225 1000 | 6,5X100            |
| A100 03 225 1050 | A110 03 225 1050 | 6,5X105            |
| A100 03 225 1100 | A110 03 225 1100 | 6,5X110            |



## 6,5 mm Cancellous Screw 32 mm Threaded



| SS               | Ti               | DIA. x LENGTH (mm) |
|------------------|------------------|--------------------|
| A100 03 226 0500 | A110 03 226 0500 | 6,5X50             |
| A100 03 226 0550 | A110 03 226 0550 | 6,5X55             |
| A100 03 226 0600 | A110 03 226 0600 | 6,5X60             |
| A100 03 226 0650 | A110 03 226 0650 | 6,5X65             |
| A100 03 226 0700 | A110 03 226 0700 | 6,5X70             |
| A100 03 226 0750 | A110 03 226 0750 | 6,5X75             |
| A100 03 226 0800 | A110 03 226 0800 | 6,5X80             |
| A100 03 226 0850 | A110 03 226 0850 | 6,5X85             |
| A100 03 226 0900 | A110 03 226 0900 | 6,5X90             |
| A100 03 226 0950 | A110 03 226 0950 | 6,5X95             |
| A100 03 226 1000 | A110 03 226 1000 | 6,5X100            |
| A100 03 226 1050 | A110 03 226 1050 | 6,5X105            |
| A100 03 226 1100 | A110 03 226 1100 | 6,5X110            |



# Large Fragment Non-Locking Screws

## 6,5 mm Cancellous Screw Full Threaded



| SS               | Ti               | DIA. x LENGTH (mm) |
|------------------|------------------|--------------------|
| A100 03 227 0500 | A110 03 227 0500 | 6,5X50             |
| A100 03 227 0550 | A110 03 227 0550 | 6,5X55             |
| A100 03 227 0600 | A110 03 227 0600 | 6,5X60             |
| A100 03 227 0650 | A110 03 227 0650 | 6,5X65             |
| A100 03 227 0700 | A110 03 227 0700 | 6,5X70             |
| A100 03 227 0750 | A110 03 227 0750 | 6,5X75             |
| A100 03 227 0800 | A110 03 227 0800 | 6,5X80             |
| A100 03 227 0850 | A110 03 227 0850 | 6,5X85             |
| A100 03 227 0900 | A110 03 227 0900 | 6,5X90             |
| A100 03 227 0950 | A110 03 227 0950 | 6,5X95             |
| A100 03 227 1000 | A110 03 227 1000 | 6,5X100            |
| A100 03 227 1050 | A110 03 227 1050 | 6,5X105            |
| A100 03 227 1100 | A110 03 227 1100 | 6,5X110            |


## 4,5 mm Malleolar Screw



| SS               | Ti               | DIA. x LENGTH (mm) |
|------------------|------------------|--------------------|
| A100 03 224 0250 | A110 03 224 0250 | 4,5X25             |
| A100 03 224 0300 | A110 03 224 0300 | 4,5X30             |
| A100 03 224 0350 | A110 03 224 0350 | 4,5X35             |
| A100 03 224 0400 | A110 03 224 0400 | 4,5X40             |
| A100 03 224 0450 | A110 03 224 0450 | 4,5X45             |
| A100 03 224 0500 | A110 03 224 0500 | 4,5X50             |
| A100 03 224 0550 | A110 03 224 0550 | 4,5X55             |
| A100 03 224 0600 | A110 03 224 0600 | 4,5X60             |
| A100 03 224 0650 | A110 03 224 0650 | 4,5X65             |
| A100 03 224 0700 | A110 03 224 0700 | 4,5X70             |

# Large Fragment Non-Locking Screws

## Washer Ø13




| SS               | DIA.<br>(mm) |
|------------------|--------------|
| A110 04 000 0014 | 13           |

## Washer Ø10



| SS               | DIA.<br>(mm) |
|------------------|--------------|
| A110 04 000 0010 | 10           |

## Washer Ø8,0



| SS        | TI        | DIA.<br>(mm) |
|-----------|-----------|--------------|
| A34004108 | A34004008 | 8            |

f aysam.orthopedics

t aysamortho

in aysam

YouTube aysam orthopaedics



"We supply WORLD's  
Orthopedic & Trauma  
Implants..."


# Mini Fragment Locking Plates & Screws



**aysam**

# Mini Fragment Locking Plates & Screws


## 1,5 mm Mini Straight Compression Locking Plate



| TI               | HOLES | LENGTH (mm) |
|------------------|-------|-------------|
| A310 02 053 2704 | 4     | 23          |
| A310 02 053 2706 | 6     | 36          |

**Intended use :** Temporary fixation of fractures, osteotomies, nonunions, malunions, pseudoarthroses, replantations, and fusions of small bones, metatarsal/ metacarpal bones and small bone fragments.


## 1,5 mm Mini Straight Reconstruction Locking Plate



| TI               | HOLES | LENGTH (mm) |
|------------------|-------|-------------|
| A310 02 053 2705 | 6     | 29          |
| A310 02 053 2712 | 12    | 59          |

**Intended use :** Temporary fixation of fractures, osteotomies, nonunions, malunions, pseudoarthroses, replantations, and fusions of small bones, metatarsal/ metacarpal bones and small bone fragments.


## 1,5 mm Mini T Locking Plate



| TI               | HOLES (SHAFT) | HOLES (HEAD) | LENGTH (mm) |
|------------------|---------------|--------------|-------------|
| A310 02 053 2709 | 8             | 3            | 44,5        |
| A310 02 053 2710 | 8             | 4            | 44,5        |

**Intended use :** Temporary fixation of fractures, osteotomies, nonunions, malunions, pseudoarthroses, replantations, and fusions of small bones, metatarsal/ metacarpal bones and small bone fragments.


## 1,5 mm Mini Strut Locking Plate



| TI               | SIDE  | LENGTH (mm) |
|------------------|-------|-------------|
| A310 02 053 2707 | RIGHT | 23          |
| A310 02 053 2708 | LEFT  | 23          |

**Intended use :** Temporary fixation of fractures, osteotomies, nonunions, malunions, pseudoarthroses, replantations, and fusions of small bones, metatarsal/ metacarpal bones and small bone fragments.

## 1,5 mm Mini Y Locking Plate (3 holes head)




| TI               | HOLES (SHAFT) | HOLES (HEAD) | LENGTH (mm) |
|------------------|---------------|--------------|-------------|
| A310 02 053 2711 | 8             | 3            | 47          |

**Intended use :** Temporary fixation of fractures, osteotomies, nonunions, malunions, pseudoarthroses, replantations, and fusions of small bones, metatarsal/ metacarpal bones and small bone fragments.

## Mini Fragment Locking Plates & Screws

### 1,5 mm Mini Condylar Locking Plate



| TI               | HOLES (SHAFT) | HOLES (HEAD) | LENGTH (mm) |
|------------------|---------------|--------------|-------------|
| A310 02 053 0206 | 6             | 2            | 36          |

**Intended use :** Temporary fixation of fractures, osteotomies, nonunions, malunions, pseudoarthroses, replantations, and fusions of small bones, metatarsal/ metacarpal bones and small bone fragments.

### 1,5 mm Mini Locking Screw



| TI       | DIA. x LENGTH (mm) | TI       | DIA. x LENGTH (mm) |
|----------|--------------------|----------|--------------------|
| A1982606 | 1,5X6              | A1982632 | 1,5X13             |
| A1982626 | 1,5X7              | A1982614 | 1,5X14             |
| A1982608 | 1,5X8              | A1982634 | 1,5X15             |
| A1982628 | 1,5X9              | A1982616 | 1,5X16             |
| A1982610 | 1,5X10             | A1982618 | 1,5X18             |
| A1982630 | 1,5X11             | A1982620 | 1,5X20             |
| A1982612 | 1,5X12             | A1982622 | 1,5X22             |
|          |                    | A1982624 | 1,5X24             |


### 1,5 mm Mini Screw



| TI               | DIA. x LENGTH (mm) |
|------------------|--------------------|
| A110 03 215 0060 | 1,5X6              |
| A110 03 215 0070 | 1,5X7              |
| A110 03 215 0080 | 1,5X8              |
| A110 03 215 0090 | 1,5X9              |
| A110 03 215 0100 | 1,5X10             |
| A110 03 215 0120 | 1,5X12             |
| A110 03 215 0140 | 1,5X14             |
| A110 03 215 0160 | 1,5X16             |
| A110 03 215 0180 | 1,5X18             |
| A110 03 215 0200 | 1,5X20             |

## Mini Fragment Locking Plates & Screws


### 2,0 mm Mini Straight Locking Plate



| TI               | HOLES | LENGTH (mm) |
|------------------|-------|-------------|
| A110 02 182 0004 | 4     | 31          |
| A110 02 182 0005 | 5     | 38          |
| A110 02 182 0006 | 6     | 45          |
| A110 02 182 0007 | 7     | 52          |
| A110 02 182 0008 | 8     | 59          |
| A110 02 182 0010 | 10    | 73          |

**Intended use :** Temporary fixation of fractures, osteotomies, nonunions, malunions, pseudoarthroses, replantations, and fusions of small bones, metatarsal/ metacarpal bones and small bone fragments.


### 2,0 mm Mini Condylar Locking Plate (2 holes head)



| TI               | HOLES (SHAFT) | HOLES (HEAD) | LENGTH (mm) |
|------------------|---------------|--------------|-------------|
| A310 02 042 0003 | 7             | 2            | 54          |

**Intended use :** Temporary fixation of fractures, osteotomies, nonunions, malunions, pseudoarthroses, replantations, and fusions of small bones, metatarsal/ metacarpal bones and small bone fragments.

### 2,0 mm Mini Rotation Correction Locking Plate (2 holes head)




| TI               | HOLES (SHAFT) | HOLES (HEAD) | LENGTH (mm) |
|------------------|---------------|--------------|-------------|
| A310 02 051 0012 | 4             | 2            | 34          |
| A310 02 051 0013 | 5             | 2            | 40          |

**Intended use :** Temporary fixation of fractures, osteotomies, nonunions, malunions, pseudoarthroses, replantations, and fusions of small bones, metatarsal/ metacarpal bones and small bone fragments.

## Mini Fragment Locking Plates & Screws


### 2,0 mm Mini T Locking Plate (3 holes head)



| TI               | HOLES (SHAFT) | HOLES (HEAD) | LENGTH (mm) |
|------------------|---------------|--------------|-------------|
| A310 02 042 0012 | 7             | 3            | 52,5        |

**Intended use :** Temporary fixation of fractures, osteotomies, nonunions, malunions, pseudoarthroses, replantations, and fusions of small bones, metatarsal/ metacarpal bones and small bone fragments.


### 2,0 mm Mini Adaptation Locking Plate



| TI               | HOLES | LENGTH (mm) |
|------------------|-------|-------------|
| A310 02 063 0003 | 12    | 81,26       |

**Intended use :** Temporary fixation of fractures, osteotomies, nonunions, malunions, pseudoarthroses, replantations, and fusions of small bones, metatarsal/ metacarpal bones and small bone fragments.

### 2,0 mm Mini Adaptation T Locking Plate (2 holes head)




| TI               | HOLES (SHAFT) | HOLES (HEAD) | LENGTH (mm) |
|------------------|---------------|--------------|-------------|
| A310 02 041 0103 | 7             | 2            | 52,5        |

**Intended use :** Temporary fixation of fractures, osteotomies, nonunions, malunions, pseudoarthroses, replantations, and fusions of small bones, metatarsal/ metacarpal bones and small bone fragments.

## Mini Fragment Locking Plates & Screws


### 2,0 mm Mini Adaptation Y Locking Plate (3 holes head)



| TI               | HOLES (SHAFT) | HOLES (HEAD) | LENGTH (mm) |
|------------------|---------------|--------------|-------------|
| A110 02 182 0009 | 7             | 3            | 54          |

**Intended use :** Temporary fixation of fractures, osteotomies, nonunions, malunions, pseudoarthroses, replantations, and fusions of small bones, metatarsal/ metacarpal bones and small bone fragments.

### 2,0 mm Mini H Locking Plate



| TI               | HOLES | LENGTH (mm) |
|------------------|-------|-------------|
| A110 02 182 0003 | 4     | 12          |

**Intended use :** Temporary fixation of fractures, osteotomies, nonunions, malunions, pseudoarthroses, replantations, and fusions of small bones, metatarsal/ metacarpal bones and small bone fragments.

### 2,0 mm Mini Locking Screw



| TI       | DIA. x LENGTH (mm) | TI       | DIA. x LENGTH (mm) |
|----------|--------------------|----------|--------------------|
| A1982005 | 2,0X6              | A1982015 | 2,0X15             |
| A1982007 | 2,0X7              | A1982025 | 2,0X16             |
| A1982017 | 2,0X8              | A1982027 | 2,0X18             |
| A1982009 | 2,0X9              | A1982636 | 2,0X20             |
| A1982019 | 2,0X10             | A1982638 | 2,0X22             |
| A1982011 | 2,0X11             | A1982640 | 2,0X24             |
| A1982021 | 2,0X12             | A1983565 | 2,0X26             |
| A1982013 | 2,0X13             | A1983570 | 2,0X28             |
| A1982023 | 2,0X14             |          |                    |


### 2,0 mm Mini Screw



| TI               | DIA. x LENGTH (mm) | TI               | DIA. x LENGTH (mm) |
|------------------|--------------------|------------------|--------------------|
| A110 03 216 0060 | 2,0X6              | A110 03 216 0160 | 2,0X16             |
| A110 03 216 0320 | 2,0X7              | A110 03 216 0180 | 2,0X18             |
| A110 03 216 0080 | 2,0X8              | A110 03 216 0200 | 2,0X20             |
| A110 03 216 0340 | 2,0X9              | A110 03 216 0220 | 2,0X22             |
| A110 03 216 0100 | 2,0X10             | A110 03 216 0240 | 2,0X24             |
| A110 03 216 0360 | 2,0X11             | A110 03 216 0260 | 2,0X26             |
| A110 03 216 0120 | 2,0X12             | A110 03 216 0280 | 2,0X28             |
| A110 03 216 0380 | 2,0X13             | A110 03 216 0300 | 2,0X30             |
| A110 03 216 0140 | 2,0X14             |                  |                    |

## Mini Fragment Locking Plates & Screws


### 2,4 mm Mini Straight Locking Plate



| TI               | HOLES | LENGTH (mm) |
|------------------|-------|-------------|
| A110 02 222 0006 | 4     | 36,0        |
| A110 02 222 0008 | 5     | 43,6        |
| A110 02 222 0010 | 6     | 51,2        |
| A110 02 222 0012 | 7     | 58,8        |
| A110 02 222 0014 | 8     | 66,4        |
| A110 02 222 0016 | 10    | 81,6        |

**Intended use :** Temporary fixation of fractures, osteotomies, nonunions, replantations, and fusions of metatarsal/ metacarpal bones and small bone fragments.


### 2,4 mm Mini T Locking Plate (3 holes head)



| TI               | HOLES (SHAFT) | HOLES (HEAD) | LENGTH (mm) |
|------------------|---------------|--------------|-------------|
| A110 02 222 0018 | 7             | 3            | 57          |

**Intended use :** Temporary fixation of fractures, osteotomies, nonunions, replantations, and fusions of metatarsal/ metacarpal bones and small bone fragments.

### 2,4 mm Mini Condylar Locking Plate




| TI               | HOLES (SHAFT) | HOLES (HEAD) | LENGTH (mm) |
|------------------|---------------|--------------|-------------|
| A110 02 222 0020 | 7             | 2            | 59          |


**Intended use :** Fixation of fractures, osteotomies, nonunions, replantations, and fusions of metatarsal/ metacarpal bones and small bone fragments.

# Mini Fragment Locking Plates & Screws

## 2,4 mm Mini Adaptation Locking Plate




| TI               | HOLES | LENGTH (mm) |
|------------------|-------|-------------|
| A110 02 222 0022 | 12    | 87,9        |




**Intended use :** Temporary fixation of fractures, osteotomies, nonunions, replantations, and fusions of metatarsal/metacarpal bones and small bone fragments.

## 2,4 mm Mini Adaptation T Locking Plate (2 holes head)




| TI               | HOLES (SHAFT) | HOLES (HEAD) | LENGTH (mm) |
|------------------|---------------|--------------|-------------|
| A110 02 222 0024 | 7             | 2            | 59          |




**Intended use :** Temporary fixation of fractures, osteotomies, nonunions, replantations, and fusions of metatarsal/metacarpal bones and small bone fragments.

## 2,4 mm Mini Adaptation Y Locking Plate (3 holes head)



| TI               | HOLES (SHAFT) | HOLES (HEAD) | LENGTH (mm) |
|------------------|---------------|--------------|-------------|
| A310 02 043 0013 | 7             | 3            | 60,5        |



**Intended use :** Temporary fixation of fractures, osteotomies, nonunions, replantations, and fusions of metatarsal/metacarpal bones and small bone fragments.

## Mini Fragment Locking Plates & Screws

### 2,4 mm Mini Locking Screw



| TI               | DIA. x LENGTH (mm) |
|------------------|--------------------|
| A110 03 213 0060 | 2,4X6              |
| A100 03216 0070  | 2,4X7              |
| A110 03 213 0080 | 2,4X8              |
| A100 03216 0090  | 2,4X9              |
| A110 03 213 0100 | 2,4X10             |
| A100 03216 0011  | 2,4X11             |
| A110 03 213 0120 | 2,4X12             |
| A100 03216 0013  | 2,4X13             |
| A110 03 213 0140 | 2,4X14             |

| TI               | DIA. x LENGTH (mm) |
|------------------|--------------------|
| A100 03216 0015  | 2,4X15             |
| A110 03 213 0160 | 2,4X16             |
| A110 03 213 0180 | 2,4X18             |
| A110 03 213 0200 | 2,4X20             |
| A110 03 213 0220 | 2,4X22             |
| A110 03 213 0240 | 2,4X24             |
| A110 03 213 0260 | 2,4X26             |
| A110 03 213 0280 | 2,4X28             |
| A110 03 213 0300 | 2,4X30             |


### 2,4 mm Mini Screw



| TI       | DIA. x LENGTH (mm) |
|----------|--------------------|
| A1982306 | 2,4X6              |
| A1982308 | 2,4X8              |
| A1982310 | 2,4X10             |
| A1982312 | 2,4X12             |
| A1982314 | 2,4X14             |
| A1982316 | 2,4X16             |
| A1982318 | 2,4X18             |
| A1982320 | 2,4X20             |
| A1982322 | 2,4X22             |

| TI              | DIA. x LENGTH (mm) |
|-----------------|--------------------|
| A1982324        | 2,4X24             |
| A1982326        | 2,4X26             |
| A1982328        | 2,4X28             |
| A1982330        | 2,4X30             |
| A100 03216 0017 | 2,4X7              |
| A100 03216 0019 | 2,4X9              |
| A100 03216 0021 | 2,4X11             |
| A100 03216 0023 | 2,4X13             |
| A100 03216 0025 | 2,4X15             |

### 2,7 mm Mini Condylar Locking Plate




| TI               | HOLES (SHAFT) | HOLES (HEAD) | LENGTH (mm) |
|------------------|---------------|--------------|-------------|
| A310 02 063 0002 | 7             | 2            | 66          |


**Intended use :** Temporary fixation of fractures, osteotomies, nonunions, replantations, and fusions of metatarsal/ metacarpal bones and small bone fragments.

## Mini Fragment Locking Plates & Screws

### 2,7 mm Mini Adaptation Locking Plate




| TI               | HOLES | LENGTH (mm) |
|------------------|-------|-------------|
| A110 02 410 0012 | 12    | 97          |



**Intended use :** Temporary fixation of fractures, osteotomies, nonunions, replantations, and fusions of metatarsal/metacarpal bones and small bone fragments.



### 2,7 mm Mini T Locking Plate (2 holes head)



| TI               | HOLES (SHAFT) | HOLES (HEAD) | LENGTH (mm) |
|------------------|---------------|--------------|-------------|
| A310 02 062 0002 | 3             | 2            | 32          |
| A310 02 062 0003 | 4             | 2            | 40          |

**Intended use :** Temporary fixation of fractures, osteotomies, nonunions, replantations, and fusions of metatarsal/metacarpal bones and small bone fragments.

### 2,7 mm Mini Straight Locking Plate




| TI               | HOLES | LENGTH (mm) |
|------------------|-------|-------------|
| A110 02 410 0004 | 4     | 40          |
| A110 02 410 0005 | 5     | 49          |
| A110 02 410 0006 | 6     | 58          |
| A110 02 410 0007 | 7     | 67          |
| A110 02 410 0008 | 8     | 76          |
| A110 02 410 0010 | 10    | 94          |

**Intended use :** Temporary fixation of fractures, osteotomies, nonunions, replantations, and fusions of metatarsal/metacarpal bones and small bone fragments.

# Mini Fragment Locking Plates & Screws


## 2,7 mm Mini L Locking Plate



| TI               | HOLES (SHAFT) | HOLES (HEAD) | LENGTH (mm) | SIDE  |
|------------------|---------------|--------------|-------------|-------|
| A310 02 040 0003 | 3             | 2            | 32          | Left  |
| A310 02 041 0003 | 3             | 2            | 32          | Right |
| A310 02 050 0003 | 4             | 2            | 40          | Left  |
| A310 02 050 0004 | 4             | 2            | 40          | Right |

**Intended use :** Temporary fixation of fractures, osteotomies, nonunions, replantations, and fusions of metatarsal/metacarpal bones and small bone fragments.

## 2,7 mm Mini Oblique L Locking Plate



| TI               | HOLES (SHAFT) | HOLES (HEAD) | LENGTH (mm) | SIDE  |
|------------------|---------------|--------------|-------------|-------|
| A310 02 060 0003 | 3             | 2            | 35          | Left  |
| A310 02 061 0003 | 3             | 2            | 35          | Right |
| A110 02 410 0002 | 4             | 2            | 42          | Left  |
| A110 02 410 0009 | 4             | 2            | 42          | Right |

**Intended use :** Temporary fixation of fractures, osteotomies, nonunions, replantations, and fusions of metatarsal/metacarpal bones and small bone fragments.

## 2,7 mm Mini Locking Screw



| TI       | DIA. x LENGTH (mm) | TI       | DIA. x LENGTH (mm) |
|----------|--------------------|----------|--------------------|
| A1982710 | 2,7X10             | A1982726 | 2,7X26             |
| A1982712 | 2,7X12             | A1982728 | 2,7X28             |
| A1982714 | 2,7X14             | A1982730 | 2,7X30             |
| A1982716 | 2,7X16             | A1982732 | 2,7X32             |
| A1982718 | 2,7X18             | A1982734 | 2,7X34             |
| A1982720 | 2,7X20             | A1982736 | 2,7X36             |
| A1982722 | 2,7X22             | A1982738 | 2,7X38             |
| A1982724 | 2,7X24             | A1982740 | 2,7X40             |

## 2,7 mm Mini Screw



| TI               | DIA. x LENGTH (mm) | TI               | DIA. x LENGTH (mm) |
|------------------|--------------------|------------------|--------------------|
| A110 03 217 0100 | 2,7X10             | A110 03 217 0260 | 2,7X26             |
| A110 03 217 0120 | 2,7X12             | A110 03 217 0280 | 2,7X28             |
| A110 03 217 0140 | 2,7X14             | A110 03 217 0300 | 2,7X30             |
| A110 03 217 0160 | 2,7X16             | A110 03 217 0320 | 2,7X32             |
| A110 03 217 0180 | 2,7X18             | A110 03 217 0340 | 2,7X34             |
| A110 03 217 0200 | 2,7X20             | A110 03 217 0360 | 2,7X36             |
| A110 03 217 0220 | 2,7X22             | A110 03 217 0380 | 2,7X38             |
| A110 03 217 0240 | 2,7X24             | A110 03 217 0400 | 2,7X40             |



■ [www.aysam.com.tr](http://www.aysam.com.tr)  
■ [www.aysamstore.com.tr](http://www.aysamstore.com.tr)  
■ [aysam@aysam.com.tr](mailto:aysam@aysam.com.tr)

f [aysam.orthopedics](#)

t [aysamortho](#)

in [aysam](#)

YouTube [aysam orthopaedics](#)

"We supply WORLD's  
Orthopedic & Trauma  
Implants..."

# Chest Wall Plates, Screws & Splints



**aysam**

# Chest Wall Plates, Screws & Splints


## 2,9 mm Sternal Locking T Plate




| TI            | HOLES |
|---------------|-------|
| A15408 128 07 | 7     |



| TI            | HOLES |
|---------------|-------|
| A15408 128 08 | 8     |



**Intended use :** Temporary fixation of fractures, malunions, nonunions, pseudoarthroses and osteotomies of the sternum.



**Intended use :** Temporary fixation of fractures, malunions, nonunions, pseudoarthroses and osteotomies of the sternum.


## 2,9 mm Sternal Locking I Plate




| TI          | HOLES |
|-------------|-------|
| A15408 2906 | 9     |



| TI          | HOLES |
|-------------|-------|
| A15408 2908 | 11    |



**Intended use :** Temporary fixation of fractures, malunions, nonunions, pseudoarthroses and osteotomies of the sternum.




**Intended use :** Temporary fixation of fractures, malunions, nonunions, pseudoarthroses and osteotomies of the sternum.

## Chest Wall Plates, Screws & Splints

### 2,9 mm Sternal Locking Straight Plate



| TI          | HOLES |
|-------------|-------|
| A15408 2910 | 8     |
| A15408 2912 | 9     |
| A15408 2914 | 10    |




**Intended use :** Fixation of fractures, malunions, nonunions, pseudoarthroses and osteotomies of the sternum.

### 2,9 mm Universal Locking Straight Plate



| TI          | HOLES |
|-------------|-------|
| A15408 2916 | 8     |
| A15408 2918 | 24    |
| A15408 2920 | 30    |




**Intended use :** Fixation of fractures, malunions, nonunions, pseudoarthroses and osteotomies of the ribs.

# Chest Wall Plates, Screws & Splints

## 2,9 mm Self-Tapping Rib Locking Screw




| TI          | DIA. x LENGTH<br>(mm) |
|-------------|-----------------------|
| A15409 2906 | 2,9X6                 |
| A15409 2908 | 2,9X8                 |
| A15409 2910 | 2,9X10                |
| A15409 2912 | 2,9X12                |
| A15409 2914 | 2,9X14                |
| A15409 2916 | 2,9X16                |
| A15409 2918 | 2,9X18                |
| A15409 2920 | 2,9X20                |



## 2,9 mm Self-Tapping Rib Screw



| TI          | DIA. x LENGTH<br>(mm) |
|-------------|-----------------------|
| A15508 2910 | 2,9X10                |
| A15508 2912 | 2,9X12                |



## Rib Intramedullary Splints



| TI               | WIDTH<br>(mm) | LENGTH<br>(mm) |
|------------------|---------------|----------------|
| A100 02 222 0003 | 3             | 50             |
| A100 02 222 0004 | 4             | 50             |
| A100 02 222 0005 | 5             | 50             |
| A100 02 222 0006 | 3             | 65             |
| A100 02 222 0007 | 4             | 65             |
| A100 02 222 0008 | 5             | 65             |
| A100 02 222 0009 | 3             | 80             |
| A100 02 222 0010 | 4             | 80             |
| A100 02 222 0011 | 5             | 80             |



**Intended use :** Fixation of fractures, malunions, nonunions, pseudoarthroses and osteotomies of the ribs.

"We supply WORLD's  
Orthopedic & Trauma  
Implants..."

# Angled Blade Plates


**& Pediatric  
Plates**



**aysam**

## Angled Blade Plates & Pediatric Plates


### 4,5 mm 135° Angled Blade Plate for operation on adult femur



**Intended use :** Temporary fixation of fractures, malunions, nonunions, pseudoarthroses and osteotomies of the proximal femur for adult patients

| SS               | HOLES | BLADE (mm) | LENGTH (mm) | SS               | HOLES | BLADE (mm) | LENGTH (mm) |
|------------------|-------|------------|-------------|------------------|-------|------------|-------------|
| A100 02 040 4050 | 4     | 50         | 76          | A100 02 140 1280 | 7     | 80         | 124         |
| A100 02 040 4060 | 4     | 60         | 76          | A100 02 140 1290 | 7     | 90         | 124         |
| A100 02 040 4070 | 4     | 70         | 76          | A100 02 140 1200 | 7     | 100        | 124         |
| A100 02 040 4080 | 4     | 80         | 76          | A100 02 141 9050 | 8     | 50         | 140         |
| A100 02 040 4090 | 4     | 90         | 76          | A100 02 141 9060 | 8     | 60         | 140         |
| A100 02 040 4100 | 4     | 100        | 76          | A100 02 141 9070 | 8     | 70         | 140         |
| A100 02 041 4050 | 5     | 50         | 92          | A100 02 141 9080 | 8     | 80         | 140         |
| A100 02 041 4060 | 5     | 60         | 92          | A100 02 141 9090 | 8     | 90         | 140         |
| A100 02 041 4070 | 5     | 70         | 92          | A100 02 141 9100 | 8     | 100        | 140         |
| A100 02 041 4080 | 5     | 80         | 92          | A100 02 041 9050 | 9     | 50         | 156         |
| A100 02 041 4090 | 5     | 90         | 92          | A100 02 041 9060 | 9     | 60         | 156         |
| A100 02 041 4100 | 5     | 100        | 92          | A100 02 041 9070 | 9     | 70         | 156         |
| A100 02 040 1250 | 6     | 50         | 108         | A100 02 041 9080 | 9     | 80         | 156         |
| A100 02 040 1260 | 6     | 60         | 108         | A100 02 041 9090 | 9     | 90         | 156         |
| A100 02 040 1270 | 6     | 70         | 108         | A100 02 041 9100 | 9     | 100        | 156         |
| A100 02 040 1280 | 6     | 80         | 108         | A100 02 041 1250 | 12    | 50         | 204         |
| A100 02 040 1290 | 6     | 90         | 108         | A100 02 041 1260 | 12    | 60         | 204         |
| A100 02 040 1200 | 6     | 100        | 108         | A100 02 041 1270 | 12    | 70         | 204         |
| A100 02 140 1250 | 7     | 50         | 124         | A100 02 041 1280 | 12    | 80         | 204         |
| A100 02 140 1260 | 7     | 60         | 124         | A100 02 041 1290 | 12    | 90         | 204         |
| A100 02 140 1270 | 7     | 70         | 124         | A100 02 041 1200 | 12    | 100        | 204         |

### 4,5 mm 95° Angled Blade Plate for operation on adult femur




**Intended use :** Temporary fixation of fractures, malunions, nonunions, pseudoarthroses and osteotomies of the proximal femur for adult patients

| SS               | HOLES | BLADE (mm) | LENGTH (mm) | SS               | HOLES | BLADE (mm) | LENGTH (mm) |
|------------------|-------|------------|-------------|------------------|-------|------------|-------------|
| A100 02 140 4050 | 4     | 50         | 86          | A100 02 140 6080 | 7     | 80         | 134         |
| A100 02 140 4060 | 4     | 60         | 86          | A100 02 140 7065 | 7     | 90         | 134         |
| A100 02 140 4070 | 4     | 70         | 86          | A100 02 140 7085 | 7     | 100        | 134         |
| A100 02 140 4080 | 4     | 80         | 86          | A100 02 140 8050 | 8     | 50         | 150         |
| A100 02 140 4090 | 4     | 90         | 86          | A100 02 140 8060 | 8     | 60         | 150         |
| A100 02 140 4100 | 4     | 100        | 86          | A100 02 140 8070 | 8     | 70         | 150         |
| A100 02 141 4050 | 5     | 50         | 102         | A100 02 140 8080 | 8     | 80         | 150         |
| A100 02 141 4060 | 5     | 60         | 102         | A100 02 141 6090 | 8     | 90         | 150         |
| A100 02 141 4070 | 5     | 70         | 102         | A100 02 141 6100 | 8     | 100        | 150         |
| A100 02 141 4080 | 5     | 80         | 102         | A100 02 140 9050 | 9     | 50         | 166         |
| A100 02 141 4090 | 5     | 90         | 102         | A100 02 140 9060 | 9     | 60         | 166         |
| A100 02 141 4100 | 5     | 100        | 102         | A100 02 140 9070 | 9     | 70         | 166         |
| A100 02 140 6051 | 6     | 50         | 118         | A100 02 140 9080 | 9     | 80         | 166         |
| A100 02 140 6061 | 6     | 60         | 118         | A100 02 140 9090 | 9     | 90         | 166         |
| A100 02 140 6071 | 6     | 70         | 118         | A100 02 140 9100 | 9     | 100        | 166         |
| A100 02 140 6081 | 6     | 80         | 118         | A100 02 141 1250 | 12    | 50         | 214         |
| A100 02 140 6090 | 6     | 90         | 118         | A100 02 141 1260 | 12    | 60         | 214         |
| A100 02 140 6100 | 6     | 100        | 118         | A100 02 141 1270 | 12    | 70         | 214         |
| A100 02 140 6050 | 7     | 50         | 134         | A100 02 141 1280 | 12    | 80         | 214         |
| A100 02 140 6060 | 7     | 60         | 134         | A100 02 141 1290 | 12    | 90         | 214         |
| A100 02 140 6070 | 7     | 70         | 134         | A100 02 141 1200 | 12    | 100        | 214         |

## Angled Blade Plates & Pediatric Plates


### 4,5 mm 90° Angled Blade Plate Adolescent



| SS               | HOLES | BLADE (mm) | LENGTH (mm) |
|------------------|-------|------------|-------------|
| A100 02 070 1040 | 3     | 40         | 78,5        |
| A100 02 070 1045 | 3     | 45         | 78,5        |
| A100 02 070 1050 | 3     | 50         | 78,5        |
| A100 02 070 1055 | 3     | 55         | 78,5        |
| A100 02 070 1060 | 3     | 60         | 78,5        |

**Intended use :** Temporary fixation of fractures, malunions, nonunions, pseudoarthroses and osteotomies of the proximal femur in children. Especially used for fixation of subtrochanteric varus osteotomies in treatment of congenital dislocation of hip in adolescents.

### 3,5 mm 90° Angled Blade Plate Pediatric




| SS               | HOLES | BLADE (mm) | LENGTH (mm) |
|------------------|-------|------------|-------------|
| A100 02 050 1225 | 3     | 25         | 51,7        |
| A100 02 050 1230 | 3     | 30         | 51,7        |
| A100 02 050 1232 | 3     | 32         | 51,7        |
| A100 02 050 1235 | 3     | 35         | 51,7        |

**Intended use :** Temporary fixation of fractures, malunions, nonunions, pseudoarthroses and osteotomies of the proximal femur in children. Especially used for fixation of subtrochanteric varus osteotomies in treatment of congenital dislocation of hip in children.

### 3,5 mm Mini Pediatric Locking Plate 100°





| SS               | TI               | HOLES | LENGTH (mm) |
|------------------|------------------|-------|-------------|
| A100 02 051 0725 | A100 02 051 0730 | 3     | 64          |



**Intended use :** Temporary fixation of fractures, malunions, nonunions, pseudoarthroses and osteotomies of the proximal femur in children. Especially used for fixation of subtrochanteric varus osteotomies in treatment of congenital dislocation of hip in children.

## Angled Blade Plates & Pediatric Plates


### 3,5 mm Pediatric Locking Plate 100°



| SS               | TI               | HOLES | LENGTH (mm) |
|------------------|------------------|-------|-------------|
| A100 02 051 0732 | A100 02 051 0735 | 3     | 76          |

**Intended use :** Temporary fixation of fractures, malunions, nonunions, pseudoarthroses and osteotomies of the proximal femur in children. Especially used for fixation of subtrochanteric varus osteotomies in treatment of congenital dislocation of hip in children


### 3,5 mm Pediatric Locking Plate 110°



| SS               | TI               | HOLES | LENGTH (mm) |
|------------------|------------------|-------|-------------|
| A100 02 051 1232 | A100 02 051 1235 | 3     | 76          |

**Intended use :** Temporary fixation of fractures, malunions, nonunions, pseudoarthroses and osteotomies of the proximal femur in children. Especially used for fixation of subtrochanteric varus osteotomies in treatment of congenital dislocation of hip in children.


### 3,5 mm Pediatric Locking Plate 150°



| SS               | TI               | HOLES | LENGTH (mm) |
|------------------|------------------|-------|-------------|
| A100 02 050 0725 | A100 02 050 0730 | 3     | 65          |

**Intended use :** Temporary fixation of fractures, malunions, nonunions, pseudoarthroses and osteotomies of the proximal femur in children. Especially used for fixation of subtrochanteric valgus osteotomies in treatment of congenital dislocation of hip in children.

### Pediatric Eight Plate (for epiphysiodesis)



| TI               | TYPE  |
|------------------|-------|
| A100 02 050 0732 | Small |
| A100 02 050 0735 | Large |

**Intended use :** Gradual correction of pediatric deformities in both the upper and lower extremities.

"We supply WORLD's  
Orthopedic & Trauma  
Implants..."

# Trochanteric Plates & Cable Systems




# Trochanteric Plates & Cable Systems

## 4,5/5,0 mm Straight Cable Grip Plate



| SS               | TYPE   | HOLES | LENGTH (mm) |
|------------------|--------|-------|-------------|
| A900 04 020 0051 | Small  | 8     | 174         |
| A900 04 020 0052 | Medium | 10    | 218         |
| A900 04 020 0053 | Large  | 12    | 262         |




**Intended use :** Temporary fixation of fractures, malunions, nonunions, pseudoarthroses and osteotomies of the proximal femur and femoral shaft, especially in periprosthetic fractures. This plate group is used together with screws, cables and crimp blocks.

## 4,5/5,0 mm Trochanteric Hook Cable Plate 1



| TI               | TYPE   | HOLES | LENGTH (mm) |
|------------------|--------|-------|-------------|
| A900 04 021 0051 | Small  | 3     | 137,5       |
| A900 04 021 0052 | Medium | 5     | 182,5       |
| A900 04 021 0053 | Large  | 7     | 227,5       |




**Intended use :** Temporary fixation of fractures, malunions, nonunions, pseudoarthroses and osteotomies of the proximal femur and femoral shaft, especially in periprosthetic fractures. This plate group is used together with screws, cables and crimp blocks.

# Trochanteric Plates & Cable Systems

## 4,5/5,0 mm Trochanteric Cable Hook Plate 2



| TI               | TYPE   | HOLES | LENGTH (mm) |
|------------------|--------|-------|-------------|
| A900 04 023 0051 | Small  | 1     | 97          |
| A900 04 023 0052 | Medium | 3     | 143         |
| A900 04 023 0053 | Large  | 6     | 189         |




**Intended use :** Temporary fixation of fractures, malunions, nonunions, pseudoarthroses and osteotomies of the proximal femur and femoral shaft, especially in periprosthetic fractures. This plate group is used together with screws, cables and crimp blocks.

## Trochanteric Cable Hook Grip Plate



| SS               | TYPE   | BLOCK | LENGTH (mm) |
|------------------|--------|-------|-------------|
| A900 04 022 0051 | Small  | 1     | 95          |
| A900 04 022 0052 | Medium | 3     | 141         |
| A900 04 022 0053 | Large  | 5     | 187         |



**Intended use :** Temporary fixation of fractures, malunions, nonunions, pseudoarthroses and osteotomies of the proximal femur. Especially used for fixation of periprosthetic fractures. This plate group is used together with cables and crimp blocks.

# Trochanteric Plates & Cable Systems

## Trochanteric Grip



SS/TI

A900 04 020 0050



**Intended use :** Temporary fixation of fractures, malunions, nonunions, pseudoarthroses and osteotomies of the proximal femur. Especially used for fixation of periprosthetic fractures. This plate group is used together with cables and crimp blocks.

## 3,5 mm Cable Locking Screw



TI

A900 03 020 0005



**Intended use :** Fixation of cables used together with 3,5 mm locking plates.

## 5,0 mm Cable Locking Screw




TI

A900 03 020 0008



**Intended use :** Fixation of cables used together with 5,0 mm locking plates.

## Cable Crimp Block



**Intended use :** Temporary fixation of fractures, malunions, nonunions, pseudoarthroses and osteotomies of the trochanteric region of proximal femur. Especially used for fixation of periprosthetic fractures. This grip is used together with cables and crimp blocks.

| SS               | TYPE               |
|------------------|--------------------|
| A900 03 015 0015 | Cable 1,5 mm       |
| A900 03 020 0020 | Cable 2,0 mm       |
| A900 04 015 0015 | Crimp Block 1,5 mm |
| A900 04 020 0020 | Crimp Block 2,0 mm |



"We supply WORLD's  
Orthopedic & Trauma  
Implants..."

# Cannulated Screws



# Cannulated Screws

## 2,7 mm Cannulated Screw Partially Threaded



| TI       | DIA. x LENGTH (mm) |
|----------|--------------------|
| A6072718 | 2,7X18             |
| A6072720 | 2,7X20             |
| A6072722 | 2,7X22             |
| A6072724 | 2,7X24             |
| A6072726 | 2,7X26             |
| A6072728 | 2,7X28             |
| A6072730 | 2,7X30             |



**Intended use :** Temporary fixation of fractures and arthrodeses of small bones & joints eg. ;  
 -fractures of forearm, hand and foot (distal radius, carpal, metacarpal, tarsal, metatarsal bones  
 -arthrodeses of forearm, hand and foot (intercarpal, intertarsal, interphalangeal, metacarpophalangeal, metatarsophalangeal Joints )

## 3,5 mm Cannulated Screw Partially Threaded



| TI       | DIA. x LENGTH (mm) | TI       | DIA. x LENGTH (mm) |
|----------|--------------------|----------|--------------------|
| A6073516 | 3,5X16             | A6073534 | 3,5X34             |
| A6073518 | 3,5X18             | A6073536 | 3,5X36             |
| A6073520 | 3,5X20             | A6073538 | 3,5X38             |
| A6073522 | 3,5X22             | A6073540 | 3,5X40             |
| A6073524 | 3,5X24             | A6073542 | 3,5X42             |
| A6073526 | 3,5X26             | A6073544 | 3,5X44             |
| A6073528 | 3,5X28             | A6073546 | 3,5X46             |
| A6073530 | 3,5X30             | A6073548 | 3,5X48             |
| A6073532 | 3,5X32             | A6073550 | 3,5X50             |



**Intended use :** Temporary fixation of fractures and arthrodeses of small bones & joints eg. ;  
 -fractures of forearm, hand and foot ( distal radius, carpal, metacarpal, tarsal, metatarsal bones  
 -arthrodeses of forearm, hand and foot (intercarpal, intertarsal, interphalangeal, metacarpophalangeal, metatarsophalangeal Joints )

## 3,5 mm Cannulated Screw Full Threaded



| TI       | DIA. x LENGTH (mm) | TI       | DIA. x LENGTH (mm) |
|----------|--------------------|----------|--------------------|
| A6083514 | 3,5X14             | A6083534 | 3,5X34             |
| A6083516 | 3,5X16             | A6083536 | 3,5X36             |
| A6083518 | 3,5X18             | A6083538 | 3,5X38             |
| A6083520 | 3,5X20             | A6083540 | 3,5X40             |
| A6083522 | 3,5X22             | A6083542 | 3,5X42             |
| A6083524 | 3,5X24             | A6083544 | 3,5X44             |
| A6083526 | 3,5X26             | A6083546 | 3,5X46             |
| A6083528 | 3,5X28             | A6083548 | 3,5X48             |
| A6083530 | 3,5X30             | A6083550 | 3,5X50             |
| A6083532 | 3,5X32             |          |                    |



**Intended use :** Temporary fixation of fractures and arthrodeses of small bones & joints eg. ;  
 -fractures of forearm, hand and foot ( distal radius, carpal, metacarpal, tarsal, metatarsal bones  
 -arthrodeses of forearm, hand and foot (intercarpal, intertarsal, interphalangeal, metacarpophalangeal, metatarsophalangeal Joints )

# Cannulated Screws

## 4,0 mm Cannulated Screw Partially Threaded



| TI       | DIA. x LENGTH (mm) |
|----------|--------------------|
| A6074020 | 4,0X20             |
| A6074022 | 4,0X22             |
| A6074024 | 4,0X24             |
| A6074026 | 4,0X26             |
| A6074028 | 4,0X28             |
| A6074030 | 4,0X30             |
| A6074032 | 4,0X32             |
| A6074034 | 4,0X34             |
| A6074036 | 4,0X36             |
| A6074038 | 4,0X38             |
| A6074040 | 4,0X40             |
| A6074042 | 4,0X42             |
| A6074044 | 4,0X44             |

| TI       | DIA. x LENGTH (mm) |
|----------|--------------------|
| A6074046 | 4,0X46             |
| A6074048 | 4,0X48             |
| A6074050 | 4,0X50             |
| A6074052 | 4,0X52             |
| A6074054 | 4,0X54             |
| A6074056 | 4,0X56             |
| A6074058 | 4,0X58             |
| A6074060 | 4,0X60             |
| A6074062 | 4,0X62             |
| A6074064 | 4,0X64             |
| A6074066 | 4,0X66             |
| A6074068 | 4,0X68             |
| A6074070 | 4,0X70             |



### Intended use :

- Temporary fixation of fractures with medium fragments, e.g.:
- tarsal and metatarsal fractures and fixation in tarsal, metatarsal osteotomies,
  - tarsometatarsal and metatarsophalangeal arthrodeses,
  - ligament fixation,
  - hallux valgus correction osteotomies

## 4,0 mm Cannulated Screw Full Threaded



| TI       | DIA. x LENGTH (mm) |
|----------|--------------------|
| A6084020 | 4,0X20             |
| A6084022 | 4,0X22             |
| A6084024 | 4,0X24             |
| A6084026 | 4,0X26             |
| A6084028 | 4,0X28             |
| A6084030 | 4,0X30             |
| A6084032 | 4,0X32             |
| A6084034 | 4,0X34             |
| A6084036 | 4,0X36             |
| A6084038 | 4,0X38             |
| A6084040 | 4,0X40             |
| A6084042 | 4,0X42             |
| A6084044 | 4,0X44             |

| TI       | DIA. x LENGTH (mm) |
|----------|--------------------|
| A6084046 | 4,0X46             |
| A6084048 | 4,0X48             |
| A6084050 | 4,0X50             |
| A6084052 | 4,0X52             |
| A6084054 | 4,0X54             |
| A6084056 | 4,0X56             |
| A6084058 | 4,0X58             |
| A6084060 | 4,0X60             |
| A6084062 | 4,0X62             |
| A6084064 | 4,0X64             |
| A6084066 | 4,0X66             |
| A6084068 | 4,0X68             |
| A6084070 | 4,0X70             |

### Intended use :

- Temporary fixation of fractures with medium fragments, e.g.:
- tarsal and metatarsal fractures and fixation in tarsal, metatarsal osteotomies,
  - tarsometatarsal and metatarsophalangeal arthrodeses,
  - ligament fixation,
  - hallux valgus correction osteotomies



# Cannulated Screws

## 4,5 mm Cannulated Screw Partially Threaded



| SS      | TI       | DIA. x LENGTH (mm) |
|---------|----------|--------------------|
| A030120 | A6074220 | 4,5X20             |
| A030122 | A6074222 | 4,5X22             |
| A030124 | A6074224 | 4,5X24             |
| A030126 | A6074226 | 4,5X26             |
| A030128 | A6074228 | 4,5X28             |
| A030130 | A6074230 | 4,5X30             |
| A030132 | A6074232 | 4,5X32             |
| A030134 | A6074234 | 4,5X34             |
| A030136 | A6074236 | 4,5X36             |
| A030138 | A6074238 | 4,5X38             |

| SS      | TI       | DIA. x LENGTH (mm) |
|---------|----------|--------------------|
| A030140 | A6074240 | 4,5X40             |
| A030142 | A6074242 | 4,5X42             |
| A030144 | A6074244 | 4,5X44             |
| A030146 | A6074246 | 4,5X46             |
| A030148 | A6074248 | 4,5X48             |
| A030150 | A6074250 | 4,5X50             |
| A030152 | A6074252 | 4,5X52             |
| A030154 | A6074254 | 4,5X54             |
| A030156 | A6074256 | 4,5X56             |
| A030158 | A6074258 | 4,5X58             |
| A030160 | A6074260 | 4,5X60             |



**Intended use :** Temporary fixation of fractures with medium fragments, e.g.:

- malleolar fractures
- pilon tibial fractures
- fractures of the calcaneus and talus
- tibial plateau fractures
- carpal and tarsal arthroses

## 4,5 mm Cannulated Screw Full Threaded



| SS      | TI       | DIA. x LENGTH (mm) |
|---------|----------|--------------------|
| A040120 | A6074520 | 4,5X20             |
| A040122 | A6074522 | 4,5X22             |
| A040124 | A6074524 | 4,5X24             |
| A040126 | A6074526 | 4,5X26             |
| A040128 | A6074528 | 4,5X28             |
| A040130 | A6074530 | 4,5X30             |
| A040132 | A6074532 | 4,5X32             |
| A040134 | A6074534 | 4,5X34             |
| A040136 | A6074536 | 4,5X36             |
| A040138 | A6074538 | 4,5X38             |
| A040140 | A6074540 | 4,5X40             |

| SS      | TI       | DIA. x LENGTH (mm) |
|---------|----------|--------------------|
| A040142 | A6074542 | 4,5X42             |
| A040144 | A6074544 | 4,5X44             |
| A040146 | A6074546 | 4,5X46             |
| A040148 | A6074548 | 4,5X48             |
| A040150 | A6074550 | 4,5X50             |
| A040152 | A6074552 | 4,5X52             |
| A040154 | A6074554 | 4,5X54             |
| A040156 | A6074556 | 4,5X56             |
| A040158 | A6074558 | 4,5X58             |
| A040160 | A6074560 | 4,5X60             |




**Intended use :** Temporary fixation of fractures with medium fragments, e.g.:

- malleolar fractures
- pilon tibial fractures
- fractures of the calcaneus and talus
- tibial plateau fractures
- carpal and tarsal arthroses

## Cannulated Screws

### 6,5 mm Cannulated Screw 16 mm Threaded




**Intended use :** Temporary fixation of fractures with large fragments, e.g.:

- femoral neck fractures,
- intercondylar femoral fractures,
- epiphysiylolysis of the femoral head,
- ankle arthrodeses,
- iliosacral dislocations,
- proximal tibia plateau fractures



| SS               | TI               | DIA. x LENGTH (mm) |
|------------------|------------------|--------------------|
| A101 04 225 0300 | A111 04 225 0300 | 6,5X30             |
| A101 04 225 0350 | A111 04 225 0350 | 6,5X35             |
| A101 04 225 0400 | A111 04 225 0400 | 6,5X40             |
| A101 04 225 0450 | A111 04 225 0450 | 6,5X45             |
| A101 04 225 0500 | A111 04 225 0500 | 6,5X50             |
| A101 04 225 0550 | A111 04 225 0550 | 6,5X55             |
| A101 04 225 0600 | A111 04 225 0600 | 6,5X60             |
| A101 04 225 0650 | A111 04 225 0650 | 6,5X65             |
| A101 04 225 0700 | A111 04 225 0700 | 6,5X70             |
| A101 04 225 0750 | A111 04 225 0750 | 6,5X75             |
| A101 04 225 0800 | A111 04 225 0800 | 6,5X80             |
| A101 04 225 0850 | A111 04 225 0850 | 6,5X85             |
| A101 04 225 0900 | A111 04 225 0900 | 6,5X90             |
| A101 04 225 0950 | A111 04 225 0950 | 6,5X95             |
| A101 04 225 1000 | A111 04 225 1000 | 6,5X100            |
| A101 04 225 1050 | A111 04 225 1050 | 6,5X105            |
| A101 04 225 1100 | A111 04 225 1100 | 6,5X110            |
| A101 04 225 1150 | A111 04 225 1150 | 6,5X115            |
| A101 04 225 1200 | A111 04 225 1200 | 6,5X120            |

### 6,5 mm Cannulated Screw 32 mm Threaded



**Intended use :** Temporary fixation of fractures with large fragments, e.g.:

- femoral neck fractures,
- intercondylar femoral fractures,
- epiphysiylolysis of the femoral head,
- ankle arthrodeses,
- iliosacral dislocations,
- proximal tibia plateau fractures



| SS               | TI               | DIA. x LENGTH (mm) |
|------------------|------------------|--------------------|
| A101 04 300 0500 | A111 04 300 0500 | 6,5X50             |
| A101 04 300 0550 | A111 04 300 0550 | 6,5X55             |
| A101 04 300 0600 | A111 04 300 0600 | 6,5X60             |
| A101 04 300 0650 | A111 04 300 0650 | 6,5X65             |
| A101 04 300 0700 | A111 04 300 0700 | 6,5X70             |
| A101 04 300 0750 | A111 04 300 0750 | 6,5X75             |
| A101 04 300 0800 | A111 04 300 0800 | 6,5X80             |
| A101 04 300 0850 | A111 04 300 0850 | 6,5X85             |
| A101 04 300 0900 | A111 04 300 0900 | 6,5X90             |
| A101 04 300 0950 | A111 04 300 0950 | 6,5X95             |
| A101 04 300 1000 | A111 04 300 1000 | 6,5X100            |
| A101 04 300 1050 | A111 04 300 1050 | 6,5X105            |
| A101 04 300 1100 | A111 04 300 1100 | 6,5X110            |
| A101 04 300 1150 | A111 04 300 1150 | 6,5X115            |
| A101 04 300 1200 | A111 04 300 1200 | 6,5X120            |

# Cannulated Screws

## 6,5 mm Cannulated Screw Full Threaded



| SS               | TI               | DIA. x LENGTH (mm) | SS               | TI               | DIA. x LENGTH (mm) |
|------------------|------------------|--------------------|------------------|------------------|--------------------|
| A101 04 330 0350 | A111 04 330 0350 | 6,5X35             | A101 04 330 0800 | A111 04 330 0800 | 6,5X80             |
| A101 04 330 0400 | A111 04 330 0400 | 6,5X40             | A101 04 330 0850 | A111 04 330 0850 | 6,5X85             |
| A101 04 330 0450 | A111 04 330 0450 | 6,5X45             | A101 04 330 0900 | A111 04 330 0900 | 6,5X90             |
| A101 04 330 0500 | A111 04 330 0500 | 6,5X50             | A101 04 330 0950 | A111 04 330 0950 | 6,5X95             |
| A101 04 330 0550 | A111 04 330 0550 | 6,5X55             | A101 04 330 1000 | A111 04 330 1000 | 6,5X100            |
| A101 04 330 0600 | A111 04 330 0600 | 6,5X60             | A101 04 330 1050 | A111 04 330 1050 | 6,5X105            |
| A101 04 330 0650 | A111 04 330 0650 | 6,5X65             | A101 04 330 1100 | A111 04 330 1100 | 6,5X110            |
| A101 04 330 0700 | A111 04 330 0700 | 6,5X70             | A101 04 330 1150 | A111 04 330 1150 | 6,5X115            |
| A101 04 330 0750 | A111 04 330 0750 | 6,5X75             | A101 04 330 1200 | A111 04 330 1200 | 6,5X120            |



**Intended use :** Temporary fixation of fractures with large fragments, e.g.:

- femoral neck fractures,
- intercondylar femoral fractures,
- epiphysiolysis of the femoral head,
- ankle arthrodeses,
- iliosacral dislocations,
- proximal tibia plateau fractures

## 7,3 mm Cannulated Screw 16 mm Threaded



| SS       | TI       | DIA. x LENGTH (mm) | SS       | TI       | DIA. x LENGTH (mm) |
|----------|----------|--------------------|----------|----------|--------------------|
| A6049050 | A6048050 | 7,3X50             | A6049090 | A6048090 | 7,3X90             |
| A6049055 | A6048055 | 7,3X55             | A6049095 | A6048095 | 7,3X95             |
| A6049060 | A6048060 | 7,3X60             | A6049100 | A6048100 | 7,3X100            |
| A6049065 | A6048065 | 7,3X65             | A6049105 | A6048105 | 7,3X105            |
| A6049070 | A6048070 | 7,3X70             | A6049110 | A6048110 | 7,3X110            |
| A6049075 | A6048075 | 7,3X75             | A6049115 | A6048115 | 7,3X115            |
| A6049080 | A6048080 | 7,3X80             | A6049120 | A6048120 | 7,3X120            |
| A6049085 | A6048085 | 7,3X85             |          |          |                    |



**Intended use :** Temporary fixation of fractures with large fragments, e.g.:

- femoral neck fractures,
- intercondylar femoral fractures,
- epiphysiolysis of the femoral head,
- ankle arthrodeses,
- iliosacral dislocations,
- proximal tibia plateau fractures

# Cannulated Screws

## 7,3 mm Cannulated Screw 32 mm Threaded




| SS       | TI       | DIA. x LENGTH (mm) | SS       | TI       | DIA. x LENGTH (mm) |
|----------|----------|--------------------|----------|----------|--------------------|
| A6059050 | A6058050 | 7,3X50             | A6059090 | A6058090 | 7,3X90             |
| A6059055 | A6058055 | 7,3X55             | A6059095 | A6058095 | 7,3X95             |
| A6059060 | A6058060 | 7,3X60             | A6059100 | A6058100 | 7,3X100            |
| A6059065 | A6058065 | 7,3X65             | A6059105 | A6058105 | 7,3X105            |
| A6059070 | A6058070 | 7,3X70             | A6059110 | A6058110 | 7,3X110            |
| A6059075 | A6058075 | 7,3X75             | A6059115 | A6058115 | 7,3X115            |
| A6059080 | A6058080 | 7,3X80             | A6059120 | A6058120 | 7,3X120            |
| A6059085 | A6058085 | 7,3X85             |          |          |                    |



**Intended use :** Temporary fixation of fractures with large fragments, e.g.:


- femoral neck fractures,
- intercondylar femoral fractures,
- epiphysiylolysis of the femoral head,
- ankle arthrodeses,
- iliosacral dislocations,
- proximal tibia plateau fractures

## Washer Ø13





| SS               | DIA. (mm) |
|------------------|-----------|
| A110 04 000 0014 | 13        |

## Washer Ø10



| SS               | DIA. (mm) |
|------------------|-----------|
| A110 04 000 0010 | 10        |

## Washer Ø8,0



| SS        | TI        | DIA. (mm) |
|-----------|-----------|-----------|
| A34004108 | A34004008 | 8         |

■ [www.aysam.com.tr](http://www.aysam.com.tr)  
■ [www.aysamstore.com.tr](http://www.aysamstore.com.tr)  
■ [aysam@aysam.com.tr](mailto:aysam@aysam.com.tr)

f [aysam.orthopedics](https://www.facebook.com/aysam.orthopedics)

t [aysamortho](https://www.tumblr.com/aysamortho)

in [aysam](https://www.instagram.com/aysam)

YouTube [aysam orthopaedics](https://www.youtube.com/aysamorthopaedics)



"We supply WORLD's  
Orthopedic & Trauma  
Implants..."

# Headless & Herbert Screws



# Headless, Herbert Screw & Interference Screws

## 3,0 mm Headless Cannulated Compression Screw



**Intended use :** Temporary fixation of fractures, malunions, nonunions, osteotomies, pseudoarthroses, arthrodesis, of the hand and foot bones(carpal, metacarpal, tarsal, metatarsal, phalanges bones), distal and proximal humerus, radius and ulna in adult and pediatric patients. Ligamentous injuries of small medium and large joints.

| TI       | DIA. x LENGTH (mm) |
|----------|--------------------|
| A6013014 | 3,0X14             |
| A6013016 | 3,0X16             |
| A6013018 | 3,0X18             |
| A6013020 | 3,0X20             |
| A6013022 | 3,0X22             |
| A6013024 | 3,0X24             |
| A6013026 | 3,0X26             |
| A6013028 | 3,0X28             |
| A6013030 | 3,0X30             |
| A6013035 | 3,0X35             |
| A6013040 | 3,0X40             |
| A6013045 | 3,0X45             |
| A6013050 | 3,0X50             |

## 3,5 mm Headless Cannulated Compression Screw



**Intended use :** Temporary fixation of fractures, malunions, nonunions, osteotomies, pseudoarthroses, arthrodesis, of the hand and foot bones(carpal, metacarpal, tarsal, metatarsal, phalanges bones), distal and proximal humerus, radius and ulna in adult and pediatric patients. Ligamentous injuries of small medium and large joints.

| TI       | DIA. x LENGTH (mm) |
|----------|--------------------|
| A6013514 | 3,5X14             |
| A6013516 | 3,5X16             |
| A6013518 | 3,5X18             |
| A6013520 | 3,5X20             |
| A6013522 | 3,5X22             |
| A6013524 | 3,5X24             |
| A6013526 | 3,5X26             |
| A6013528 | 3,5X28             |
| A6013530 | 3,5X30             |
| A6013535 | 3,5X35             |
| A6013540 | 3,5X40             |
| A6013545 | 3,5X45             |
| A6013550 | 3,5X50             |

## 4,0 mm Headless Cannulated Compression Screw



**Intended use :** Temporary fixation of fractures, malunions, nonunions, osteotomies, pseudoarthroses, arthrodesis, of the distal and proximal humerus, radius and ulna in adult and pediatric patients. ligamentous injuries of small medium and large joints. Fixation of fractures with medium fragments, e.g.:

- malleolar fractures
- tibia pilon fractures
- fractures of the calcaneus and talus
- tibial plateau fractures
- carpal and tarsal arthrodesis

| TI       | DIA. x LENGTH (mm) |
|----------|--------------------|
| A6014015 | 4,0X15             |
| A6014020 | 4,0X20             |
| A6014025 | 4,0X25             |
| A6014030 | 4,0X30             |
| A6014035 | 4,0X35             |
| A6014040 | 4,0X40             |
| A6014045 | 4,0X45             |
| A6014050 | 4,0X50             |
| A6014055 | 4,0X55             |
| A6014060 | 4,0X60             |
| A6014065 | 4,0X65             |
| A6014070 | 4,0X70             |
| A6014075 | 4,0X75             |

# Headless, Herbert Screw & Interference Screws

## 5,0 mm Headless Cannulated Compression Screw



**Intended use :** Temporary fixation of fractures, malunions, nonunions, osteotomies, pseudoarthroses, arthrodesis, of the, distal and proximal humerus, radius and ulna in adult and pediatric patients. ligamentous injuries of small medium and large joints. Fixation of fractures with medium fragments, e.g.:

- malleolar fractures
- tibia pilon fractures
- fractures of the calcaneus and talus
- tibial plateau fractures
- carpal and tarsal arthroses



| TI       | DIA. x LENGTH (mm) |
|----------|--------------------|
| A6015030 | 5,0X30             |
| A6015035 | 5,0X35             |
| A6015040 | 5,0X40             |
| A6015045 | 5,0X45             |
| A6015050 | 5,0X50             |
| A6015055 | 5,0X55             |
| A6015060 | 5,0X60             |
| A6015065 | 5,0X65             |
| A6015070 | 5,0X70             |
| A6015075 | 5,0X75             |
| A6015080 | 5,0X80             |

## 6,5 mm Headless Cannulated Compression Screw



**Intended use :** Temporary fixation of fractures with large fragments, e.g.:

- femoral neck fractures
- intercondylar femoral fractures
- epiphyseolysis of the femoral head
- ankle arthroses
- iliosacral dislocations



| TI       | DIA. x LENGTH (mm) |
|----------|--------------------|
| A6016540 | 6,5X40             |
| A6016545 | 6,5X45             |
| A6016550 | 6,5X50             |
| A6016555 | 6,5X55             |
| A6016560 | 6,5X60             |
| A6016565 | 6,5X65             |
| A6016570 | 6,5X70             |
| A6016575 | 6,5X75             |
| A6016580 | 6,5X80             |
| A6016585 | 6,5X85             |
| A6016590 | 6,5X90             |

## 3,0/4,0 mm Herbert Screw




**Intended use :** Temporary fixation of fractures, malunions, nonunions, osteotomies, pseudoarthroses, arthrodesis, of the, hand and foot bones(carpal, metacarpal, tarsal, metatarsal, phalanges bones), distal and proximal humerus, radius and ulna in adult and pediatric patients. Ligamentous injuries of small medium and large joints.



| TI      | LENGTH (mm) |
|---------|-------------|
| A34.114 | 14          |
| A34.116 | 16          |
| A34.118 | 18          |
| A34.120 | 20          |
| A34.122 | 22          |
| A34.124 | 24          |
| A34.126 | 26          |
| A34.128 | 28          |
| A34.130 | 30          |
| A34.132 | 32          |

# Headless, Herbert Screw & Interference Screws

## Interference Screw



| TI       | Diameter<br>(mm) | LENGTH<br>(mm) | LENGTH<br>(mm) | LENGTH<br>(mm) |
|----------|------------------|----------------|----------------|----------------|
| A6072710 | 7                | 20             | 25             | 30             |
| A6072712 | 8                | 20             | 25             | 30             |
| A6072714 | 9                | 20             | 25             | 30             |
| A6072716 | 10               | 20             | 25             | 30             |
| A6072732 | 11               | 20             | 25             | 30             |



**Intended use :** Temporary fixation of both bone-patellar tendon-bone (BTB) and soft tissue grafts in ACL and PCL reconstructions.




"We supply WORLD's  
Orthopedic & Trauma  
Implants..."

# Intramedullary Nails



# Intramedullary Nails/Humerus


## Humerus Intramedullary Nail



**Intended use :** Temporary fixation of fractures, osteotomies, nonunions, malunions, pseudoarthroses of humerus. Especially used in humeral shaft fractures.

| SS               | TI               | DIA x LENGTH (mm) |
|------------------|------------------|-------------------|
| A100 03 036 0718 | A100 03 035 0718 | 7X180             |
| A100 03 036 0720 | A100 03 035 0720 | 7X200             |
| A100 03 036 0722 | A100 03 035 0722 | 7X220             |
| A100 03 036 0724 | A100 03 035 0724 | 7X240             |
| A100 03 036 0726 | A100 03 035 0726 | 7X260             |
| A100 03 036 0728 | A100 03 035 0728 | 7X280             |
| A100 03 036 0730 | A100 03 035 0730 | 7X300             |
| A100 03 036 0732 | A100 03 035 0732 | 7X320             |
| A100 03 036 0818 | A100 03 035 0818 | 8X180             |
| A100 03 036 0820 | A100 03 035 0820 | 8X200             |
| A100 03 036 0822 | A100 03 035 0822 | 8X220             |
| A100 03 036 0824 | A100 03 035 0824 | 8X240             |
| A100 03 036 0826 | A100 03 035 0826 | 8X260             |
| A100 03 036 0828 | A100 03 035 0828 | 8X280             |
| A100 03 036 0830 | A100 03 035 0830 | 8X300             |
| A100 03 036 0832 | A100 03 035 0832 | 8X320             |
| A100 03 036 1018 | A100 03 035 1018 | 8,5X180           |
| A100 03 036 1020 | A100 03 035 1020 | 8,5X200           |
| A100 03 036 1022 | A100 03 035 1022 | 8,5X220           |
| A100 03 036 1024 | A100 03 035 1024 | 8,5X240           |
| A100 03 036 1026 | A100 03 035 1026 | 8,5X260           |
| A100 03 036 1028 | A100 03 035 1028 | 8,5X280           |
| A100 03 036 1030 | A100 03 035 1030 | 8,5X300           |
| A100 03 036 1032 | A100 03 035 1032 | 8,5X320           |
| A100 03 036 0918 | A100 03 035 0918 | 9X180             |
| A100 03 036 0920 | A100 03 035 0920 | 9X200             |
| A100 03 036 0922 | A100 03 035 0922 | 9X220             |
| A100 03 036 0924 | A100 03 035 0924 | 9X240             |
| A100 03 036 0926 | A100 03 035 0926 | 9X260             |
| A100 03 036 0928 | A100 03 035 0928 | 9X280             |
| A100 03 036 0930 | A100 03 035 0930 | 9X300             |
| A100 03 036 0932 | A100 03 035 0932 | 9X320             |

## 3,5 mm Intramedullary Nail Locking Screw



| SS               | TI               | DIA x LENGTH (mm) |
|------------------|------------------|-------------------|
| A101 03 036 4020 | A101 03 035 4020 | 3,5X20            |
| A101 03 036 4025 | A101 03 035 4025 | 3,5X25            |
| A101 03 036 4030 | A101 03 035 4030 | 3,5X30            |
| A101 03 036 4035 | A101 03 035 4035 | 3,5X35            |
| A101 03 036 4040 | A101 03 035 4040 | 3,5X40            |
| A101 03 036 4045 | A101 03 035 4045 | 3,5X45            |
| A101 03 036 4050 | A101 03 035 4050 | 3,5X50            |


## Humerus Intramedullary Nail End Cap



| SS               | TI               |
|------------------|------------------|
| A100 03 201 0000 | A100 03 202 0000 |

# Intramedullary Nails/Femur

## Femur Intramedullary Nail (ANTEGRADE)




### Intended use :

Temporary fixation of fractures, osteotomies, nonunions, malunions, pseudoarthroses of femur. Especially used in femoral shaft fractures.

| SS       | TI       | DIA x LENGTH (mm) | SS       | TI       | DIA x LENGTH (mm) |
|----------|----------|-------------------|----------|----------|-------------------|
| A4099320 | A4089320 | 9X320             | A4090320 | A4080320 | 10X320            |
| A4099340 | A4089340 | 9X340             | A4090340 | A4080340 | 10X340            |
| A4099360 | A4089360 | 9X360             | A4090360 | A4080360 | 10X360            |
| A4099380 | A4089380 | 9X380             | A4090380 | A4080380 | 10X380            |
| A4099400 | A4089400 | 9X400             | A4090400 | A4080400 | 10X400            |
| A4099420 | A4089420 | 9X420             | A4090420 | A4080420 | 10X420            |
| A4099440 | A4089440 | 9X440             | A4090440 | A4080440 | 10X440            |
| A4099460 | A4089460 | 9X460             | A4090460 | A4080460 | 10X460            |
| A4091320 | A4081320 | 11X320            | A4092320 | A4082320 | 12X320            |
| A4091340 | A4081340 | 11X340            | A4092340 | A4082340 | 12X340            |
| A4091360 | A4081360 | 11X360            | A4092360 | A4082360 | 12X360            |
| A4091380 | A4081380 | 11X380            | A4092380 | A4082380 | 12X380            |
| A4091400 | A4081400 | 11X400            | A4092400 | A4082400 | 12X400            |
| A4091420 | A4081420 | 11X420            | A4092420 | A4082420 | 12X420            |
| A4091440 | A4081440 | 11X440            | A4092440 | A4082440 | 12X440            |
| A4091460 | A4081460 | 11X460            | A4092460 | A4082460 | 12X460            |
| A4093320 | A4083320 | 13X320            | A4094320 | A4084320 | 14X320            |
| A4093340 | A4083340 | 13X340            | A4094340 | A4084340 | 14X340            |
| A4093360 | A4083360 | 13X360            | A4094360 | A4084360 | 14X360            |
| A4093380 | A4083380 | 13X380            | A4094380 | A4084380 | 14X380            |
| A4093400 | A4083400 | 13X400            | A4094400 | A4084400 | 14X400            |
| A4093420 | A4083420 | 13X420            | A4094420 | A4084420 | 14X420            |
| A4093440 | A4083440 | 13X440            | A4094440 | A4084440 | 14X440            |
| A4093460 | A4083460 | 13X460            | A4094460 | A4084460 | 14X460            |
| A4095320 | A4085320 | 15X320            | A4096320 | A4086320 | 16X320            |
| A4095340 | A4085340 | 15X340            | A4096340 | A4086340 | 16X340            |
| A4095360 | A4085360 | 15X360            | A4096360 | A4086360 | 16X360            |
| A4095380 | A4085380 | 15X380            | A4096380 | A4086380 | 16X380            |
| A4095400 | A4085400 | 15X400            | A4096400 | A4086400 | 16X400            |
| A4095420 | A4085420 | 15X420            | A4096420 | A4086420 | 16X420            |
| A4095440 | A4085440 | 15X440            | A4096440 | A4086440 | 16X440            |
| A4095460 | A4085460 | 15X460            | A4096460 | A4086460 | 16X460            |

## 4,5 mm Intramedullary Nail Locking Screw



| SS       | TI       | DIA x LENGTH (mm) |
|----------|----------|-------------------|
| A4194020 | A4184020 | 4,5X20            |
| A4194025 | A4184025 | 4,5X25            |
| A4194030 | A4184030 | 4,5X30            |
| A4194035 | A4184035 | 4,5X35            |
| A4194040 | A4184040 | 4,5X40            |
| A4194045 | A4184045 | 4,5X45            |
| A4194050 | A4184050 | 4,5X50            |
| A4194055 | A4184055 | 4,5X55            |
| A4194060 | A4184060 | 4,5X60            |
| A4194065 | A4184065 | 4,5X65            |
| A4194070 | A4184070 | 4,5X70            |
| A4194075 | A4184075 | 4,5X75            |


## Femur Intramedullary Nail End Cap



| SS               | TI               |
|------------------|------------------|
| A100 03 040 4900 | A100 03 040 5000 |

# Intramedullary Nails/Femur


## Distal Femur Intramedullary Nail (DFN / RETROGRADE)



**Intended use :**  
Temporary fixation of fractures, osteotomies, nonunions, malunions, pseudoarthroses of distal femur.


| SS       | TI       | DIA. x LENGTH (mm) |
|----------|----------|--------------------|
| A4499200 | A4489200 | 9X200              |
| A4499220 | A4489220 | 9X220              |
| A4499240 | A4489240 | 9X240              |
| A4499260 | A4489260 | 9X260              |
| A4490200 | A4480200 | 10X200             |
| A4490220 | A4480220 | 10X220             |
| A4490240 | A4480240 | 10X240             |
| A4490260 | A4480260 | 10X260             |
| A4491200 | A4481200 | 11X200             |
| A4491220 | A4481220 | 11X220             |
| A4491240 | A4481240 | 11X240             |
| A4491260 | A4481260 | 11X260             |
| A4490300 | A4480300 | 12X200             |
| A4490320 | A4480320 | 12X220             |
| A4490340 | A4480340 | 12X240             |
| A4490360 | A4480360 | 12X260             |
| A4491300 | A4481300 | 13X200             |
| A4491320 | A4481320 | 13X220             |
| A4491340 | A4481340 | 13X240             |
| A4491360 | A4481360 | 13X260             |

## 4,5 mm Intramedullary Nail Locking Screw




| SS       | TI       | DIA. x LENGTH (mm) |
|----------|----------|--------------------|
| A4194020 | A4184020 | 4,5X20             |
| A4194025 | A4184025 | 4,5X25             |
| A4194030 | A4184030 | 4,5X30             |
| A4194035 | A4184035 | 4,5X35             |
| A4194040 | A4184040 | 4,5X40             |
| A4194045 | A4184045 | 4,5X45             |
| A4194050 | A4184050 | 4,5X50             |
| A4194055 | A4184055 | 4,5X55             |
| A4194060 | A4184060 | 4,5X60             |
| A4194065 | A4184065 | 4,5X65             |
| A4194070 | A4184070 | 4,5X70             |
| A4194075 | A4184075 | 4,5X75             |

## Femur Short Reconstruction Intramedullary Nail ( PFN )



| TI       | DIA. x LENGTH (mm) |
|----------|--------------------|
| A4289240 | 9X240              |
| A4280240 | 10X240             |
| A4281240 | 11X240             |
| A3282240 | 12X240             |



**Intended use :** Temporary fixation of fractures, osteotomies, nonunions, malunions, pseudoarthroses of proximal femur. Especially used in ipsilateral pertrochanteric, intertrochanteric, subtrochanteric, femoral neck fractures.

## Distal Femur Intramedullary Nail (DFN / RETROGRADE) End Cap



| SS               | TI               |
|------------------|------------------|
| A100 03 040 4700 | A100 03 040 4800 |

# Intramedullary Nails/**Femur**

## Femur Long Reconstruction Intramedullary Nail



| TI       | DIA. x LENGTH (mm) | SIDE | TI       | DIA. x LENGTH (mm) | SIDE  |
|----------|--------------------|------|----------|--------------------|-------|
| A4389340 | 9X340              | Left | A4399340 | 9X340              | Right |
| A4389360 | 9X360              | Left | A4399360 | 9X360              | Right |
| A4389380 | 9X380              | Left | A4399380 | 9X380              | Right |
| A4389400 | 9X400              | Left | A4399400 | 9X400              | Right |
| A4389420 | 9X420              | Left | A4399420 | 9X420              | Right |
| A4380340 | 10X340             | Left | A4390340 | 10X340             | Right |
| A4380360 | 10X360             | Left | A4390360 | 10X360             | Right |
| A4380380 | 10X380             | Left | A4390380 | 10X380             | Right |
| A4380400 | 10X400             | Left | A4390400 | 10X400             | Right |
| A4380420 | 10X420             | Left | A4390420 | 10X420             | Right |
| A4381340 | 11X340             | Left | A4391340 | 11X340             | Right |
| A4381360 | 11X360             | Left | A4391360 | 11X360             | Right |
| A4381380 | 11X380             | Left | A4391380 | 11X380             | Right |
| A4381400 | 11X400             | Left | A4391400 | 11X400             | Right |
| A4381420 | 11X420             | Left | A4391420 | 11X420             | Right |
| A4382340 | 12X340             | Left | A4392340 | 12X340             | Right |
| A4382360 | 12X360             | Left | A4392360 | 12X360             | Right |
| A4382380 | 12X380             | Left | A4392380 | 12X380             | Right |
| A4382400 | 12X400             | Left | A4392400 | 12X400             | Right |
| A4382420 | 12X420             | Left | A4392420 | 12X420             | Right |
| A4383340 | 13X340             | Left | A4393340 | 13X340             | Right |
| A4383360 | 13X360             | Left | A4393360 | 13X360             | Right |
| A4383380 | 13X380             | Left | A4393380 | 13X380             | Right |
| A4383400 | 13X400             | Left | A4393400 | 13X400             | Right |
| A4383420 | 13X420             | Left | A4393420 | 13X420             | Right |

**Intended use :** Temporary fixation of fractures, osteotomies, nonunions, malunions, pseudoarthroses of proximal femur. Especially used in ipsilateral pertrochanteric, intertrochanteric, subtrochanteric, femoral neck fractures and femoral shaft fractures.

## Femur Long Reconstruction Intramedullary Nail Lag Screw



| TI       | DIA. x LENGTH (mm) | TI       | DIA. x LENGTH (mm) |
|----------|--------------------|----------|--------------------|
| A4385065 | 6,4X65             | A4385095 | 6,4X95             |
| A4385070 | 6,4X70             | A4385100 | 6,4X100            |
| A4385075 | 6,4X75             | A4385105 | 6,4X105            |
| A4385080 | 6,4X80             | A4385110 | 6,4X110            |
| A4385085 | 6,4X85             | A4385115 | 6,4X115            |
| A4385090 | 6,4X90             | A4385120 | 6,4X120            |

# Intramedullary Nails/Femur

## 4,5 mm Intramedullary Nail Locking Screw




| TI       | DIA. x LENGTH (mm) |
|----------|--------------------|
| A4184020 | 4,5X20             |
| A4184025 | 4,5X25             |
| A4184030 | 4,5X30             |
| A4184035 | 4,5X35             |
| A4184040 | 4,5X40             |
| A4184045 | 4,5X45             |
| A4184050 | 4,5X50             |
| A4184055 | 4,5X55             |
| A4184060 | 4,5X60             |
| A4184065 | 4,5X65             |
| A4184070 | 4,5X70             |
| A4184075 | 4,5X75             |

## Femur Long Reconstruction Intramedullary Nail End Cap



| TI       |
|----------|
| A3184030 |

## APFNA Femur Short Intramedullary Nail



| TI       | DIA. x LENGTH (mm) |
|----------|--------------------|
| A3080320 | 9X180              |
| A3080340 | 9X200              |
| A3080360 | 9X240              |
| A3080380 | 10X180             |
| A3080400 | 10X200             |
| A3080420 | 10X240             |
| A3080440 | 11X180             |
| A3080460 | 11X200             |
| A3081320 | 11X240             |
| A3081340 | 12X180             |
| A3081360 | 12X200             |
| A3081380 | 12X240             |
| A3081400 | 13X180             |
| A3081420 | 13X200             |
| A3081440 | 13X240             |
| A3081460 | 14X180             |
| A3082320 | 14X200             |
| A3082340 | 14X240             |

**Intended use :** Temporary fixation of fractures, osteotomies, nonunions, malunions, pseudoarthroses of proximal femur. Especially used in pertrochanteric, intertrochanteric, subtrochanteric, femoral neck fractures

## Intramedullary Nails/**Femur**

### APFNA Intramedullary Nail End Cap




| TI       | DIA. x LENGTH (mm) |
|----------|--------------------|
| A3184020 | 26                 |
| A3184025 | 40                 |

### APFNA Intramedullary Nail Blade Screw



| TI       | DIA. x LENGTH (mm) |
|----------|--------------------|
| A3395075 | 10,4X75            |
| A3395080 | 10,4X80            |
| A3395085 | 10,4X85            |
| A3395090 | 10,4 X90           |
| A3395095 | 10,4X95            |
| A3395100 | 10,4X100           |
| A3395105 | 10,4X105           |
| A3395110 | 10,4X110           |
| A3395115 | 10,4X115           |
| A3395120 | 10,4X120           |

### AROCOM Proximal Femur Anti-Rotation Compression Intramedullary Nail



| TI       | DIA. x LENGTH (mm) |
|----------|--------------------|
| A3089320 | 10X220             |
| A3089340 | 11X220             |
| A3089360 | 12X220             |
| A3089380 | 13X220             |

**Intended use :** Temporary fixation of fractures, osteotomies, nonunions, malunions, pseudoarthroses of proximal femur. Especially used in pertrochanteric, intertrochanteric, subtrochanteric, femoral neck fractures.

### AROCOM Proximal Femur Anti-Rotation Intramedullary Nail End Cap



| TI       |
|----------|
| A3184035 |



| TI       | LENGTH |
|----------|--------|
| A3184040 | 35     |

## Intramedullary Nails/**Femur**

### AROCOM Proximal Femur Anti-Rotation Intramedullary Nail Lag Screw



| TI       | DIA x LENGTH<br>(mm) | TI       | DIA. x LENGTH<br>(mm) |
|----------|----------------------|----------|-----------------------|
| A3082360 | 10X75                | A3083320 | 10X105                |
| A3082380 | 10X80                | A3083340 | 10X110                |
| A3082400 | 10X85                | A3083360 | 10X115                |
| A3082420 | 10X90                | A3083380 | 10X120                |
| A3082440 | 10X95                | A3083400 | 10X125                |
| A3082460 | 10X100               | A3083420 | 10X130                |

### AROCOM Proximal Femur Anti-Rotation Intramedullary Nail Anti-Rotation Screw



| TI       | DIA x LENGTH<br>(mm) | TI       | DIA. x LENGTH<br>(mm) |
|----------|----------------------|----------|-----------------------|
| A3083460 | 6X65                 | A3084440 | 6X100                 |
| A3084320 | 6X70                 | A3084460 | 6X105                 |
| A3084340 | 6X75                 | A3085320 | 6X110                 |
| A3084360 | 6X80                 | A3085340 | 6X115                 |
| A3084380 | 6X85                 | A3085360 | 6X120                 |
| A3084400 | 6X90                 | A3085380 | 6X125                 |
| A3084420 | 6X95                 | A3085400 | 6X130                 |

### AROCOM Proximal Femur Anti-Rotation Intramedullary Nail Connection Sleeve



| TI       | DIA. x LENGTH<br>(mm) |
|----------|-----------------------|
| A3083440 | 12X52                 |

### AROCOM Proximal Femur Anti-Rotation Intramedullary Nail Locking Screw



| TI       | DIA. x LENGTH<br>(mm) | TI       | DIA. x LENGTH<br>(mm) |
|----------|-----------------------|----------|-----------------------|
| A3085420 | 5.5X26                | A3086400 | 5,5X40                |
| A3085440 | 5.5X28                | A3086420 | 5,5X42                |
| A3085460 | 5.5X30                | A3086440 | 5,5X44                |
| A3086320 | 5.5X32                | A3086460 | 5,5X46                |
| A3086340 | 5.5X34                | A4283240 | 5,5X48                |
| A3086360 | 5.5X36                | A4284240 | 5,5X50                |
| A3086380 | 5.5X38                |          |                       |

## Intramedullary Nails/**Femur**

### ANTIRO Proximal Femur Anti-Rotation Intramedullary Nail



| TI       | DIA. x LENGTH (mm) |
|----------|--------------------|
| A3089400 | 10X220             |
| A3089420 | 11X220             |
| A3089440 | 12X220             |
| A3089460 | 13X220             |



**Intended use :** Temporary fixation of fractures, osteotomies, nonunions, malunions, pseudoarthroses of proximal femur. Especially used in pertrochanteric, intertrochanteric, subtrochanteric, femoral neck fractures.

### ANTIRO Proximal Femur Anti-Rotation Intramedullary Nail End Cap



| TI       |
|----------|
| A3184045 |

### ANTIRO Proximal Femur Anti-Rotation Intramedullary Nail Compression Screw



| TI       |
|----------|
| A3184050 |

### ANTIRO Proximal Femur Anti-Rotation Intramedullary Nail Lag Screw



| TI       | DIA. x LENGTH (mm) |
|----------|--------------------|
| A4384320 | 11,5X75            |
| A4384340 | 11,5X80            |
| A4384360 | 11,5X85            |
| A4384380 | 11,5X90            |
| A4384400 | 11,5X95            |
| A4384420 | 11,5X100           |
| A4384440 | 11,5X105           |
| A4384460 | 11,5X110           |
| A4385320 | 11,5X115           |
| A4385340 | 11,5X120           |

# Intramedullary Nails/Femur

## ANTIRO Proximal Femur Anti-Rotation Intramedullary Nail Pin



| TI       | DIA. x LENGTH (mm) |
|----------|--------------------|
| A4385360 | 5X75               |
| A4385380 | 5X80               |
| A4385400 | 5X85               |
| A4385420 | 5X90               |
| A4385440 | 5X95               |

| TI       | DIA. x LENGTH (mm) |
|----------|--------------------|
| A4385460 | 5X100              |
| A4386320 | 5X105              |
| A4386340 | 5X110              |
| A4386360 | 5X115              |
| A4386380 | 5X120              |

## ANTIRO Proximal Femur Anti-Rotation Intramedullary Nail Locking Screw



| TI        | DIA. x LENGTH (mm) |
|-----------|--------------------|
| A43879320 | 5X20               |
| A43879340 | 5X25               |
| A43879360 | 5X30               |
| A43879380 | 5X35               |
| A43879400 | 5X40               |
| A43879420 | 5X45               |



| TI        | DIA. x LENGTH (mm) |
|-----------|--------------------|
| A43879440 | 5X50               |
| A43879460 | 5X55               |
| A43870320 | 5X60               |
| A43870340 | 5X65               |
| A43870360 | 5X70               |
| A43870380 | 5X75               |



# Intramedullary Nails/Tibia

## Tibia Intramedullary Nail

| SS       | TI       | DIA. x LENGTH<br>(mm) | SS       | TI       | DIA. x LENGTH<br>(mm) |
|----------|----------|-----------------------|----------|----------|-----------------------|
| A4199260 | A4189260 | 8X255                 | A4192260 | A4182260 | 11X255                |
| A4199280 | A4189280 | 8X270                 | A4192280 | A4182280 | 11X270                |
| A4199300 | A4189300 | 8X285                 | A4192300 | A4182300 | 11X285                |
| A4199320 | A4189320 | 8X300                 | A4192320 | A4182320 | 11X300                |
| A4199340 | A4189340 | 8X315                 | A4192340 | A4182340 | 11X315                |
| A4199360 | A4189360 | 8X330                 | A4192360 | A4182360 | 11X330                |
| A4199380 | A4189380 | 8X345                 | A4192380 | A4182380 | 11X345                |
| A4199400 | A4189400 | 8X360                 | A4192400 | A4182400 | 11X360                |
| A3199340 | A3189340 | 8X375                 | A3192340 | A3182340 | 11X375                |
| A3199360 | A3189360 | 8X390                 | A3192360 | A3182360 | 11X390                |
| A3199380 | A3189380 | 8X405                 | A3192380 | A3182380 | 11X405                |
| A3199400 | A3189400 | 8X420                 | A3192400 | A3182400 | 11X420                |
| A4190260 | A4180260 | 9X255                 | A4193260 | A4183260 | 12X255                |
| A4190280 | A4180280 | 9X270                 | A4193280 | A4183280 | 12X270                |
| A4190300 | A4180300 | 9X285                 | A4193300 | A4183300 | 12X285                |
| A4190320 | A4180320 | 9X300                 | A4193320 | A4183320 | 12X300                |
| A4190340 | A4180340 | 9X315                 | A4193340 | A4183340 | 12X315                |
| A4190360 | A4180360 | 9X330                 | A4193360 | A4183360 | 12X330                |
| A4190380 | A4180380 | 9X345                 | A4193380 | A4183380 | 12X345                |
| A4190400 | A4180400 | 9X360                 | A4193400 | A4183400 | 12X360                |
| A3190340 | A3180340 | 9X375                 | A3193340 | A3183340 | 12X375                |
| A3190360 | A3180360 | 9X390                 | A3193360 | A3183360 | 12X390                |
| A3190380 | A3180380 | 9X405                 | A3193380 | A3183380 | 12X405                |
| A3190400 | A3180400 | 9X420                 | A3193400 | A3183400 | 12X420                |
| A4191260 | A4181260 | 10X255                | A4194260 | A4184260 | 13X255                |
| A4191280 | A4181280 | 10X270                | A4194280 | A4184280 | 13X270                |
| A4191300 | A4181300 | 10X285                | A4194300 | A4184300 | 13X285                |
| A4191320 | A4181320 | 10X300                | A4194320 | A4184320 | 13X300                |
| A4191340 | A4181340 | 10X315                | A4194340 | A4184340 | 13X315                |
| A4191360 | A4181360 | 10X330                | A4194360 | A4184360 | 13X330                |
| A4191380 | A4181380 | 10X345                | A4194380 | A4184380 | 13X345                |
| A4191400 | A4181400 | 10X360                | A4194400 | A4184400 | 13X360                |
| A3191340 | A3181340 | 10X375                | A3194340 | A3184340 | 13X375                |
| A3191360 | A3181360 | 10X390                | A3194360 | A3184360 | 13X390                |
| A3191380 | A3181380 | 10X405                | A3194380 | A3184380 | 13X405                |
| A3191400 | A3181400 | 10X420                | A3194400 | A3184400 | 13X420                |




**Intended use :** Temporary fixation of fractures, osteotomies, nonunions, malunions, pseudoarthroses of tibia. Espceially used for tibia shaft fractures.

# Intramedullary Nails/Tibia


## ANT Tibia Intramedullary Nail

| TI               | DIA. x LENGTH<br>(mm) | TI               | DIA. x LENGTH<br>(mm) |
|------------------|-----------------------|------------------|-----------------------|
| A100 03 037 2708 | 8X270                 | A100 03 037 2709 | 9X270                 |
| A100 03 037 2808 | 8X285                 | A100 03 037 2809 | 9X285                 |
| A100 03 037 3008 | 8X300                 | A100 03 037 3009 | 9X300                 |
| A100 03 037 3108 | 8X315                 | A100 03 037 3109 | 9X315                 |
| A100 03 037 3308 | 8X330                 | A100 03 037 3309 | 9X330                 |
| A100 03 037 3408 | 8X345                 | A100 03 037 3409 | 9X345                 |
| A100 03 037 3608 | 8X360                 | A100 03 037 3609 | 9X360                 |
| A100 03 037 3808 | 8X375                 | A100 03 037 3809 | 9X375                 |
| A100 03 037 4008 | 8X390                 | A100 03 037 4009 | 9X390                 |
| A100 03 037 4208 | 8X405                 | A100 03 037 4209 | 9X405                 |
| A100 03 037 2710 | 10X270                | A100 03 037 4409 | 9X420                 |
| A100 03 037 2810 | 10X285                | A100 03 037 2711 | 11X270                |
| A100 03 037 3010 | 10X300                | A100 03 037 2811 | 11X285                |
| A100 03 037 3110 | 10X315                | A100 03 037 3011 | 11X300                |
| A100 03 037 3310 | 10X330                | A100 03 037 3111 | 11X315                |
| A100 03 037 3410 | 10X345                | A100 03 037 3311 | 11X330                |
| A100 03 037 3610 | 10X360                | A100 03 037 3411 | 11X345                |
| A100 03 037 3810 | 10X375                | A100 03 037 3611 | 11X360                |
| A100 03 037 4010 | 10X390                | A100 03 037 3811 | 11X375                |
| A100 03 037 4210 | 10X405                | A100 03 037 4011 | 11X390                |
| A100 03 037 4410 | 10X420                | A100 03 037 4211 | 11X405                |
| A100 03 037 2712 | 12X270                | A100 03 037 4411 | 11X420                |
| A100 03 037 2812 | 12X285                | A100 03 037 2813 | 13X285                |
| A100 03 037 3012 | 12X300                | A100 03 037 3013 | 13X300                |
| A100 03 037 3112 | 12X315                | A100 03 037 3113 | 13X315                |
| A100 03 037 3312 | 12X330                | A100 03 037 3313 | 13X330                |
| A100 03 037 3412 | 12X345                | A100 03 037 3413 | 13X345                |
| A100 03 037 3612 | 12X360                | A100 03 037 3613 | 13X360                |
| A100 03 037 3812 | 12X375                | A100 03 037 3813 | 13X375                |
| A100 03 037 4012 | 12X390                | A100 03 037 4013 | 13X390                |
| A100 03 037 4212 | 12X405                | A100 03 037 4213 | 13X405                |
| A100 03 037 4412 | 12X420                | A100 03 037 4413 | 13X420                |



Temporary fixation of fractures, osteotomies, nonunions, malunions, pseudoarthroses of tibia. especially used for fixation of ipsilateral tibia plateau and shaft fractures.

## 4,5 mm Intramedullary Nail Locking Screw



| TI       | DIA x LENGTH<br>(mm) | TI       | DIA x LENGTH<br>(mm) |
|----------|----------------------|----------|----------------------|
| A4184020 | 4,5X20               | A4184050 | 4,5X50               |
| A4184025 | 4,5X25               | A4184055 | 4,5X55               |
| A4184030 | 4,5X30               | A4184060 | 4,5X60               |
| A4184035 | 4,5X35               | A4184065 | 4,5X65               |
| A4184040 | 4,5X40               | A4184070 | 4,5X70               |
| A4184045 | 4,5X45               | A4184075 | 4,5X75               |

## ANT Tibia Intramedullary Nail End Cap



TI  
A100 03 040 4600

"We supply WORLD's  
Orthopedic & Trauma  
Implants..."


# Elastic Nails



aysam

# Elastic Nails

## Elastic Nail




| TI        | DIAMETER (mm) | LENGTH (mm) |
|-----------|---------------|-------------|
| A45015440 | 1,5           | 440         |
| A45020440 | 2,0           | 440         |
| A45025440 | 2,5           | 440         |
| A45030440 | 3,0           | 440         |
| A45035440 | 3,5           | 440         |
| A45040440 | 4,0           | 440         |
| A45045440 | 4,5           | 440         |
| A45050440 | 5,0           | 440         |
| A45055440 | 5,5           | 440         |



**Intended use :** Temporary fixation of fractures, osteotomies, malunions, nonunions, pseduarthroses of long bones (humerus, ulna, radius, femur, tibia, fibula). Especially used in pediatric and adolescent patients.

## Elastic Nail End Cap



| TI     | DIAMETER (mm) |
|--------|---------------|
| A45015 | 1,5           |
| A45020 | 2,0           |
| A45025 | 2,5           |
| A45030 | 3,0           |
| A45035 | 3,5           |
| A45040 | 4,0           |
| A45045 | 4,5           |
| A45050 | 5,0           |
| A45055 | 5,5           |

"We supply WORLD's  
Orthopedic & Trauma  
Implants..."

# Spinal Systems



## Asino Fixed Angle Screw



**Intended use :** Temporary fixation of fractures, pseudoarthrosis, spondylolisthesis of vertebrae.



| TI       | DIA x LENGTH (mm) |
|----------|-------------------|
| A8124025 | 4,0X25            |
| A8124030 | 4,0X30            |
| A8124035 | 4,0X35            |
| A8124040 | 4,0X40            |
| A8124525 | 4,5X25            |
| A8124530 | 4,5X30            |
| A8124535 | 4,5X35            |
| A8124540 | 4,5X40            |
| A8124545 | 4,5X45            |
| A8125030 | 5,0X30            |
| A8125035 | 5,0X35            |
| A8125040 | 5,0X40            |
| A8125045 | 5,0X45            |
| A8125050 | 5,0X50            |
| A8125530 | 5,5X30            |
| A8125535 | 5,5X35            |
| A8125540 | 5,5X40            |
| A8125545 | 5,5X45            |
| A8125550 | 5,5X50            |
| A8126030 | 6,0X30            |
| A8126035 | 6,0X35            |
| A8126040 | 6,0X40            |
| A8126045 | 6,0X45            |
| A8126050 | 6,0X50            |
| A8126055 | 6,0X55            |

| TI       | DIA x LENGTH (mm) |
|----------|-------------------|
| A8126530 | 6,5X30            |
| A8126535 | 6,5X35            |
| A8126540 | 6,5X40            |
| A8126545 | 6,5X45            |
| A8126550 | 6,5X50            |
| A8126555 | 6,5X55            |
| A8126560 | 6,5X60            |
| A8126565 | 6,5X65            |
| A8127030 | 7,0X30            |
| A8127035 | 7,0X35            |
| A8127040 | 7,0X40            |
| A8127045 | 7,0X45            |
| A8127050 | 7,0X50            |
| A8127055 | 7,0X55            |
| A8127060 | 7,0X60            |
| A8127065 | 7,0X65            |
| A8127530 | 7,5X30            |
| A8127535 | 7,5X35            |
| A8127540 | 7,5X40            |
| A8127545 | 7,5X45            |
| A8127550 | 7,5X50            |
| A8127555 | 7,5X55            |
| A8127560 | 7,5X60            |
| A8127565 | 7,5X65mm          |

## Asino Fixed-Angle Reduction Screw



**Intended use :** Temporary fixation of fractures, pseudoarthrosis, spondylolisthesis of vertebrae. Especially used for reduction of displaced vertebrae.



| TI       | DIA x LENGTH (mm) |
|----------|-------------------|
| A8164025 | 4,0X25            |
| A8164030 | 4,0X30            |
| A8164035 | 4,0X35            |
| A8164040 | 4,0X40            |
| A8164530 | 4,5X30            |
| A8164535 | 4,5X35            |
| A8164540 | 4,5X40            |
| A8164545 | 4,5X45            |
| A8165030 | 5,0X30            |
| A8165035 | 5,0X35            |
| A8165040 | 5,0X40            |
| A8165045 | 5,0X45            |
| A8165530 | 5,5X30            |
| A8165535 | 5,5X35            |
| A8165540 | 5,5X40            |
| A8165545 | 5,5X45            |
| A8165550 | 5,5X50            |

| TI       | DIA x LENGTH (mm) |
|----------|-------------------|
| A8166030 | 6,0X30            |
| A8166035 | 6,0X35            |
| A8166040 | 6,0X40            |
| A8166045 | 6,0X45            |
| A8166050 | 6,0X50            |
| A8166530 | 6,5X30            |
| A8166535 | 6,5X35            |
| A8166540 | 6,5X40            |
| A8166545 | 6,5X45            |
| A8166550 | 6,5X50            |
| A8166555 | 6,5X55            |
| A8167030 | 7,0X30            |
| A8167035 | 7,0X35            |
| A8167040 | 7,0X40            |
| A8167045 | 7,0X45            |
| A8167050 | 7,0X50            |

# Spinal Systems

## Asino Multi-Axial Screw



**Intended use :** Temporary fixation of fractures, pseudoarthrosis, spondylolisthesis of vertebrae.

| TI       | DIA x LENGTH (mm) |
|----------|-------------------|
| A8104025 | 4,0X25            |
| A8104030 | 4,0X30            |
| A8104035 | 4,0X35            |
| A8104040 | 4,0X40            |
| A8104530 | 4,5X30            |
| A8104535 | 4,5X35            |
| A8104540 | 4,5X40            |
| A8104545 | 4,5X45            |
| A8105030 | 5,0X30            |
| A8105035 | 5,0X35            |
| A8105040 | 5,0X40            |
| A8105045 | 5,0X45            |
| A8105530 | 5,5X30            |
| A8105535 | 5,5X35            |
| A8105540 | 5,5X40            |
| A8105545 | 5,5X45            |
| A8105550 | 5,5X50            |
| A8106030 | 6,0X30            |
| A8106035 | 6,0X35            |
| A8106040 | 6,0X40            |
| A8106045 | 6,0X45            |
| A8106050 | 6,0X50            |
| A8106055 | 6,0X55            |

| TI       | DIA x LENGTH (mm) |
|----------|-------------------|
| A8106530 | 6,5X30            |
| A8106535 | 6,5X35            |
| A8106540 | 6,5X40            |
| A8106545 | 6,5X45            |
| A8106550 | 6,5X50            |
| A8106555 | 6,5X55            |
| A8107030 | 7,0X30            |
| A8107035 | 7,0X35            |
| A8107040 | 7,0X40            |
| A8107045 | 7,0X45            |
| A8107050 | 7,0X50            |
| A8107530 | 7,5X30            |
| A8107535 | 7,5X35            |
| A8107540 | 7,5X40            |
| A8107545 | 7,5X45            |
| A8107550 | 7,5X50            |
| A8107555 | 7,5X55            |
| A8107560 | 7,5X60            |
| A8107565 | 7,5X65            |

## Asino Multi-Axial Reduction Screw



**Intended use :** Temporary fixation of fractures, pseudoarthrosis, spondylolisthesis of vertebrae. especially used for reduction of displaced vertebrae.


| TI       | DIA x LENGTH (mm) |
|----------|-------------------|
| A8144025 | 4,0X25            |
| A8144030 | 4,0X30            |
| A8144035 | 4,0X35            |
| A8144040 | 4,0X40            |
| A8144530 | 4,5X30            |
| A8144535 | 4,5X35            |
| A8144540 | 4,5X40            |
| A8144545 | 4,5X45            |
| A8145030 | 5,0X30            |
| A8145035 | 5,0X35            |
| A8145040 | 5,0X40            |
| A8145045 | 5,0X45            |
| A8145050 | 5,0X50            |
| A8145530 | 5,5X30            |
| A8145535 | 5,5X35            |
| A8145540 | 5,5X40            |
| A8145545 | 5,5X45            |

| TI       | DIA x LENGTH (mm) |
|----------|-------------------|
| A8146030 | 6,0X30            |
| A8146035 | 6,0X35            |
| A8146040 | 6,0X40            |
| A8146045 | 6,0X45            |
| A8146050 | 6,0X50            |
| A8146530 | 6,5X30            |
| A8146535 | 6,5X35            |
| A8146540 | 6,5X40            |
| A8146545 | 6,5X45            |
| A8146550 | 6,5X50            |
| A8147035 | 7,0X35            |
| A8147040 | 7,0X40            |
| A8147045 | 7,0X45            |

Asino Cross Link Plate



| TI       | DIA x LENGTH<br>(mm) |
|----------|----------------------|
| A8482001 | 6,0X28               |
| A8482002 | 6,0X48               |
| A8482003 | 6,0X68               |
| A8482004 | 6,0X88               |



**Intended use :** Temporary fixation of fractures, pseudoarthrosis, spondylolisthesis of vertebrae. Used for purpose of cross linking of two rods.

6,0 mm Hook



| TI       | TYPE                    |
|----------|-------------------------|
| A8460001 | Pedicle Hook            |
| A8460002 | Laminar Hook            |
| A8460003 | Transverse Process Hook |




**Intended use :** Temporary fixation of fractures, pseudoarthrosis, spondylolisthesis of vertebrae.

Set Screw



| TI       |
|----------|
| A8420600 |




# Spinal Systems

## Lateral Connector 6,0 mm



TI

A8426006




## Connector 6,0 mm





TI

A8440810




## Rod 6,0 mm



| TI       | DIA x LENGTH (mm) |
|----------|-------------------|
| A8001030 | 6,0X30            |
| A8001035 | 6,0X35            |
| A8001040 | 6,0X40            |
| A8001050 | 6,0X50            |
| A8001060 | 6,0X60            |
| A8001070 | 6,0X70            |
| A8001080 | 6,0X80            |
| A8001090 | 6,0X90            |
| A8001100 | 6,0X100           |
| A8001110 | 6,0X110           |
| A8001250 | 6,0X250           |
| A8001500 | 6,0X500           |

# Spinal Systems

## Pediatric Rod 5,0 mm




| TI               | DIA x LENGTH (mm) |
|------------------|-------------------|
| A425 20 060 0060 | 5,0X60            |
| A425 20 060 0080 | 5,0X80            |
| A425 20 060 0100 | 5,0X100           |
| A425 20 060 0110 | 5,0X110           |
| A425 20 060 0120 | 5,0X120           |
| A425 20 060 0130 | 5,0X130           |
| A425 20 060 0140 | 5,0X140           |
| A425 20 060 0150 | 5,0X150           |
| A425 20 060 0160 | 5,0X160           |
| A425 20 060 0180 | 5,0X180           |


| TI               | DIA x LENGTH (mm) |
|------------------|-------------------|
| A425 20 060 0200 | 5,0X200           |
| A425 20 060 0240 | 5,0X240           |
| A425 20 060 0260 | 5,0X260           |
| A425 20 060 0280 | 5,0X280           |
| A425 20 060 0300 | 5,0X300           |
| A425 20 060 0320 | 5,0X320           |
| A425 20 060 0340 | 5,0X340           |
| A425 20 060 0360 | 5,0X360           |
| A425 20 060 0380 | 5,0X380           |



## Anterior Cervical Plate 1



| TI               | LENGTH (mm) |
|------------------|-------------|
| A433 02 011 0020 | 26          |
| A433 02 011 0030 | 30          |
| A433 02 011 0040 | 33          |
| A433 02 011 0050 | 36          |
| A433 02 011 0060 | 39          |
| A433 02 011 0070 | 42          |
| A433 02 011 0080 | 45          |
| A433 02 013 0020 | 48          |
| A433 02 013 0030 | 51          |
| A433 02 013 0040 | 54          |
| A433 02 013 0050 | 57          |
| A433 02 013 0060 | 61          |
| A433 02 013 0070 | 64          |
| A433 02 013 0080 | 70          |



**Intended use :** Temporary fixation of fractures, pseudoarthrosis, spondylolisthesis, fusion of cervical vertebrae.

## Anterior Cervical Screw-Plate Locking Apparatus




| TI               |
|------------------|
| A442 00 000 0001 |

## Cervical Screw



| TI       | DIA x LENGTH (mm) |
|----------|-------------------|
| A6133528 | 3,5X12            |
| A6133530 | 3,5X14            |
| A6133532 | 3,5X16            |
| A6133534 | 3,5X18            |
| A6133536 | 3,5X20            |

## Anterior Cervical Plate 2



**Intended use :** Temporary fixation of fractures, pseudoarthrosis, spondylolisthesis, fusion of cervical vertebrae.

| TI               | LENGTH (mm) |
|------------------|-------------|
| A433 02 015 0020 | 25          |
| A433 02 015 0030 | 27,5        |
| A433 02 015 0040 | 30          |
| A433 02 015 0050 | 32,5        |
| A433 02 015 0060 | 35          |
| A433 02 015 0070 | 37,5        |
| A433 02 015 0080 | 40          |
| A438 02 002 0070 | 45          |
| A438 02 002 0075 | 50          |
| A438 02 002 0080 | 55          |
| A438 02 002 0085 | 60          |
| A438 02 002 0090 | 65          |
| A438 02 002 0095 | 70          |
| A438 02 002 0100 | 75          |

## Cervical Expansion Head Locking Screw



**Cervical Screw Head Expansion Aparatus**




| TI               | DIA x LENGTH (mm) |
|------------------|-------------------|
| A400 02 051 0030 | 4,0X13            |
| A400 02 051 0034 | 4,0X14            |
| A400 02 051 0036 | 4,0X15            |
| A400 02 051 0040 | 4,0X16            |
| A400 02 051 0044 | 4,0X17            |
| A400 02 052 0025 | 4,0X18            |

Anterior Thoracic and Lumbar Plate



| TI               | LENGTH<br>(mm) |
|------------------|----------------|
| A400 02 440 0040 | 60             |
| A400 02 440 0050 | 65             |
| A400 02 440 0060 | 70             |
| A400 02 440 0070 | 75             |
| A400 02 440 0080 | 80             |
| A400 02 440 0090 | 85             |
| A400 02 440 0100 | 90             |
| A400 02 440 0110 | 100            |



**Intended use :** Temporary fixation of fractures, pseudoarthrosis, spondylolisthesis, fusion of thoracic and lumbar vertebrae.

Anterior Thoracic and Lumbar Bolt



| TI       | DIA x LENGTH<br>(mm) |
|----------|----------------------|
| A6134510 | 6,5X30               |
| A6134512 | 6,5X35               |
| A6134514 | 6,5X40               |
| A6134516 | 6,5X45               |
| A6134518 | 6,5X50               |



Anterior Thoracic and Lumbar Screw




| TI       | DIA x LENGTH<br>(mm) |
|----------|----------------------|
| A6144014 | 6,5X30               |
| A6144016 | 6,5X35               |
| A6144018 | 6,5X40               |
| A6144020 | 6,5X45               |
| A6144022 | 6,5X50               |




"We supply WORLD's  
Orthopedic & Trauma  
Implants..."

# External Fixators



# Ayfix Dynamic Axial External Fixators

## Ayfix Dynamic Axial External Fixator Femur-Tibia




| CAT. NO          | TYPE   | DISTRACTION (mm) | LENGTH (mm) |
|------------------|--------|------------------|-------------|
| A600 11 100 0001 | Long   | 60               | 270         |
| A600 11 100 0002 | Medium | 40               | 250         |
| A600 11 100 0003 | Short  | 20               | 230         |

**Intended use :** Temporary fixation of fractures, nonunions, malunions, osteotomies, pseduarthrosis, of femur and tibia, arthrodesis of hip, knee, ankle joints. Especially used for fixation of open fractures of femur, tibia and for distraction and compression purposes on fracture site.

| available as aluminium/stainless steel/peek.

## Ayfix Dynamic Axial External Fixator Femur-Tibia T-Clamp



| CAT. NO          | TYPE   | DISTRACTION (mm) | LENGTH (mm) |
|------------------|--------|------------------|-------------|
| A600 16 100 0006 | Long   | 60               | 240         |
| A600 16 100 0007 | Medium | 40               | 220         |
| A600 16 100 0008 | Short  | 20               | 200         |

**Intended use :** Temporary fixation of fractures, nonunions, malunions, osteotomies, pseduarthrosis, of femur and tibia, arthrodesis of hip, knee, ankle joints. Especially used for fixation of open fractures of proximal and distal femur-tibia and for distraction and compression purposes on fracture site.

| available as aluminium/stainless steel/peek.


# Ayfix Dynamic Axial External Fixators

## Ayfix Dynamic Axial External Fixator Trochanteric Clamp (135°)



**Half-Pin (Schanz Screw) to be used**  
Ø (5,0/5,5/6,0 mm)

**Allen Wrench to be used**  
Ø (6,0 mm)



| CAT. NO          | TYPE   | DISTRACTION (mm) | LENGTH (mm) |
|------------------|--------|------------------|-------------|
| A600 11 100 0040 | Long   | 60               | 270         |
| A600 11 100 0041 | Medium | 40               | 250         |
| A600 11 100 0042 | Short  | 20               | 230         |

**Intended use :** Temporary fixation of fractures, nonunions, malunions, osteotomies, pseduarthrosis, of femur. Especially used for fixation of fractures of proximal femur (intertrochanteric, pertrochanteric) in elderly patients for early mobilization, for distraction and compression purposes on fracture site.


| available as aluminium/stainless steel/peek.

## Ayfix Dynamic Axial External Fixator Square Clamp Femur-Tibia



**Half-Pin (Schanz Screw) to be used**  
Ø (5,0/5,5/6,0 mm)

**Allen Wrench to be used**  
Ø (6,0 mm)



| CAT. NO          | TYPE   | DISTRACTION (mm) | LENGTH (mm) |
|------------------|--------|------------------|-------------|
| A600 11 100 0140 | Long   | 60               | 270         |
| A600 11 100 0141 | Medium | 40               | 250         |
| A600 11 100 0142 | Short  | 20               | 230         |

**Intended use :** Temporary fixation of fractures, nonunions, malunions, osteotomies, pseduarthrosis, of femur and tibia, arthrodesis of hip, knee, ankle joints. Especially used for fixation of open fractures of femur, tibia and for distraction and compression purposes on fracture site.

| available as aluminium/stainless steel/peek.


# Ayfix Dynamic Axial External Fixators

## Ayfix Limb Lengthening External Fixator Femur-Tibia



Half-Pin (Schanz Screw) to be used  
Ø (5,0/5,5/6,0 mm)

Allen Wrench to be used  
Ø (6,0 mm)



| CAT. NO          | TYPE   | DISTRACTION (mm) | LENGTH (mm) |
|------------------|--------|------------------|-------------|
| A600 11 100 0020 | Long   | 97               | 250         |
| A600 11 100 0021 | Medium | 77               | 230         |
| A600 11 100 0022 | Short  | 57               | 210         |

**Intended use :** Temporary fixation of fractures, nonunions, malunions, osteotomies, pseduarthrosis, of femur and tibia, arthrodesis of hip knee ankle joints. Especially used for fixation of open fractures of femur, tibia and for distraction and compression purposes on fracture site also used for lengthening purposes of femur and tibia bones.


| available as aluminium/stainless steel/peek.

## Ayfix Dynamic Axial Ankle Fixator



Half-Pin (Schanz Screw) to be used  
Ø (5,0/5,5/6,0 mm)

Allen Wrench to be used  
Ø (6,0 mm)



| CAT. NO          | TYPE           | DISTRACTION (mm) | LENGTH (mm) |
|------------------|----------------|------------------|-------------|
| A600 11 100 0004 | Left - Long    | 80               | 275         |
| A600 11 100 0005 | Left - Medium  | 60               | 255         |
| A600 11 100 0006 | Left - Short   | 40               | 235         |
| A600 11 101 0004 | Right - Long   | 80               | 275         |
| A600 11 101 0005 | Right - Medium | 60               | 255         |
| A600 11 101 0006 | Right - Short  | 40               | 235         |

**Intended use :** Temporary fixation of fractures, nonunions, malunions, osteotomies, pseduarthrosis, of distal tibia, arthrodesis of ankle joints. Especially used for fixation of open fractures of distal tibia and for distraction and compression purposes on fracture site.

| available as aluminium/stainless steel/peek.

# Ayfix Dynamic Axial External Fixators


## Ayfix Dynamic Axial External Fixator Humerus



**Half-Pin (Schanz Screw) to be used**  
Ø (4,0/4,5/5,0 mm)

**Allen Wrench to be used**  
Ø (5,0 mm)

| CAT. NO          | TYPE   | DISTRACTION (mm) | LENGTH (mm) |
|------------------|--------|------------------|-------------|
| A600 16 100 0001 | Long   | 51               | 235         |
| A600 16 100 0002 | Medium | 31               | 215         |
| A600 16 100 0003 | Short  | 11               | 195         |



**Intended use :** Temporary fixation of fractures, nonunions, malunions, osteotomies, pseduarthrosis, of humerus, arthrodesis of elbow joint. Especially used for fixation of open fractures of humerus and for distraction and compression purposes on fracture site.

| available as aluminium/stainless steel/peek.


## Ayfix Dynamic Axial External Fixator T-Clamp Humerus



**Half-Pin (Schanz Screw) to be used**  
Ø (4,0/4,5/5,0 mm)

**Allen Wrench to be used**  
Ø (5,0 mm)

| CAT. NO          | TYPE   | DISTRACTION (mm) | LENGTH (mm) |
|------------------|--------|------------------|-------------|
| A600 16 101 0001 | Long   | 51               | 207         |
| A600 16 101 0002 | Medium | 31               | 187         |
| A600 16 101 0003 | Short  | 11               | 167         |



**Intended use :** Temporary fixation of fractures, nonunions, malunions, osteotomies, pseduarthrosis, of proximal humerus, arthrodesis of shoulder joint. Especially used for fixation of open fractures of proximal humerus and for distraction and compression purposes on fracture

| available as aluminium/stainless steel/peek.

# Ayfix Dynamic Axial External Fixators


## Ayfix Limb Lengthening External Fixator Humerus



Half-Pin (Schanz Screw) to be used  
Ø (4,0/4,5/5,0 mm)

Allen Wrench to be used  
Ø (5,0 mm)

| CAT. NO          | TYPE   | DISTRACTION (mm) | LENGTH (mm) |
|------------------|--------|------------------|-------------|
| A600 16 100 0010 | Long   | 81               | 215         |
| A600 16 100 0011 | Medium | 61               | 195         |
| A600 16 100 0012 | Short  | 41               | 175         |



**Intended use :** Temporary fixation of fractures, nonunions, malunions, osteotomies, pseduarthrosis, of humerus, arthrodesis of elbow joint. Especially used for fixation of open fractures of humerus and for distraction and compression purposes on fracture site also used for lengthening purposes of humerus bone.


| available as aluminium/stainless steel/peek.

## Ayfix Dynamic External Fixator Forearm

Allen Wrench to be used  
Ø (5,0 mm)



Half-Pin (Schanz Screw) to be used  
Ø (3,5/4,0/4,5 mm)



| CAT. NO          | TYPE   | DISTRACTION (mm) | LENGTH (mm) |
|------------------|--------|------------------|-------------|
| A600 17 100 0001 | Long   | 42               | 200         |
| A600 17 100 0002 | Medium | 22               | 180         |
| A600 17 100 0003 | Short  | 14               | 168         |

**Intended use :** Temporary fixation of fractures, nonunions, malunions, osteotomies, pseduarthrosis, of ulna and radius, arthrodesis of wrist joint. Especially used for fixation of open fractures of ulna and radius and for distraction and compression purposes on fracture site. Also used for lengthening purposes ulna and radius. It is used for temporary fixation of PEDIATRIC FEMUR & TIBIA FRACTURES .

| available as aluminium/stainless steel/peek.

# Ayfix Dynamic Axial External Fixators


## Ayfix Colles Type Wrist Fixator



**Allen Wrench to be used**  
Ø (3,0/4,0/5,0 mm)

**Half-Pin (Schanz Screw) to be used**  
Ø (2,5/3,0/3,5 mm)

| CAT. NO.         | LENGTH (mm) |
|------------------|-------------|
| A600 12 110 0010 | 218         |



**Intended use :** Temporary fixation of fractures, nonunions, malunions, osteotomies, pseduarthrosis, of distal radius, arthrodesis of wrist joint. Especially used for fixation of open and/or communited intra-articular fractures of distal radius and for distraction and compression purposes on fracture site.

| available as aluminium/stainless steel/peek.

## Wrist Fixator



**Half-Pin (Schanz Screw) to be used**  
Ø (2,5/3,0/3,5 mm)

**Allen Wrench to be used**  
Ø (4,0 mm)

| CAT. NO.         | LENGTH (mm) |
|------------------|-------------|
| A600 12 110 0011 | 210         |




**Intended use :** Temporary fixation of fractures, nonunions, malunions, osteotomies, pseduarthrosis, of distal radius, arthrodesis of wrist joint. Especially used for fixation of open and/or communited intra-articular fractures of distal radius and for distraction and compression purposes on fracture site.

## Penning Type Wrist Fixator



**Allen Wrench to be used**  
Ø (3,0/4,0 mm)

**Half-Pin (Schanz Screw) to be used**  
Ø (2,5/3,0/3,5 mm)



**Intended use :** Temporary fixation of fractures, nonunions, malunions, osteotomies, pseduarthrosis, of distal radius, arthrodesis of wrist joint. Especially used for fixation of open and/or communited intra-articular fractures of distal radius and for distraction and compression purposes on fracture site.

| CAT. NO          | TYPE  | LENGTH (mm) |
|------------------|-------|-------------|
| A600 12 110 0012 | SHORT | 190         |
| A600 12 110 0013 | LONG  | 220         |


| available as aluminium/stainless steel/peek.

# Ayfix Dynamic Axial External Fixators

## Ayfix Pelvis Fixator



| CAT. NO          | TYPE     | DISTRACTION (mm) | LENGTH (mm) |
|------------------|----------|------------------|-------------|
| A600 25 000 0000 | LONG     | 60               | 373         |
| A600 25 000 0001 | STANDARD | 40               | 333         |




**Intended use :** Temporary fixation of fractures, osteotomies of pelvis in adults. Especially used for fixation of open and/or communitated fractures of pelvis.

| available as aluminium/stainless steel/peek.

## Ayfix Pediatric Pelvis Fixator



| CAT. NO          | DISTRACTION (mm) | LENGTH (mm) |
|------------------|------------------|-------------|
| A600 25 000 0110 | 25               | 275         |



**Intended use :** Temporary fixation of fractures, osteotomies of pelvis in adolescents.. Especially used for fixation of open and/or communitated fractures of pelvis.

| available as aluminium/stainless steel/peek.

# Ayfix Dynamic Axial External Fixators


## Dynamic Axial Knee Fixator



Allen Wrench to be used  
Ø (6,0 mm)

Half-Pin (Schanz Screw) to be used  
Ø (5,0/5,5/6,0 mm)

| CAT. NO          | DISTRACTION (mm) | LENGTH (mm) |
|------------------|------------------|-------------|
| A600 15 100 0002 | 71,5             | 350         |



**Intended use :** Temporary fixation of fractures, nonunions, malunions, osteotomies, pseduarthrosis, of distal femur and proximal tibia, arthrodesis of knee joint. Especially used for fixation of open and communited fractures of distal femur and proximal tibia and for distraction and compression purposes on knee joint.

| available as aluminium/stainless steel/peek.


## Dynamic Axial Elbow Fixator



Allen Wrench to be used  
Ø (5,0 mm)

Half-Pin (Schanz Screw) to be used  
Ø (3,5/4,0/4,5 mm)

| CAT. NO          | DISTRACTION (mm) | LENGTH (mm) |
|------------------|------------------|-------------|
| A600 15 100 0001 | 59               | 298         |



**Intended use :** Temporary fixation of fractures, nonunions, malunions, osteotomies, pseduarthrosis, of distal humerus and proximal radius/ulna, arthrodesis of elbow joint. Especially used for fixation of open and communited fractures of distal humerus and proximal radius/ulna and for distraction and compression purposes on elbow joint.


| available as aluminium/stainless steel/peek.

# Ayfix Dynamic Axial External Fixators

## Ayfix Limb Reconstruction Fixator Femur-Tibia




| CAT. NO          | TYPE   | DISTRACTION (mm) | LENGTH (mm) |
|------------------|--------|------------------|-------------|
| A600 11 100 0010 | Long   | 140              | 400         |
| A600 11 100 0013 | Medium | 90               | 300         |
| A600 11 100 0012 | Short  | 40               | 250         |




| available as aluminium/stainless steel/peek.

**Intended use :** Temporary fixation of fractures, nonunions, malunions, osteotomies, pseduarthrosis, of femur and tibia, arthrodesis of knee and ankle joints. Especially used for the treatment of reconstruction of femoral and tibial bone defects by sliding distal or proximal osteotomized bone with elongation of fibrous callus at osteotomy site. (new bone formation/ reconstruction of femoral and tibial bone defects. For adult patients.

## Ayfix Limb Reconstruction Pediatric Fixator



| CAT. NO          | TYPE   | DISTRACTION (mm) | LENGTH (mm) |
|------------------|--------|------------------|-------------|
| A600 11 101 0010 | Long   | 90               | 250         |
| A600 11 101 0013 | Medium | 90               | 200         |
| A600 11 101 0012 | Short  | 45               | 150         |




| available as aluminium/stainless steel/peek.


**Intended use :** Temporary fixation of fractures, nonunions, malunions, osteotomies, pseduarthrosis, of femur and tibia, arthrodesis of knee and ankle joints. Especially used for the treatment of reconstruction of femoral and tibial bone defects by sliding distal or proximal osteotomized bone with elongation of fibrous callus at osteotomy site. (new bone formation/ reconstruction of femoral and tibial bone defects. For pediatric patients.

# Ayfix Dynamic Axial External Fixators

## Polyaxial Femur-Tibia External Fixator




| CAT. NO          | LENGTH (mm) |
|------------------|-------------|
| A600 11 115 0000 | 245         |




**Intended use :** Temporary fixation of fractures, nonunions, malunions, osteotomies, pseduarthrosis, of femur/tibia. Especially used for fixation of open fractures of femur/ tibia and for distraction and compression purposes on fracture site.

| available as aluminium/stainless steel/peek.

## Polyaxial Humerus External Fixator



| CAT. NO          | LENGTH (mm) |
|------------------|-------------|
| A600 11 120 0000 | 170         |




**Intended use :** Temporary fixation of fractures, nonunions, malunions, osteotomies, pseduarthrosis, of humerus,Especially used for fixation of open fractures of humerus and for distraction and compression purposes on fracture site.

| available as aluminium/stainless steel/peek.

# Ayfix Dynamic Axial External Fixators


## Ayfix Pertrochanteric Fixator



Allen Wrench to be used  
Ø (6,0 mm)


| CAT. NO.         | LENGTH (mm) |
|------------------|-------------|
| A600 12 120 0003 | 117,4       |

Half-Pin (Schanz Screw) to be used  
Ø (5,0/5,5/6,0 mm)



**Intended use :** Temporary fixation of fractures, nonunions, malunions, osteotomies, pseduarthrosis, of proximal femur. Especially used for fixation of open fractures of proximal femur. Especially used for fixation of fractures of proximal femur(intertrochanteric, pertrochanteric) in elderly patients for early mobilization, for distraction and compression purposes on fracture site.

## Hybrid External Fixator Femur-Tibia




Allen Wrench to be used  
Ø (6,0 mm)

| CAT. NO.         | LENGTH (mm) |
|------------------|-------------|
| A600 11 014 0005 | 230,5       |

Rings to be used : Ø (160 mm)


Half-Pin (Schanz Screw) to be used  
: Ø (4,5/5,0/5,5/6,0 mm)




**Intended use :** Temporary fixation of fractures, nonunions, malunions, osteotomies, pseduarthrosis, of proximal femur. Especially used for fixation of open fractures of proximal femur. Especially used for fixation of fractures of proximal femur(intertrochanteric, pertrochanteric) in elderly patients for early mobilization. and for distraction and compression purposes on fracture site.

# Ayfix Dynamic Axial External Fixators

## Hybrid External Fixator Wrist



| CAT. NO.         | LENGTH (mm) |
|------------------|-------------|
| A600 25 001 0025 | 174,3       |




**Intended use :** Temporary fixation of fractures, nonunions, malunions, osteotomies, pseduarthrosis, of distal radius and ulna, arthrodesis of wrist joint. Especially used for fixation of open fractures of distal radius and ulna and for distraction and compression purposes on fracture site. This fixator provides multiplanar and stronger fixation when needed.

## Hybrid External Fixator Finger




| CAT. NO.         | LENGTH (mm) |
|------------------|-------------|
| A600 25 001 0026 | 57,4        |



**Intended use :** Temporary fixation of fractures, nonunions, malunions, osteotomies, pseduarthrosis, of metacarpal, metatarsal and phalangeal bones, arthrodesis of metacarpophalangeal and interphalangeal joints. Especially used for fixation of open fractures of for distraction and compression purposes on fracture site and joints. This fixator provides multiplanar and stronger fixation when needed.

## Ayfix Finger Fixator



| CAT. NO.         | LENGTH (mm) |
|------------------|-------------|
| A600 12 120 0001 | 72,7        |




**Intended use :** Temporary fixation of fractures, nonunions, malunions, osteotomies, pseduarthrosis, of metacarpal, metatarsal, phalangeal bones. Arthrodesis of metacarpophalangeal and interphalangeal joints. Especially used for fixation of open fracture, for lengthening purposes of metacarpal metatarsal, phalangeal bones.

# Mini Rail External Fixators


## Mini Rail Fixator Short

Allen Wrench to be used  
Ø (3,0 mm)



Half-Pin (Schanz Screw) to be used  
Ø (1,5/2,0/2,5 mm)

| CAT. NO          | DISTRACTION (mm) | LENGTH (mm) |
|------------------|------------------|-------------|
| A600 12 120 0005 | 19,5             | 71,5        |




**Intended use :** Temporary fixation of fractures, nonunions, malunions, osteotomies, pseduarthrosis, of metacarpal & metatarsal bones. Arthrodesis of metametacarpophalangeal & metatarsophalangeal joints. Especially used for fixation of open fracture, for lengthening purposes of metacarpal & metatarsal bones.

| available as aluminium/stainless steel/peek.


## Mini Rail Fixator Medium

Allen Wrench to be used  
Ø (3,0 mm)



Half-Pin (Schanz Screw) to be used  
Ø (1,5/2,0/2,5 mm)

| CAT. NO          | DISTRACTION (mm) | LENGTH (mm) |
|------------------|------------------|-------------|
| A600 12 120 0006 | 29,5             | 81,5        |




**Intended use :** Temporary fixation of fractures, nonunions, malunions, osteotomies, pseduarthrosis, of metacarpal & metatarsal bones. Arthrodesis of metametacarpophalangeal & metatarsophalangeal joints. Especially used for fixation of open fracture, for lengthening purposes of metacarpal & metatarsal bones.

| available as aluminium/stainless steel/peek.

## Mini Rail Fixator Long

Allen Wrench to be used  
Ø (3,0 mm)



Half-Pin (Schanz Screw) to be used  
Ø (1,5/2,0/2,5 mm)

| CAT. NO          | DISTRACTION (mm) | LENGTH (mm) |
|------------------|------------------|-------------|
| A600 12 120 0007 | 39,5             | 91,5        |



**Intended use :** Temporary fixation of fractures, nonunions, malunions, osteotomies, pseduarthrosis, of metacarpal, metatarsal, occasionally ulna, radius bones. Arthrodesis of metacarpophalangeal, matatarsophalangeal & occasionally wrist joints. Especially used for fixation of open fracture, for lengthening purposes of matacarpal, metatarsal, occasionally ulna, radius bones.

| available as aluminium/stainless steel/peek.


## Mini Rail Fixator X Long

Allen Wrench to be used  
Ø (3,0 mm)



Half-Pin (Schanz Screw) to be used  
Ø (1,5/2,0/2,5 mm)

| CAT. NO          | DISTRACTION (mm) | LENGTH (mm) |
|------------------|------------------|-------------|
| A600 12 120 0002 | 69,5             | 121,5       |




**Intended use :** Temporary fixation of fractures, nonunions, malunions, osteotomies, pseduarthrosis, of metacarpal, metatarsal, occasionally ulna, radius bones. Arthrodesis of metacarpophalangeal, matatarsophalangeal & occasionally wrist joints. Especially used for fixation of open fracture, for lengthening purposes of matacarpal, metatarsal, occasionally ulna, radius bones.

| available as aluminium/stainless steel/peek.

# Mini Rail External Fixators


## Mini Rail Fixator XX Long



**Intended use :** Temporary fixation of fractures, nonunions, malunions, osteotomies, pseduarthrosis, of metacarpal, metatarsal, occasionally ulna, radius bones. Arthrodesis of metacarpophalangeal, matatarsophalangeal & occasionally wrist joints. Especially used for fixation of open fracture, for lengthening purposes of matacarpal, metatarsal, occasionally ulna, radius bones.

| available as aluminium/stainless steel/peek.


## Mini Rail Fixator Vertical Axis



**Intended use :** Temporary fixation of fractures, nonunions, malunions, osteotomies, pseduarthrosis, of metacarpal & metatarsal bones. Arthrodesis of metametacarpophalangeal & metatarsophalangeal joints. Especially used for fixation of open fracture, for lengthening purposes of metacarpal & metatarsal bones.

| available as aluminium/stainless steel/peek.

## Mini Rail Fixator Horizontal Axis



**Intended use :** Temporary fixation of fractures, nonunions, malunions, osteotomies, pseduarthrosis, of metacarpal & metatarsal bones. Arthrodesis of metametacarpophalangeal & metatarsophalangeal joints. Especially used for fixation of open fracture, for lengthening purposes of metacarpal & metatarsal bones.

| available as aluminium/stainless steel/peek.

## Mini Rail Fixator Horizontal Axis 90°



**Intended use :** Temporary fixation of fractures, nonunions, malunions, osteotomies, pseduarthrosis, of metacarpal & metatarsal bones. Arthrodesis of metametacarpophalangeal & metatarsophalangeal joints. Especially used for fixation of open fracture, for lengthening purposes of metacarpal & metatarsal bones.

| available as aluminium/stainless steel/peek.

# Tubular External Fixators

## Tubular Fixator Rod-Rod Clamp



| ALUMINIUM        | ROD DIAMETER (mm) |
|------------------|-------------------|
| A600 11 006 0001 | 5,0               |
| A600 11 014 0001 | 8,0               |
| A600 11 010 0001 | 10,0              |

| available as aluminium/stainless steel/peek.

## Tubular Fixator Rod-Pin Clamp



| ALUMINIUM        | ROD DIAMETER (mm) | SCREW DIAMETER (mm) |
|------------------|-------------------|---------------------|
| A600 11 006 0002 | 5,0               | 2,5-3,0-3,5         |
| A600 11 014 0002 | 8,0               | 3,5-4,0-4,5         |
| A600 11 010 0002 | 10,0              | 5,0-5,5-6,0         |

| available as aluminium/stainless steel/peek.

## Tubular External Fixators

### Tubular Fixator Carbon Rod



| CARBON           | DIAMETER x LENGTH<br>(mm) |
|------------------|---------------------------|
| A600 11 005 0050 | 5,0X50                    |
| A600 11 005 0100 | 5,0X100                   |
| A600 11 005 0150 | 5,0X150                   |
| A600 11 005 0200 | 5,0X200                   |
| A600 11 005 0250 | 5,0X250                   |
| A600 11 005 0300 | 5,0X300                   |
| A600 11 005 0350 | 5,0X350                   |
| A600 11 006 0080 | 6,0X80                    |
| A600 11 006 0100 | 6,0X100                   |
| A600 11 006 0150 | 6,0X150                   |
| A600 11 006 0200 | 6,0X200                   |
| A600 11 006 0250 | 6,0X250                   |
| A600 11 006 0300 | 6,0X300                   |
| A600 11 006 0350 | 6,0X350                   |
| A600 11 008 0050 | 8,0X50                    |
| A600 11 008 0100 | 8,0X100                   |
| A600 11 008 0150 | 8,0X150                   |

| CARBON           | DIAMETER x LENGTH<br>(mm) |
|------------------|---------------------------|
| A600 11 008 0200 | 8,0X200                   |
| A600 11 008 0250 | 8,0X250                   |
| A600 11 008 0300 | 8,0X300                   |
| A600 11 008 0350 | 8,0X350                   |
| A600 11 010 0100 | 10,0X100                  |
| A600 11 010 0150 | 10,0X150                  |
| A600 11 010 0200 | 10,0X200                  |
| A600 11 010 0250 | 10,0X250                  |
| A600 11 010 0300 | 10,0X300                  |
| A600 11 010 0350 | 10,0X350                  |
| A600 11 014 0100 | 14,0X100                  |
| A600 11 014 0150 | 14,0X150                  |
| A600 11 014 0200 | 14,0X200                  |
| A600 11 014 0250 | 14,0X250                  |
| A600 11 014 0300 | 14,0X300                  |
| A600 11 014 0350 | 14,0X350                  |

### Tubular Fixator Multiple Pin Clamp



| ALUMINIUM        | DIAMETER<br>(mm) |
|------------------|------------------|
| A600 11 006 0003 | 5,0              |
| A600 11 014 0003 | 8,0              |
| A600 11 010 0003 | 10,0             |

| available as aluminium/stainless steel/peek.

# Tubular External Fixators

## Tubular Fixator Straight Support



| SS               | DIAMETER (mm) |
|------------------|---------------|
| A600 11 005 0004 | 5,0           |
| A600 11 008 0004 | 8,0           |
| A600 11 010 0007 | 10,0          |

## Tubular Fixator 30° Angled Support



| SS               | DIAMETER (mm) |
|------------------|---------------|
| A600 29 121 0120 | 5,0           |
| A600 29 121 0080 | 8,0           |
| A600 29 121 0100 | 10,0          |

## Tubular Fixator Multiple Pin T Clamp



| ALUMINIUM        | DIAMETER (mm) |
|------------------|---------------|
| A600 11 006 0004 | 5,0           |
| A600 11 014 0004 | 8,0           |
| A600 11 010 0004 | 10,0          |

| available as aluminium/stainless steel/peek.

## Tubular Fixator Semi Circular Rod



| SS               | DIAMETER x LENGTH (mm) |
|------------------|------------------------|
| A600 10 151 0160 | 8,0X160                |
| A600 10 151 0180 | 8,0X180                |
| A600 10 151 0200 | 8,0X200                |
| A600 10 121 0160 | 10,0X160               |
| A600 10 121 0180 | 10,0X180               |
| A600 10 121 0200 | 10,0X200               |

## Tubular Fixator Compression Distraction Tube



| ALUMINIUM        | DIAMETER x LENGTH (mm) |
|------------------|------------------------|
| A600 11 007 0022 | 8,0X200                |

\* DIA. = DIAMETER (mm) \*\* LENGTH = LENGTH (mm)

## Tubular External Fixators

### Tubular Fixator Ankle Hinge



| SS               | DIAMETER x LENGTH<br>(mm) |
|------------------|---------------------------|
| A600 11 007 0020 | 8,0x200                   |

### Tubular Fixator Elbow-Wrist Hinge



| SS               | DIAMETER x LENGTH<br>(mm) |
|------------------|---------------------------|
| A600 11 007 0019 | 5,0x200                   |



# Hybrid External Fixators

## Hybrid Fixator Aluminium 2/3 Ring



| CAT. NO          | DIAMETER (mm) |
|------------------|---------------|
| A600 11 007 0014 | 155           |
| A600 11 007 0016 | 185           |
| A600 11 007 0018 | 215           |

available as aluminium/stainless steel/peek.

## Hybrid Fixator Aluminium 1/3 Ring



| CAT. NO          | DIAMETER (mm) |
|------------------|---------------|
| A600 11 008 0014 | 155           |
| A600 11 008 0016 | 185           |
| A600 11 008 0018 | 215           |

available as aluminium/stainless steel/peek.

## Hybrid Fixator Three Hole Sliding Wire Clamp



| CAT. NO.         |
|------------------|
| A600 11 008 0026 |

available as aluminium/stainless steel/peek.

## Washer



| SS               | DIAMETER (mm) |
|------------------|---------------|
| A600 10 030 0012 | 12            |
| A600 10 030 0013 | 14            |

## Hybrid Fixator Two Hole Wire Securing Pin



| SS               |
|------------------|
| A600 10 153 0001 |

# Hybrid External Fixators

## Olive Wire



| CAT. NO          | DIA. x LENGTH (mm) |
|------------------|--------------------|
| A600 10 188 2030 | 2,0X300            |
| A600 10 188 2040 | 2,0X400            |

## Trocar Pointed Wire



| SS               | DIAMETER (mm) | LENGTH (mm) |
|------------------|---------------|-------------|
| A600 10 100 2030 | 2,0           | 300         |
| A600 10 100 2040 | 2,0           | 400         |

## Hybrid Fixator Independent Wire Clamp



| SS               |
|------------------|
| A600 11 008 0024 |

## Hybrid Fixator Supplementary Pin Holder Clamp



| CAT. NO          |
|------------------|
| A600 11 008 0025 |

| available as aluminium/stainless steel/peek.

## Hybrid Fixator Carbon Rod



| CAT. NO          |
|------------------|
| A600 11 006 0060 |

# Hybrid External Fixators

## Wire Fixation Bolt



| SS               |
|------------------|
| A600 10 090 0000 |

## Hybrid Fixator Straight Posts



| SS               |
|------------------|
| A600 11 008 0021 |

## Hybrid Fixator 30° Angled Posts



| SS               |
|------------------|
| A600 11 008 0022 |

## Half Pin Universal Retaining Bolt



| SS               |
|------------------|
| A600 10 099 0001 |

## Connection Bolt (M6)



| SS               | LENGTH (mm) |
|------------------|-------------|
| A600 10 020 0025 | 25          |
| A600 10 020 0030 | 30          |

## Spiked Washer



| SS               |
|------------------|
| A600 10 032 0000 |

## Nut M6



| SS               |
|------------------|
| A600 10 021 0000 |

## Hybrid Fixator Single Pin Clamp



| SS               |
|------------------|
| A600 11 007 0012 |

## Hybrid Fixator Hinge Clamp



| available as aluminium/stainless steel/peek.

| CAT. NO          |
|------------------|
| A600 11 008 0019 |

\* DIA. = DIAMETER (mm) \*\* LENGTH = LENGTH (mm)

# Hybrid External Fixators

## Hybrid Fixator Kit



CAT. NO

A600 11 007 0023

| available as aluminium/stainless steel/peek.

## Hybrid Fixator Translation Rotation Hinge



CAT. NO

A600 11 008 0020

| available as aluminium/stainless steel/peek.

## Hybrid Fixator Telescopic Reduction Unit



CAT. NO

A600 11 007 0021

| available as aluminium/stainless steel/peek.

## Square Socket



| SS               | HOLES |
|------------------|-------|
| A600 14 002 0001 | 1     |
| A600 14 002 0002 | 2     |
| A600 14 002 0003 | 3     |
| A600 14 002 0004 | 4     |

| available as aluminium/stainless steel/peek.

## Hybrid Fixator Extension Plate




SS

A600 11 050 2251

# Hybrid External Fixators

## Threaded Rod (M6)



| SS               | LENGTH (mm) |
|------------------|-------------|
| A600 10 050 0050 | 50          |
| A600 10 050 0075 | 75          |
| A600 10 050 0100 | 100         |
| A600 10 050 0115 | 115         |
| A600 10 050 0125 | 125         |
| A600 10 050 0150 | 150         |
| A600 10 050 0175 | 175         |
| A600 10 050 0200 | 200         |
| A600 10 050 0250 | 250         |
| A600 10 050 0300 | 300         |



Circular External Fixators

Half Ring



| ALUMINIUM        | DIAMETER (mm) |
|------------------|---------------|
| A600 10 120 0080 | 80            |
| A600 10 120 0100 | 100           |
| A600 10 120 0120 | 120           |
| A600 10 120 0130 | 130           |
| A600 10 120 0140 | 140           |
| A600 10 120 0150 | 150           |
| A600 10 120 0160 | 160           |
| A600 10 120 0170 | 170           |
| A600 10 120 0180 | 180           |
| A600 10 120 0200 | 200           |
| A600 10 120 0220 | 220           |

5/8 C Ring



| ALUMINIUM        | DIAMETER (mm) |
|------------------|---------------|
| A600 10 150 0080 | 80            |
| A600 10 150 0100 | 100           |
| A600 10 150 0120 | 120           |
| A600 10 150 0140 | 140           |
| A600 10 150 0150 | 150           |
| A600 10 150 0160 | 160           |
| A600 10 150 0170 | 170           |
| A600 10 150 0180 | 180           |
| A600 10 150 0200 | 200           |
| A600 10 150 0220 | 220           |

3/4 Foot Ring



| ALUMINIUM        | DIAMETER (mm) |
|------------------|---------------|
| A600 29 120 0220 | 80            |
| A600 29 120 0100 | 100           |
| A600 29 120 0120 | 120           |
| A600 29 120 0140 | 140           |
| A600 29 120 0160 | 160           |

## Circular External Fixators

### Femur Arch



| ALUMINIUM        | SS               | ANGLE/DIAMETER (mm) |
|------------------|------------------|---------------------|
| A601 10 160 0090 | A601 11 160 0090 | 90° / Ø 160         |
| A601 10 180 0090 | A601 11 180 0090 | 90° / Ø 180         |
| A601 10 200 0090 | A601 11 200 0090 | 90° / Ø 200         |
| A601 10 160 0120 | A601 11 160 0120 | 120° / Ø 160        |
| A601 10 180 0120 | A601 11 180 0120 | 120° / Ø 180        |
| A601 10 200 0120 | A601 11 200 0120 | 120° / Ø 200        |

### Multiple Pin Fixation Clamp



| SS               |
|------------------|
| A601 10 013 0000 |
| A601 10 013 0001 |
| A601 10 013 0002 |
| A601 10 013 0003 |

### Oblique Support Connector



| SS               | LENGTH (mm) |
|------------------|-------------|
| A600 10 080 0000 | 30          |
| A600 10 080 0001 | 58          |
| A600 10 080 0002 | 82          |

### Adjustable Clamp



| SS               |
|------------------|
| A600 10 001 0001 |

### Telescopic Rod



| SS               | LENGTH (mm) |
|------------------|-------------|
| A600 10 071 0100 | 100         |
| A600 10 071 0125 | 125         |
| A600 10 071 0150 | 150         |
| A600 10 071 0175 | 175         |
| A600 10 071 0215 | 215         |

### Double Female Hinge



| SS               |
|------------------|
| A601 10 099 0010 |

# Circular External Fixators

## Male Half Hinge



SS

A601 10 099 0015

## Female Half Hinge



SS

A601 10 099 0020

## Cannulated Bolt



SS

A600 10 090 0001

## Slotted Bolt



SS

A600 10 091 0000

## Half Pin Universal Retaining Bolt



SS

A600 10 099 0001

## Single Pin Fixation Bolt



SS

A601 10 010 0000

## Female Support



| SS               | HOLES | LENGTH (mm) |
|------------------|-------|-------------|
| A600 10 065 0000 | 1     | 23          |
| A600 10 065 0002 | 2     | 32          |
| A600 10 065 0003 | 3     | 41          |


## Male Support



| SS               | HOLES | LENGTH (mm) |
|------------------|-------|-------------|
| A600 10 060 0112 | 1     | 33          |
| A600 10 060 0222 | 2     | 42          |
| A600 10 060 0332 | 3     | 51          |

# Circular External Fixators

## Connection Bolt (M6)



| SS               | LENGTH (mm) |
|------------------|-------------|
| A600 10 020 0016 | 16          |
| A600 10 020 0020 | 20          |
| A600 10 020 0025 | 25          |
| A600 10 020 0030 | 30          |

## Connecting Plate



| SS               | HOLES | LENGTH (mm) |
|------------------|-------|-------------|
| A600 10 050 2035 | 2     | 41          |
| A600 10 050 3050 | 3     | 51          |
| A600 10 050 4060 | 4     | 61          |
| A600 10 050 5070 | 5     | 71          |
| A600 10 050 6085 | 6     | 81          |
| A600 10 050 7100 | 7     | 91          |
| A600 10 050 8110 | 8     | 101         |

## Connecting Plate with Threaded End



| SS               | HOLES | LENGTH (mm) |
|------------------|-------|-------------|
| A600 11 050 2131 | 5     | 131         |
| A600 11 050 2171 | 7     | 171         |
| A600 11 050 2211 | 9     | 211         |

## 90 ° Angled Connection Plate



| SS               | HOLES | LENGTH (mm) |
|------------------|-------|-------------|
| A600 14 000 1150 | 2     | 50          |
| A600 14 000 1175 | 4     | 75          |

## Twisted Plate



| SS               | HOLES | LENGTH (mm) |
|------------------|-------|-------------|
| A600 10 054 0046 | 2     | 50          |
| A600 10 054 0067 | 3     | 70          |
| A600 10 054 0088 | 4     | 90          |

## Universal Hinge



| SS               |
|------------------|
| A601 15 000 0000 |

## Curved Plate



| ALUMINIUM        |
|------------------|
| A601 15 000 0020 |

## Circular External Fixators

### Square Socket



| SS               | HOLES |
|------------------|-------|
| A600 14 002 0001 | 1     |
| A600 14 002 0002 | 2     |
| A600 14 002 0003 | 3     |
| A600 14 002 0004 | 4     |

### Square Socket Centering Sleeve



| SS               | DIAMETER (mm) |
|------------------|---------------|
| A600 14 002 0605 | 6             |

### Washer



| SS               | DIAMETER (mm) |
|------------------|---------------|
| A600 10 030 0012 | 12            |
| A600 10 030 0013 | 14            |

### Slotted Washer



| SS               |
|------------------|
| A600 10 034 0000 |

### Oval Washer



| SS               | DIAMETER (mm) |
|------------------|---------------|
| A600 10 031 0012 | 12            |

### Conical Washer



| SS               |
|------------------|
| A600 10 110 0001 |

### Star Washer



| SS               |
|------------------|
| A600 10 032 0000 |

### Nut M6



| SS               |
|------------------|
| A600 10 021 0000 |

# Circular External Fixators

## Threaded Rod (M6)



| SS               | LENGTH (mm) |
|------------------|-------------|
| A600 10 050 0050 | 50          |
| A600 10 050 0075 | 75          |
| A600 10 050 0100 | 100         |
| A600 10 050 0115 | 115         |
| A600 10 050 0125 | 125         |
| A600 10 050 0150 | 150         |
| A600 10 050 0175 | 175         |
| A600 10 050 0200 | 200         |
| A600 10 050 0250 | 250         |
| A600 10 050 0300 | 300         |
| A600 10 050 0350 | 350         |
| A600 10 050 0400 | 400         |

## Lengthening Rod (M6)



| CAT. NO          | LENGTH (mm) |
|------------------|-------------|
| A600 10 011 0060 | 60          |
| A600 10 011 0080 | 80          |
| A600 10 011 0100 | 100         |
| A600 10 011 0120 | 120         |
| A600 10 011 0150 | 150         |
| A600 10 011 0200 | 200         |

## Slotted Rod (M6)



| SS               | LENGTH (mm) |
|------------------|-------------|
| A600 08 011 0040 | 40          |
| A600 08 011 0050 | 50          |
| A600 08 011 0060 | 60          |
| A600 08 011 0070 | 70          |
| A600 08 011 0100 | 100         |

## Trocar Pointed Wire



| SS               | DIA x LENGTH (mm) |
|------------------|-------------------|
| A600 10 100 1830 | 1,8X300           |
| A600 10 100 1840 | 1,8X400           |
| A600 10 100 2030 | 2,0X300           |
| A600 10 100 2040 | 2,0X400           |

## Circular External Fixators

### Bayonet Pointed Wire



| SS               | DIA x LENGTH (mm) |
|------------------|-------------------|
| A600 10 139 1830 | 1,8X300           |
| A600 10 139 1840 | 1,8X400           |
| A600 10 139 2030 | 2,0X300           |
| A600 10 139 2040 | 2,0X400           |

### Olive Wire



| SS               | DIA x LENGTH (mm) |
|------------------|-------------------|
| A600 10 188 1830 | 1,8X300           |
| A600 10 188 1840 | 1,8X400           |
| A600 10 188 2030 | 2,0X300           |
| A600 10 188 2040 | 2,0X400           |

### Hybrid External Fixator (For adaptation to External Circular Fixation System)



| CAT.NO           | TYPE   | LENGTH (mm) |
|------------------|--------|-------------|
| A600 25 001 0025 | Medium | 230         |
| A600 11 014 0005 | Long   | 250         |

## Pediatric Circular External Fixators

### Pediatric Half Ring



| ALUMINIUM        | DIAMETER (mm) |
|------------------|---------------|
| A600 11 121 0070 | 70            |
| A600 11 121 0080 | 80            |
| A600 11 121 0090 | 90            |
| A600 11 121 0100 | 100           |
| A600 11 121 0110 | 110           |
| A600 11 121 0120 | 120           |

### Pediatric 5/8 C Ring



| ALUMINIUM        | DIAMETER (mm) |
|------------------|---------------|
| A600 11 151 0070 | 70            |
| A600 11 151 0080 | 80            |
| A600 11 151 0090 | 90            |
| A600 11 151 0100 | 100           |
| A600 11 151 0110 | 110           |
| A600 11 151 0120 | 120           |

### Pediatric Femur Arch



| ALUMINIUM        | ANGLE | DIAMETER (mm) |
|------------------|-------|---------------|
| A601 11 090 0090 | 90°   | 90            |
| A601 11 110 0090 | 90°   | 110           |
| A601 11 090 0120 | 120°  | 90            |
| A601 11 110 0120 | 120°  | 110           |

### Pediatric Adjustable Clamp



**SS**  
A600 11 002 0002

### Pediatric Multiple Pin Fixation Clamp



**SS**  
A601 11 014 0000

Pediatric Circular External Fixators

Pediatric Slotted Bolt



| SS               |
|------------------|
| A600 11 091 0000 |

Pediatric Oblique Support Connector



| SS               | LENGTH (mm) |
|------------------|-------------|
| A600 11 080 0000 | 25          |
| A600 11 080 0001 | 50          |
| A600 11 080 0002 | 75          |

Pediatric Cannulated Bolt



| SS               |
|------------------|
| A600 11 090 0000 |

Pediatric Telescopic Rod



| SS               |
|------------------|
| A600 11 072 0001 |

Pediatric Half Pin Universal Retaining Bolt



| SS               |
|------------------|
| A600 11 099 0001 |

Pediatric Male Half Hinge



| SS               |
|------------------|
| A601 11 099 0015 |

Pediatric Female Half Hinge



| SS               |
|------------------|
| A601 11 099 0020 |

Pediatric Lengthening Rod



| CAT. NO.         |
|------------------|
| A601 11 001 0101 |

## Pediatric Circular External Fixators

### Pediatric Double Female Hinge



**SS**

A601 11 099 0010

### Pediatric Connecting Plate



| SS               | HOLES | LENGTH (mm) |
|------------------|-------|-------------|
| A600 11 050 2035 | 2     | 25          |
| A600 11 050 3050 | 3     | 45          |
| A600 11 050 4060 | 4     | 55          |
| A600 11 050 5070 | 5     | 65          |

### Pediatric Universal Hinge



**SS**

A601 16 000 0000

### Pediatric Washer



**SS**

A600 11 030 0101

### Pediatric Star Washer



**SS**

A600 11 032 0000

### Pediatric Curved Plate



**ALUMINIUM**

A601 16 000 0020

### Pediatric Nut (M5)



**SS**

A600 11 021 0000

### Pediatric Slotted Washer



**SS**

A600 11 034 0000

### Pediatric Oval Washer



**SS**

A600 11 031 0101

## Pediatric Circular External Fixators

### Pediatric Square Socket



| SS               | HOLES |
|------------------|-------|
| A600 16 002 0001 | 1     |
| A600 16 002 0002 | 2     |
| A600 16 002 0003 | 3     |
| A600 16 002 0004 | 4     |

### Pediatric Square Socket Centering Sleeve



| SS               | DIAMETER (mm) |
|------------------|---------------|
| A600 14 002 0305 | 3             |
| A600 14 002 0405 | 4             |
| A600 14 002 0505 | 5             |

### Pediatric Male Support



| SS               | HOLES | LENGTH (mm) |
|------------------|-------|-------------|
| A600 11 060 0112 | 1     | 31,5        |
| A600 11 060 0222 | 2     | 40,5        |
| A600 11 060 0332 | 3     | 49,5        |

### Pediatric Connection Bolt (M5)



| SS               | LENGTH (mm) |
|------------------|-------------|
| A600 11 020 0012 | 12          |
| A600 11 020 0016 | 16          |

### Pediatric 90° Angled Connection Plate



| SS               | HOLES | LENGTH (mm) |
|------------------|-------|-------------|
| A600 14 001 1150 | 2     | 41          |
| A600 14 001 1175 | 4     | 65          |

### Twisted Plate



| SS               | HOLES | LENGTH (mm) |
|------------------|-------|-------------|
| A600 11 054 0046 | 2     | 45          |
| A600 11 054 0067 | 3     | 63          |
| A600 11 054 0088 | 4     | 81          |

### Pediatric Female Support



| SS               | HOLES | LENGTH (mm) |
|------------------|-------|-------------|
| A600 11 065 0000 | 1     | 22,5        |
| A600 11 065 0002 | 2     | 31,5        |
| A600 11 065 0003 | 3     | 40,5        |

### Pediatric Threaded Socket



| SS               | LENGTH (mm) |
|------------------|-------------|
| A600 16 001 1115 | 15          |
| A600 16 001 1125 | 25          |

## Pediatric Circular External Fixators

### Pediatric Trocar Pointed Wire



| SS               | DIA x LENGTH<br>(mm) |
|------------------|----------------------|
| A600 11 100 1520 | 1,5X200              |
| A600 11 100 1525 | 1,5X250              |
| A600 11 100 1820 | 1,8X200              |
| A600 11 100 1825 | 1,8X250              |

### Pediatric Bayonet Pointed Wire




| SS               | DIA x LENGTH<br>(mm) |
|------------------|----------------------|
| A600 11 135 1520 | 1,5X200              |
| A600 11 135 1525 | 1,5X250              |
| A600 11 135 1820 | 1,8X200              |
| A600 11 135 1825 | 1,8X250              |

### Pediatric Olive Wire



| SS               | DIA x LENGTH<br>(mm) |
|------------------|----------------------|
| A600 11 188 1520 | 1,5X200              |
| A600 11 188 1525 | 1,5X250              |
| A600 11 188 1820 | 1,8X200              |
| A600 11 188 1825 | 1,8X250              |

### Pediatric Threaded Rod (M5)



| SS               | LENGTH<br>(mm) |
|------------------|----------------|
| A600 11 050 0040 | 40             |
| A600 11 050 0050 | 50             |
| A600 11 050 0060 | 60             |
| A600 11 050 0080 | 80             |
| A600 11 050 0090 | 90             |
| A600 11 050 0100 | 100            |
| A600 11 050 0120 | 120            |
| A600 11 050 0150 | 150            |
| A600 11 050 0180 | 180            |

"We supply WORLD's  
Orthopedic & Trauma  
Implants..."

# Wires, Pins, Staples, & Drills Bits



**aysam**

## Wires, Pins, Staples & Drill Bits

### Kirschner Wires with Trocar Point



| SS               | DIAMETER (mm) | SS               | DIAMETER (mm) |
|------------------|---------------|------------------|---------------|
| A100 03 070 0620 | 0,6           | A100 03 070 2028 | 2,0           |
| A100 03 070 0825 | 0,8           | A100 03 070 2225 | 2,2           |
| A100 03 070 0925 | 0,9           | A100 03 070 2510 | 2,5           |
| A100 03 070 1025 | 1,0           | A100 03 070 2810 | 2,8           |
| A100 03 070 1210 | 1,2           | A100 03 070 3010 | 3,0           |
| A100 03 070 1510 | 1,4           | A100 03 070 3225 | 3,2           |
| A100 03 070 1525 | 1,5           | A100 03 070 3510 | 3,5           |
| A100 03 070 1625 | 1,6           | A100 03 070 4010 | 4,0           |
| A100 03 070 1810 | 1,8           | A100 03 070 4529 | 4,5           |
|                  |               | A100 03 070 5025 | 5,0           |

### Kirschner Wires with Double Trocar Point



| SS               | DIAMETER (mm) | SS               | DIAMETER (mm) |
|------------------|---------------|------------------|---------------|
| A100 03 073 0625 | 0,6           | A100 03 073 2012 | 2,0           |
| A100 03 073 0812 | 0,8           | A100 03 073 2212 | 2,2           |
| A100 03 073 0925 | 0,9           | A100 03 073 2510 | 2,5           |
| A100 03 073 1012 | 1,0           | A100 03 073 3015 | 2,8           |
| A100 03 073 1212 | 1,2           | A100 03 073 3010 | 3,0           |
| A100 03 073 1429 | 1,4           | A100 03 073 3229 | 3,2           |
| A100 03 073 1512 | 1,5           | A100 03 073 3510 | 3,5           |
| A100 03 073 1629 | 1,6           | A100 03 073 4010 | 4,0           |
| A100 03 073 1812 | 1,8           | A100 03 073 6032 | 4,5           |
|                  |               | A100 03 073 5020 | 5,0           |

## Wires, Pins, Staples & Drill Bits

### Threaded Kirschner Wires with Trocar Point



| SS               | DIAMETER (mm) | SS               | DIAMETER (mm) |
|------------------|---------------|------------------|---------------|
| A100 03 071 0020 | 2,0           | A100 03 071 0032 | 3,2           |
| A100 03 071 0025 | 2,5           | A100 03 071 0035 | 3,5           |
| A100 03 071 0030 | 3,0           | A100 03 071 0040 | 4,0           |

### Cerclage Wire



| SS               | DIA x LENGTH (mm) | SS               | DIA x LENGTH (mm) |
|------------------|-------------------|------------------|-------------------|
| A400 03 011 0329 | 0,3X5000          | A400 03 011 1019 | 1,0X5000          |
| A400 03 011 0427 | 0,4X5000          | A400 03 011 1218 | 1,2X5000          |
| A400 03 011 0525 | 0,5X5000          | A400 03 011 1517 | 1,5X5000          |
| A400 03 011 0623 | 0,6X5000          | A400 03 011 1826 | 1,8X5000          |
| A400 03 011 0722 | 0,7X5000          | A400 03 011 2027 | 2,0X5000          |
| A400 03 011 0821 | 0,8X5000          | A400 03 011 2328 | 2,3X5000          |
| A400 03 011 0920 | 0,9X5000          | A400 03 011 2529 | 2,5X5000          |

### Steinmann Pin with Trocar Point



| SS               | DIAMETER (mm) | SS               | DIAMETER (mm) |
|------------------|---------------|------------------|---------------|
| A100 03 050 1829 | 1,8           | A100 03 050 4012 | 4,0           |
| A100 03 050 2029 | 2,0           | A100 03 050 4512 | 4,5           |
| A100 03 050 2529 | 2,5           | A100 03 050 5012 | 5,0           |
| A100 03 050 3012 | 3,0           | A100 03 050 6012 | 6,0           |
| A100 03 050 3512 | 3,5           |                  |               |

## Wires, Pins, Staples & Drill Bits

### Steinmann Pins with Central Thread



| SS               | DIAMETER (mm) | THREAD LENGTH (mm) | SS               | DIAMETER (mm) | THREAD LENGTH (mm) |
|------------------|---------------|--------------------|------------------|---------------|--------------------|
| A100 03 051 3230 | 3,0           | 60                 | A100 03 051 5230 | 5,0           | 60                 |
| A100 03 051 4230 | 4,0           | 60                 | A100 03 051 5524 | 5,5           | 60                 |
| A100 03 051 4523 | 4,5           | 60                 | A100 03 050 6018 | 6,0           | 60                 |

### Schanz Screw



| SS               | DIAMETER (mm) | SS               | DIAMETER (mm) |
|------------------|---------------|------------------|---------------|
| A100 03 090 5012 | 1,5           | A100 03 090 4011 | 4,0           |
| A100 03 090 6026 | 2,0           | A100 03 090 4517 | 4,5           |
| A100 03 091 6020 | 2,5           | A100 03 090 5013 | 5,0           |
| A100 03 090 3115 | 3,0           | A100 03 090 5515 | 5,5           |
| A100 03 090 6030 | 3,5           | A100 03 090 6012 | 6,0           |

### Schanz Screw Cortical



| SS               | DIA.x LENGTH (mm) |
|------------------|-------------------|
| A600 11 110 0130 | 6,0X130           |
| A600 11 110 0115 | 6,0X150           |
| A600 11 110 0160 | 6,0X160           |
| A600 11 110 0180 | 6,0X180           |
| A600 11 110 0200 | 6,0X200           |
| A600 11 110 0220 | 6,0X220           |

## Wires, Pins, Staples & Drill Bits

### Schanz Screw Cancellous



| SS               | DIA.x LENGTH<br>(mm) |
|------------------|----------------------|
| A600 11 111 0150 | 6,0X150              |
| A600 11 111 0160 | 6,0X160              |
| A600 11 111 0180 | 6,0X180              |
| A600 11 111 0200 | 6,0X200              |
| A600 11 111 0220 | 6,0X220              |

### Calibrated Self Tapping Schanz Screw



| SS               | DIAMETER<br>(mm) |
|------------------|------------------|
| A600 13 110 4018 | 4,0              |
| A600 13 110 5013 | 5,0              |
| A600 13 110 0160 | 6,0              |

### Rush Intramedullary Pin



| SS               | DIAMETER<br>(mm) | SS               | DIAMETER<br>(mm) |
|------------------|------------------|------------------|------------------|
| A100 04 051 2016 | 2,0              | A100 04 051 4514 | 4,5              |
| A100 04 051 2518 | 2,5              | A100 04 051 5016 | 5,0              |
| A100 04 051 3016 | 3,0              | A100 04 051 5522 | 5,5              |
| A100 04 051 3218 | 3,2              | A100 04 051 6020 | 6,0              |
| A100 04 051 3220 | 3,5              | A100 04 051 7025 | 7,0              |
| A100 04 051 4016 | 4,0              |                  |                  |


### Pelvic Traction Screw



| SS               | DIAMETER<br>(mm) |
|------------------|------------------|
| A602 00 010 3232 | 3,2              |
| A602 00 010 4545 | 4,5              |
| A602 00 010 8016 | 8,0              |
| A602 00 010 8020 | 8,0              |


# Wires, Pins, Staples & Drill Bits

## Epiphyseal Staple Stainless Steel




| SS               | A LENGTH (mm) | B LENGTH (mm) | DIAMETER (mm) |
|------------------|---------------|---------------|---------------|
| A100 03 110 0015 | 19            | 16            | 1,5           |
| A100 03 110 0020 | 19            | 22            | 2,0           |
| A100 03 110 0025 | 19            | 26            | 2,5           |

## Coventry Staple Stainless Steel




| SS               | A LENGTH (mm) | B LENGTH (mm) | C LENGTH (mm) | DIAMETER (mm) |
|------------------|---------------|---------------|---------------|---------------|
| A100 03 111 0140 | 33            | 25            | 14            | 2,5           |
| A100 03 111 0050 | 33            | 25            | 10            | 2,5           |
| A100 03 111 0100 | 33            | 25            | 5             | 2,5           |

## Thorny Epiphyseal Staple Titanium



| TI               | A LENGTH (mm) | B LENGTH (mm) |
|------------------|---------------|---------------|
| A100 03 115 0650 | 20            | 6,5           |
| A100 03 115 0800 | 20            | 8             |
| A100 03 115 0115 | 20            | 10            |
| A100 03 115 0150 | 20            | 15            |

## Thorny Coventry Staple Titanium



| TI               | A LENGTH (mm) | B LENGTH (mm) | C LENGTH (mm) |
|------------------|---------------|---------------|---------------|
| A100 03 111 0050 | 33            | 25            | 5             |
| A100 03 111 0100 | 33            | 25            | 10            |
| A100 03 111 0140 | 33            | 25            | 14            |
| A100 03 111 0200 | 33            | 25            | 16            |

# Wires, Pins, Staples & Drill Bits

## Thorny Washer



| SS               |
|------------------|
| A600 12 032 0000 |

## Tibial Bolt



| SS               | DIAMETER (mm) |
|------------------|---------------|
| A100 03 135 1000 | 3             |
| A100 03 130 0800 | 4             |

## Knowless Hip Pin



| SS               | DIAMETER (mm) |
|------------------|---------------|
| A100 14 051 4070 | 3,0           |
| A100 14 051 8075 | 3,2           |
| A100 14 051 5060 | 4,0           |
| A100 14 051 4585 | 4,5           |
| A100 14 051 5016 | 5,0           |

## Wires, Pins, Staples & Drill Bits

### Drill Bit



| SS               | DIAMETER (mm) |
|------------------|---------------|
| A101 03 090 0020 | 1,2           |
| A101 03 091 0027 | 1,5           |
| A101 03 091 0028 | 1,8           |
| A101 03 090 2012 | 2,0           |
| A101 03 090 2212 | 2,2           |
| A101 03 091 0022 | 2,5           |
| A101 03 090 2712 | 2,7           |
| A101 03 092 0030 | 3,0           |
| A101 03 090 3212 | 3,2           |
| A101 03 090 3512 | 3,5           |
| A101 03 090 0045 | 3,7           |
| A101 03 092 0040 | 4,0           |
| A101 03 090 4512 | 4,5           |
| A101 03 090 4028 | 5,0           |
| A101 03 090 0032 | 5,5           |
| A101 03 091 0045 | 6,0           |

### Gigli Saw



| SS               | LENGTH (mm) |
|------------------|-------------|
| A400 09 111 1030 | 300         |
| A400 09 111 1040 | 400         |
| A400 09 111 1050 | 500         |
| A400 09 111 1060 | 600         |




**STERILIZATION CONTAINERS  
AND ACCESSORIES**



**ORTHOPEDIC SURGICAL  
INSTRUMENTS**



**ORTHOPEDIC  
INSTRUMENT CONTAINERS**



**ORTHOPEDIC SURGERY  
INSTRUMENT SETS**



**PERFORATED  
SHEET BASKETS**



ORTHOPAEDICS & MEDICAL DEVICES



### AYSAM | Bafra

Doğankaya Mahallesi, Organize Sanayi Bölgesi,  
1.Cadde, No:51/1 Bafra, SAMSUN/TÜRKİYE



### AYSAM | Tekkeköy

Organize Sanayi Sitesi Erdoğan Cebeci Bulvarı No. 34  
Kutlukent Tekkeköy, SAMSUN/TÜRKİYE



### AYSAM | Istanbul

Molla Gürani Mah. Molla Gürani Cad. No:59/9  
Fındıkzade Fatih, ISTANBUL/TÜRKİYE



### AYSAM | Warehouse

Kirazlık Mah 1034 Cd. No:14  
Tekkeköy, SAMSUN/TÜRKİYE



Manufacturer of Orthopaedics & Medical Devices  
quality & trust



[aysam.com.tr](http://aysam.com.tr)